VILLACRYL HARD INSTRUCTIONS FOR USE - EN

I. Product description

Villacryl Hard: Hard denture relining material.

Processing parameters: necessary parameters for the proper processing of the materials are given in further section of IFU.

Indications for use

Villacryl Hard: relining of acrylic removable complete and partial dentures directly in the patient's mouth or indirectly in the dental laboratory based on the impression on the model.

Clinical Benefits:

Restoring lost biting and chewing functions to users of dental prostheses that have lost their retention to the denture foundation.

Target group or target groups of patients and intended users of the devices

Users: Professionals in the dental technology industry.

Qualified dentists

Qualified dental technicians producing dental prostheses.

Patients: People of different ages and sexes, denture users that have lost their retention to the denture foundation.

Contraindications

Do not use the products by people diagnosed with a lack of resistance to allergenic factors, including the ingredients of the product. In such situations, use a different material suitable for the patient.

For people sensitive to taste or smell sensations, with gag reflexes, it is recommended to perform indirect relining in the dental laboratory. Do not use the product by pregnant women and nursing mothers in the group of dental technicians.

Do not reline in case of mucositis.

Product form

Separate: a mixture of powder in plastic containers, liquid and etchant in bottles.

The product is available in the form of a set: powder + liquid + etchant + hardener

Chemical composition

Powder: acrylic polymer, barbituric acid derivative, pigments, benzoyl peroxide, silica

Fluid: acrylic monomers, N,N dimethyl-p-toluidine (DMPT), hydroquinone

Etching agent: acetone, ethyl acetate

Product compatibility

For use with acrylic dentures.

Safety instructions

Warnings

H225 Highly flammable liquid and vapor.

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation

H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled

H335 May cause respiratory irritation.

H336 May cause drowsiness or dizziness

H411 Toxic to aquatic life with long-lasting effects

EUH066: Repeated exposure may cause skin dryness or cracking

The liquid is sensitive to temperature and UV radiation. An uncontrolled polymerization reaction may occur.

Precautions

In case of any allergic reaction, stop work immediately and consult a dermatologist.

It is absolutely recommended that the patient rinse his mouth with water after relining, directly in the mouth.

In the case of patients who are sensitive to the components of the product, special care should be taken, do not reline directly in the patient's mouth, and in special cases use a different material.

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303 + P361 + P353 IF ON SKIN: Remove/ Take off immediately all contaminated clothing and rinse skin with water/shower.

P304 + P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Contact an ophthalmologist if advisable in case of irritation.

P312 Call a POISON CENTER or doctor/physician if you feel unwell.

P403 + P233 Store in a well-ventilated place. Keep container tightly closed.

P235 Keep cool.

P501 Dispose of the contents/container to an authorized consignee of the waste.

It is recommended to use gloves made of HDPE polyethylene for hand protection, a filtering fluidcompatible face mask for respiratory protection, goggles and workwear.

Do not eat or drink while working with the material. In the event of accidental ingestion of an ingredient of the product, consult a doctor immediately.

Do not pour the liquid down the sanitation system. Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

Information to be provided to patients / prosthesis users:

To reduce the residual monomer, it is advisable during the first three days, soak the relined prosthesis in water when not in use.

If any symptoms of an allergic reaction or irritation occur, discontinue using the prosthesis immediately and consult a doctor.

The color of the relined layer after a few months of use of the denture may change due to the consumption of food and drinks, pharmaceutical drugs or the use of chemical agents for disinfection and cleaning of the prosthesis. It is a natural phenomenon; it does not pose a threat to the patient and does not reduce the usefulness of the prosthesis.

Strong oxidizing compounds contained in cleaning agents and mouthwashes may change the color of the prosthesis after polymerization.

After each meal, the prosthesis should be washed with a soft brush under cool, running water.

After each meal denture should be cleaned with soft brush under cool tap water. For cleaning can be used unscented soap, various preparations for cleaning and caring for dental prostheses that do not contain abrasives. It is recommended to wash the dentures from time to time in a vibration washer filled with a special cleaning liquid for dentures.

Report any damage of the denture to the dentist. Repairs may be only performed by a dental technician. If a piece of damaged dentures is swallowed, consult a doctor immediately. The material after polymerization is not visible on X-ray pictures.

Side effects

Dental technicians/dentists

Irritation, redness or symptoms of hypersensitivity, including the occurrence of anaphylactic shock, in case of allergy to any of the product ingredients.

Patients

Direct relining technique:

A slight burning sensation in the mouth.

Irritation, inflammation and allergic reactions caused by the contact of the chemical components of the product with the mucosa.

In people sensitive to taste and smell sensations, gag reflexes may occur.

Material after polymerization:

Oral irritation, inflammation and allergic reactions caused by residual monomer washed out of the relined denture during use.

Irritation and inflammation of the oral cavity caused by fungi and bacteria adhering to the prosthesis underneath the material.

Storage conditions

Product is sensitive to temperature and UV radiation. Protect against light and UV radiation.

Store the product components in their original, tightly closed containers, in a dark and well-ventilated room, at a temperature of 5-25 °C , out of the reach of third parties.

Waste disposal

Dispose in accordance with applicable regulations of economic activity waste management.

Cardboard packages, labels, instructions: non-hazardous segregated municipal waste - paper.

Powder: non-hazardous mixed municipal waste

Powder packaging: non-hazardous segregated municipal waste – plastic

The liquid and its packaging – hazardous waste, comply with the regulations on hazardous waste management (storage, transport and disposal) and hand it over to a professional hazardous waste consignee.

How to use (Step by step instructions)

Initial steps:

Read the instructions carefully before using the product.

Make sure the patient for whom the prosthesis is reline has not been diagnosed with the lack of resistance to allergenic factors.

For people sensitive to taste or small sensations, with gag reflexes, it is recommended to perform indirect relining in the dental laboratory.

Before reusing the product, check if the fluid retains its properties. In the event of visible traces of density change, cloudiness or discoloration, stop using the product.

Preparation:

Cleaning and disinfecting dentures before repair

Remove any plaque and scale from the surface of the denture. Use protective gloves and masks. **Denture Cleaner by Everall7** is recommended for cleaning the denture. Wash the cleaned denture with warm water.

Product preparation:

Shake the powder and liquid containers 3 times as shown on the package and product labels to help disperse the powder in the monomer. Open the powder and liquid containers as shown on the packaging and product labels.

I. Making direct relining in dental office

Preparation: Remove the mucosal part of the denture with a bur. When a thick base layer is required, channels should be made to improve the distribution of the material on the surface of the prosthesis. Using a wick on prosthesis with Villacryl Hard Etchant and wait 2-3 minutes.

Measuring and mixing: Pour 1.5 ml (1.30 g) of Villacryl Hard liquid into the attached vessel, then add 1 scoop of powder (2.35 g). Mix gently for about 10 seconds.

Application: Protect the mucosa in the patient's mouth with petroleum jelly. Immediately after mixing the ingredients, apply the mass to the prepared areas on the denture. Instruct patient to moisten their lips with saliva to prevent Villacryl Hard from sticking to the tissues when inserting the denture into the mouth.

Placement in the patient's mouth: Place the prosthesis in the mouth maximum 90 seconds after mixing. The patient should perform the Herbst test movements as if taking a functional impression. Keep the prosthesis in this position for 2 minutes. During this time, patient may feel a slight burning sensation in the mouth. Carefully remove the denture from the mouth and remove any excess material with a scalpel. After making the relining, the patient should rinse the mouth thoroughly with water.

Polymerization: After removing the relined denture from the mouth, place it in a vessel filled with water at approx. 50 °C for approx. 30 minutes, adding the Hardener solution to the vessel with water (0.5 scoops of approx. 2 g per glass of water approx. 200 ml). Better results are obtained by using the above procedure by placing the vessel with water and the prosthesis in a pressure polymer at 2 bar.

Processing: After the material is completely cured, the denture should be processed and polished with burs and rubbers, intended for acrylic materials.

Reduction of the residual monomer level: After relining, place the denture back in water at 50 ° C for 30 minutes. It is advisable for the patient to soak the relined prosthesis in water during the first three days when not in use.

II. Making indirect relining in dental laboratory

Preparation: After taking the functional impression on the denture with silicone mass, use **Stodent III** to cast the stone model. Fix the obtained model with technical silicone in the fixator (verticulator) in order to maintain the height of the central occlusion set by the dentist during the procedure. Remove the impression material and prepare the base so that the future layer of Villacryl Hard of the same thickness. Apply with a brush the Villacryl Hard Etching fluid twice to the mucosal part of denture and wait until it evaporates. With a brush cover the scaled flask **twice with Iso-sol** - isolating liquid. Wait for the creation of isolation membrane.

Measuring and mixing: Pour 1.5 ml (1.30 g) of Villacryl Hard liquid into the attached vessel, then add 1 scoop of powder (2.35 g). Mix gently for about 10 seconds.

Application: Immediately apply the acrylic mass on the prepared surface of the prosthesis. Place the denture on the isolated model, press down and close the fixator (verticulator). Remove excess mass with a scalpel.

Polymerization: Place the denture with the model in a pressure polymerizer for 30 minutes at approx. 50 °C by adding the Hardener solution to a vessel with water (0.5 scoops, approx. 2 g per glass of water, approx. 200 ml) and pressurize to 2 bar.

Final processing: Use standard method of treatment of denture's border with burs and **Everall's Pumeks** and **Pasta polerska**.

Reduction of the residual monomer level: After relining, place the denture again in water at 50-60 ° C for 30 minutes.

It is recommended to soak the relined denture in water for a minimum of 24 hours after polymerization before handing it over to the patient. It is advisable for the patient to soak the relined prosthesis in water during the first three days when not in use.

Hygiene, sterilization, disinfection

Product

The product is intended for multiple use. It is not supplied sterile and does not require sterilization or disinfection prior to use. To protect the product against contamination, apply the principles of good laboratory practices in dental laboratories.

Relined prosthesis

There is a risk of transmission of infection.

After removing the denture from the patient's mouth for relining, wash and disinfect.

After relining, the prosthesis should be washed and disinfected again before being placed in the patient's mouth.

For the disinfection of dentures made of materials, it is recommended to use disinfectants for disinfecting dental prostheses, following the procedures recommended by the manufacturer.

It is also possible to use agents based on at least 65% ethyl or propyl alcohol. The denture should be sprayed with a disinfectant solution, and then gently wiped to dryness with a paper tissue. You can also use ready-to-use disposable wipes soaked in a disinfecting solution. Additionally, the denture should be washed with warm soapy water.

Batch number and expiry date - see imprint on product labels and packaging.

The expiry date refers to properly stored unopened packaging. Do not use the product after the expiry date.

Additional information

Reporting of serious incidents: Any serious incidents related to the products must be reported to the manufacturer and the competent authority in accordance with local regulations.

Data necessary to identify the product and the contents of the package

Reference number, product name, color, package size, batch number and expiry date - imprint on the product packaging.

Package content - see table below.

Assortment*

Color shades

Villacryl Hard	
pink	

Packaging

	Powder + liquid + etchant + hardener kit
Villacryl Hard	60 g powder + 40 ml liquid +10 ml etchant + 40 g hardener

Problem prevention (1)

See website: www.everall7.pl

MD

Please send any comments regarding our products to our address.

Information provided in any way, even during the demonstration, does not invalidate the instructions for use. Users are responsible for checking that the product is suitable for the intended use. The manufacturer is not liable for damages, including to third parties, resulting from non-compliance with the instructions or non-use. The manufacturer's liability is in any case limited to the value of the products delivered.

MD	The symbol means a medical device as defined in EU Regulation 2017/745.
-----------	---

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