Master all your challenges

initial Zirconia Coloring Liquid


- EN Intro Kit Instruction For Use
- DE Intro Kit Gebrauchsanweisung
- FR Intro Kit Mode d'emploi
- IT Intro Kit Istruzioni per l'uso
- ES Intro Kit Instrucciones de uso
- PL Intro Kit Instrukcja obsługi

JUNE 2024

Since 1921
100 years of Quality in Dental

EN Intended Use

Intended for use by professional dental technicians.

Water-based coloring liquids for the complete or partial coloring of pre-sintered zirconia restorations made from zirconia milling blanks for the manufacturing of permanent dental restorations for adults.

Indications

Individualization of anatomically reduced and fully anatomical (monolithic) crowns and bridges made of zirconium oxide before sintering.

Contraindications

Intolerance to contained components.

Adverse effects

We are not aware of any risks or adverse effects.

Warnings

1. Causes serious eye damage.
2. May cause an allergic skin reaction.
3. Processing of GC Initial Zirconia Disk produces dust which can lead to eye and skin irritation and may cause lung damage.
Always ensure that the extraction unit of the machine is in good working order. Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
4. Do not use autoclave to sterilize the framework.
5. Dispose of contents/container in accordance with local/regional/national/international regulations.

Storage

Liquid can only be stored in the original bottle. Keep it separated from edibles and ensure the bottle is tightly sealed. Avoid storing it with alkalis (lyes). Close all bottles immediately after use. For reproducible results, shake each liquid before use.

Note:

The "Initial Zirconia Coloring Liquid" are sensitive to a slight color shift within their shelf life during storage or use. This however does not affect the color of the final sintered restoration.

Packaging

GC Initial Zirconia Coloring Liquid, SET, 750g

- 4 Body colors: Body A / Body B / Body C / Body D
- 4 Cervical colors: Cervical A / Cervical B / Cervical C / Cervical D
- 3 Effect colors: Blue, Grey, Pink
- 1 Neutral
- 4 Marking Liquids: Red, Blue, Yellow, Green

GC Initial Zirconia Coloring Liquid, INTRO KIT, 250g

1x Body A, 1x Cervical A, 1x Effect Grey, 1x Neutral, 1x Marking Liquid Red, 1x Marking Liquid Blue

EN Instructions For Use

One Liquid – two techniques

1. Painting Technique – individual & detailed procedure

- Application of the different Body – Cervical – Effect colors by painting
Shade-assignment tables mixing information or immersion timings, please refer to the dedicated Technical Manual - GC Initial Zirconia Coloring Liquid under www.gc.dental/europe
- Drying
Dry the colored zirconia restoration before the sintering process at 80 °C for 60 minutes in a dedicated small oven with a heating lamp. To ensure a uniform color result, the dental restoration should be dried immediately after the color application. Position the crowns and bridges on their occlusal surface while drying.
- Sintering Process
Sinter the colored and dried zirconia restoration according to the zirconia manufacturer's recommendation.
- Final aesthetics can be achieved with the GC Initial IQ ONE SQIN concept.

2. Immersion Technique – fast & easy procedure

- Immersion in coloring liquid
Shade-assignment tables mixing information or immersion timings, please refer to the dedicated Technical Manual - GC Initial Zirconia Coloring Liquid under www.gc.dental/europe
- Drying
Dry the colored zirconia restoration before the sintering process at 80 °C for 60 minutes in a dedicated small oven with a heating lamp. To ensure a uniform color result, the dental restoration should be dried immediately after the color application. Position the crowns and bridges on their occlusal surface while drying.
- Sintering Process
Sinter the colored and dried zirconia restoration according to the zirconia manufacturer's recommendation.
- Final aesthetics can be achieved with the GC Initial IQ ONE SQIN concept.

For more tips and tricks please consult the Initial Zirconia Coloring Liquid Technical Manual on our dedicated website www.gc.dental/europe

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EN Medical Device Regulation Information

Composition

Erbiumchloride Hexahydrate < 8, Ironchloride Hexahydrate < 6, Praseodymiumacetate Hydrate < 6, Manganesulfate Monohydrate < 0,3, Neodymiumchloride < 8, Copper (II) sulfate Pentahydrate < 0,1, Chromium acetate < 0,1, Polyethylenglycol 6-9, Hydrochloric Acid 0.2-0.4, Water > 80

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gceurope.com> or for The Americas: <http://www.gcamerica.com>. They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at klema.qm@klema.at.

Undesired effects - Reporting: If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: klema.qm@klema.at.

In this way you will contribute to improve the safety of this product.

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EU: KLEMA Dentalprodukte GmbH
ALSIP, IL 60803 U.S.A.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
GC ASIA DENTAL PTE. LTD.
GC AUSTRALASIA DENTAL PTY.LTD