**PRINT SPECIFICATION PAGE ONLY.  
DO NOT PRINT THIS PAGE IN THE FINAL ARTWORK.**

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2512256PE32  
2020-03-10

ETCHING GEL  
Etching Gel  
for Enamel and Dentin

**DEVICE DESCRIPTION**  
Etching Gel is a 37% phosphoric acid gel used for etching enamel and dentin to prepare the tooth surface prior to the application of a bonding agent. The non-slumping consistency of the blue-coloured Etching Gel enables precise application.

Etching Gel is available in syringes with single-use application tips.

#### INDICATIONS FOR USE

Etching of enamel surface or dentin before the preparation of a composite or bonded amalgam filling, before the preparation of a sandwich composite filling with a glass ionomer cement; before adhesive fixation of an aesthetic facet, either ceramic or composite, before the fixation of an inlay or onlay, the root pin, a crown or a bridge with a composite cement.

#### CONTRAINDICATIONS

None Known

#### WARNINGS

1. Never use excessive force on the syringe plunger to dispense etch, which may lead to over dispensing of the material. If resistance is felt, replace the current application tip with a new tip and check flow extrorally.  
2. Contains phosphoric acid. Avoid contact with the skin, eyes and soft tissue. In case of contact, flush immediately with water. Get medical attention for the eyes. Do not take internally.  
3. This product shall not be placed in direct contact with the pulp tissue.

#### PRECAUTIONS

1. The use of personal protective equipment and a contact-free technique is recommended when handling and working with the Etching Gel.  
2. The use of a rubber dam is recommended to protect the surrounding soft tissues.  
3. Application tips are for single patient use only to prevent cross-contamination between patients.  
4. This product is intended to be used by qualified dental professionals on a general patient population.

#### ADVERSE REACTIONS

None Known

#### INSTRUCTIONS FOR USE

1. Use a rubber dam to isolate the oral cavity.  
2. Securely attach the application tip to the syringe and check the flow extrorally.  
3. The etching gel is applied to the enamel edges directly with the tip and left to work for 20 seconds.  
4. If dentin is to be etched ("total-etch"), the etching gel is applied to the dentin 15-20 seconds later so that the overall working time does not exceed 30 seconds.  
5. The etching gel is removed with a sufficient stream of water (for 15 seconds) and the cavity surface is dried with a stream of air which is free of oil particles, so that the enamel edges become white and the dentin is not overried, i.e., does not lose its natural moisture.  
6. If the etched surface is contaminated, e.g., with saliva, the etching must be repeated for 10 seconds and the surface washed in water and dried in the usual manner.  
7. The selected adhesive system is applied according to the manufacturer's instructions.

#### CLEANING AND DISINFECTION INSTRUCTIONS

The used application tips should be removed and disposed. After closing the syringe with a stopper, disinfect its surface with a solution of a common surface disinfectant.

#### STORAGE, SHELF LIFE AND DISPOSAL

To be stored in a dry place, at a temperature of 5 - 25 °C, away from heat sources. Refer to the outer package for the expiration date. Do not use after the expiration date. For proper disposal, always follow your local and national regulations.

#### VIGILANCE

If an adverse medical event occurs, report it to the manufacturer and to the competent authority for that country.

#### DISCLAIMER

The instructions provided above have been validated by SpofaDental as being CAPABLE of preparing a medical device, including reusable ones. It remains the responsibility of the dental professional / user to ensure

that the reprocessing actually performed by equipment, materials and personnel is provided for the instructions provided above in order to achieve the desired result. Any deviation by the dental professional / user from the instructions provided above will be done at your own risk; the SpofaDental will not be able to address any requests for refund or exchanges under guarantee of any products that have not been handled or reprocessed in accordance with the above instructions.

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