# Bluephase® EasyCure

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# Appendix

Bluephase EasyCure is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 3th (see the following tables). Bluephase EasyCure is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

## Guidance and manufacturer's declaration – electromagnetic emission

The following tables are guidelines according to the medical standard IEC 60601-1-2.

Bluephase EasyCure is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase EasyCure should assure that it is used in such an environment.

| Emission test  | Compliance | Electromagnetic environment - guidance  |
|--|------------|---|
| RF emissions CISPR 11                                  | Group 1    | Bluephase EasyCure uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.   |
| RF emissions CISPR 11                                  | Class B    | Bluephase EasyCure is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2                       | Class A    | N/A   |
| Voltage fluctuations / flicker emissions IEC 61000-3-3 | Complies   | N/A   |

Table: According to IEC 60601-1-2 3th Edition

## Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase EasyCure is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase EasyCure should assure that it is used in such an environment.

| Immunity test   | IEC 60601 test level  | Compliance level  | Electromagnetic environment – guidance   |
|---|---|---|--|
| Electrostatic discharge (ESD) IEC 61000-4-2   | ± 6 kV contact<br>± 8 kV air  | ± 6 kV contact<br>± 8 kV air  | Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.  |
| Electrical fast transient/burst IEC 61000-4-4   | ± 2 kV for power supply lines<br>± 1 kV for input/output lines  | ± 2 kV for power supply lines<br>± 1 kV for input/output lines  | Mains power quality should be that of typical commercial or dental environment.  |
| Surge IEC 61000-4-5   | ± 1 kV line - line<br>± 2 kV line - earth   | ± 1 kV line - line<br>no prot. earth  | Mains power quality should be that of typical commercial or dental environment.  |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | <5% $U_T$<br>(>95% dip in $U_T$ ) for 0.5 cycle<br>40% $U_T$<br>(60% dip in $U_T$ ) for 5 cycles<br>70% $U_T$<br>(30% dip in $U_T$ ) for 25 cycles<br><5% $U_T$<br>(>95% dip in $U_T$ ) for 5 sec | <5% $U_T$<br>(>95% dip in $U_T$ ) for 0.5 cycle<br>40% $U_T$<br>(60% dip in $U_T$ ) for 5 cycles<br>70% $U_T$<br>(30% dip in $U_T$ ) for 25 cycles<br><5% $U_T$<br>(>95% dip in $U_T$ ) for 5 sec | Mains power quality should be that of typical commercial or dental environment. If the user of Bluephase EasyCure requires continued operation during power mains interruptions, it is recommended that Bluephase EasyCure be powered from an uninterruptible power supply or battery. |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8   | 3 A/m   | 30 A/m  | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or dental environment.  |

Table: According to IEC60601-1-2 3th Edition

NOTE:  $U_T$  is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase EasyCure is intended for use in the electromagnetic environment specified below. The customer or the user of Bluephase EasyCure should assure that it is used in such an environment.

| Immunity test                 | IEC 60601 test level                    | Compliance level | Electromagnetic environment – guidance   |
|-------------------------------|---|------------------|--|
|                               |   |                  | Portable and mobile RF communications equipment should not be used closer to any part of Bluephase EasyCure, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.   |
|                               |   |                  | Recommended separation distance  |
| Conducted RF<br>IEC 61000-4-6 | 3 V <sub>rms</sub><br>150 kHz to 80 MHz | 10 V             | $d = 0.35 \sqrt{P}$  |
| Radiated RF<br>IEC 61000-4-3  | 3 V/m<br>80 MHz to 2.5 GHz              | 10 V/m           | $d = 0.35 \sqrt{P}$ 80 MHz to 800 MHz<br>$d = 0.70 \sqrt{P}$ 800 MHz to 2.5 GHz  |
|                               |   |                  | Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).<br><br>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup><br><br>Interference may occur in the vicinity of equipment marked with the following symbol:<br><br> |

Table: According to IEC60601-1-2 3th Edition

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Bluephase EasyCure is used exceeds the applicable RF compliance level above, Bluephase EasyCure should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Bluephase EasyCure.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strength should be less than 10 V/m. NOTE: U<sub>i</sub> is the a.c. mains voltage prior to application of the test level.

## Recommended separation distances between portable and mobile RF communications equipment and "Bluephase EasyCure"

Bluephase EasyCure is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of Bluephase EasyCure can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Bluephase EasyCure as recommended below, according to the maximum out-put power of the communication equipment.

| Rated maximum output power of transmitter (W) | Separation distance according to frequency of transmitter (m) |  |  |
|---|---|--|--|
|   | 150 kHz to 80 MHz<br>$d = 0.35 \sqrt{P}$                      | 80 MHz to 800 MHz<br>$d = 0.35 \sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 0.7 \sqrt{P}$ |
| 0.01  | 0.035   | 0.12                                     | 0.23                                     |
| 0.1   | 0.11  | 0.38                                     | 0.73                                     |
| 1   | 0.35  | 1.2                                      | 2.3                                      |
| 10  | 1.1   | 3.8                                      | 7.3                                      |
| 100   | 3.5   | 12                                       | 23                                       |

Bluephase EasyCure is EMC-tested in conformity with the requirements of 60601-1-2:2014+AMD1:2020 (see the following tables). Bluephase EasyCure is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

### Guidance and manufacturer's declaration – electromagnetic emission

Bluephase EasyCure is intended for use in a professional healthcare environment. The customer or the user of Bluephase EasyCure should assure that it is used in such an electromagnetic environment specified below.

| Emission test  | Compliance | Electromagnetic environment - guidance  |
|--|------------|---|
| RF emissions CISPR 11                                  | Group 1    | Bluephase EasyCure uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.   |
| RF emissions CISPR 11                                  | Class B    | Bluephase EasyCure is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2                       | Class A    |   |
| Voltage fluctuations / flicker emissions IEC 61000-3-3 | Complies   |   |

Table: According to IEC60601-1-2:2014+AMD1:2020, table 2

### Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase EasyCure is intended for use in a professional healthcare environment. The customer or the user of Bluephase EasyCure should assure that it is used in such an electromagnetic environment specified below.

| Immunity test  | IEC60601-1-2 ed. 4.1 test level   | Compliance level   | Electromagnetic environment – guidance  |
|--|---|--|---|
| Electrostatic discharge (ESD) IEC 61000-4-2                              | ± 8 kV contact<br>± 15 kV air   | ± 8 kV contact<br>± 15 kV air  | Floors should be concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%. |
| Radiated RF EM fields IEC 61000-4-3                                      | 3 V/m<br>80 MHz – 2.7 GHz<br>80% AM at 1 kHz  | 10 V/m<br>80 MHz – 2.7 GHz<br>80% AM at 1 kHz  | Bluephase EasyCure fulfills the test levels for a home healthcare environment.  |
| Proximity fields from RF wireless communications equipment IEC 61000-4-3 | 27 V/m, 385 MHz, Pulse mod. 18 Hz<br>28 V/m, 450 MHz, Pulse mod. 18 Hz<br>9 V/m, 710/745/780 MHz, Pulse mod. 217 Hz<br>28 V/m, 810/870/930 MHz, Pulse mod. 18 Hz<br>28 V/m, 1720/1845/1970 MHz, Pulse mod. 217 Hz<br>28 V/m, 2450 MHz, Pulse mod. 217 Hz<br>9 V/m, 5240/5500/5785MHz, Pulse mod. 217 Hz | 27 V/m, 385 MHz Pulse mod. 18 Hz<br>28 V/m, 450 MHz Pulse mod. 18 Hz<br>9 V/m, 710/745/780 MHz Pulse mod. 217 Hz<br>28 V/m, 810/870/930 MHz Pulse mod. 18 Hz<br>28 V/m, 1720/1845/1970 MHz Pulse mod. 217 Hz<br>28 V/m, 2450 MHz Pulse mod. 217 Hz<br>9 V/m, 5240/5500/5785MHz Pulse mod. 217 Hz |   |

|  |   |   |   |
|--|---|---|---|
| Rated power frequency magnetic field IEC 61000-4-8 | 30 A/m  | 100 A/m   |   |
| Proximity magnetic fields IEC 61000-4-39           | 8 A/m, 30 kHz<br>CW<br>65 A/m, 134.2 kHz<br>Pulse mod. 2.1 kHz<br>7.5 A/m, 13.56 MHz<br>Pulse mod. 50 kHz | 8 A/m, 30 kHz<br>CW<br>65 A/m, 134.2 kHz<br>Pulse mod. 2.1 kHz<br>7.5 A/m, 13.56 MHz<br>Pulse mod. 50 kHz | For information only:<br>Bluephase EasyCure does not contain magnetic sensitive parts or circuitry (acc. 8.11a) |

Table: According to IEC60601-1-2:2014+AMD1:2020, table 4

### Warning

Portable wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkie etc. can affect the Bluephase EasyCure and should be kept at least a distance of 30 cm (12 inches) to any part of Bluephase EasyCure.

### Guidance and manufacturer's declaration – electromagnetic immunity

Bluephase EasyCure is intended for use in a professional healthcare environment. The customer or the user of Bluephase EasyCure should assure that it is used in such an electromagnetic environment specified below.

| Immunity test   | IEC60601-1-2 ed. 4.1 test level   | Compliance level  | Electromagnetic environment – guidance   |
|---|---|---|--|
| Electrical fast transients /bursts IEC 61000-4-4          | ± 2 kV<br>100 kHz repetition frequency  | ± 2 kV<br>100 kHz repetition frequency  |  |
| Surges<br>Line-to-line<br>IEC 61000-4-5                   | ± 0.5 kV, ± 1 kV  | ± 0.5 kV, ± 1 kV  |  |
| Surges<br>Line-to-ground<br>IEC 61000-4-5                 | ± 0.5 kV, ± 1 kV, ± 2 kV  | no protection earth   |  |
| Conducted disturbances induced by RF fields IEC 61000-4-6 | 3 V<br>0.15 MHz – 80 MHz<br>6 V in ISM bands<br>between 0.15 MHz and 80 MHz<br>80 % AM at 1 kHz   | 10 V<br>0.15 MHz – 80 MHz<br>10 V in ISM and amateur bands<br>between 0.15 MHz and 80 MHz<br>80 % AM at 1 kHz   | Bluephase EasyCure fulfills the test levels for a home healthcare environment. |
| Voltage dips IEC 61000-4-11                               | 0 % U <sub>T</sub> : 0.5 cycles<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br><br>0 % U <sub>T</sub> : 1 cycle<br>and<br>70 % U <sub>T</sub> : 25/30 cycles<br>Single phase: at 0° | 0 % U <sub>T</sub> : 0.5 cycles<br>At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br><br>0 % U <sub>T</sub> : 1 cycle<br>and<br>70 % U <sub>T</sub> : 25/30 cycles<br>Single phase: at 0° |  |
| Voltage interruptions IEC 61000-4-11                      | 0 % U <sub>T</sub> : 250/300 cycles   | 0 % U <sub>T</sub> : 250/300 cycles   |  |

Table: According to IEC 60601-1-2:2014+AMD1:2020, table 5

NOTE: U<sub>T</sub> is the a.c. mains voltage prior to application of the test level.

[ivoclar.com](http://ivoclar.com)

**Making People Smile**