Manufactured By: Pentron Clinical
1717 West Collins Avenue Orange, CA 92867 USA
800-551-0283 / 714-516-7557 / www.pentron.com

Universal C&B Resin Cement

Distribution Facility:
Kerr Australia Pty. Ltd., Unit 10, 112-118 Talavera Road, North Ryde 2113, Australia
+61-29-870-7500

CE [REP] SpofaDental a.s.
Markova 238
CZ-506 01 Jičín
Phone: +420 493 583 111

CE
0086
8°C (46°F)
2°C (36°F)
Temperature Limit

DMC1639_Rev.02

Rx ONLY

INDICATIONS FOR USE
Cement-It is a luting cement bonding system for porcelain crowns, inlays, and onlays, and other dental restorations, including ceramics, metal alloys and composite restoratives.

CONTRAINDICATIONS

None known

WARNINGS

- Uncured methacrylate resin may cause contact dermatitis and harm the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.
- *37% Phosphoric Acid Etching Gel contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water; get medical attention for the eyes. Do not take internally.
- Bond-1™ One-Component Total Etch Bonding System contains ethanol and is flammable and volatile. Keep away from sources of ignition and excessive heat.

PRECAUTIONS

- Tips and applicator brushes are for single patient use only to prevent cross-contamination between patients.
- If using phosphoric acid to etch the enamel and dentin, never use excessive force on the syringe plunger to dispense etchant as it may lead to over dispensing of the material. If resistance is felt, replace the tip with a new one and check the flow extra-orally.
- Cement-It requires refrigeration and storage at 2-8°C (36-46°F).

ADVERSE REACTIONS

None known

INSTRUCTIONS FOR USE

Cementing Ceramic or Indirect Composite Crowns, Bridges, Inlays & Onlays

- Apply 37% Phosphoric Acid Etching Gel™ to the prepared enamel and dentin surfaces for 15 seconds. Rinse thoroughly for 15 seconds with water. It is important not to splash etchant into the eyes, onto soft tissue or skin.
- Dry with clean, oil-free air without desiccating the dentin. Do not allow the etched surface to become contaminated with saliva or blood during the bonding procedure.
- Apply two consecutive coats of Bond-1 within 20 seconds. Dry with a gentle stream of air to volatilize any residual ethanol solvent for approximately 10-15 seconds. The tooth should have a shiny appearance. If not, apply an additional coat and gently air dry for 10 seconds. Light cure for 10 seconds. **Optional: Etch the inner surface of the sandblasted metal restoration. Rinse and air dry.**
- Apply one coat of Bond-1 to the inside of the restoration. Remove any excess and residual ethanol with a jet of air for approximately 10-15 seconds.
- Place the Cement-It mixture

ENGLISH

Cement-It™

Universal C&B Resin Cement

DEVICE DESCRIPTION

Cement-It is the ideal choice for all your indirect bonding needs. This comprehensive, dual-cure cementation system delivers high bond strength, reduced sensitivity, and versatility to any bonded restoration. Cement-It is available in auto-mix syringes in these shades: A2, Opaqueous and Translucent.

INDICATIONS FOR USE
Cement-It is a luting cement bonding system for procelain crowns, inlays, and onlays, and other dental restorations, including ceramics, metal alloys and composite restoratives.

CONTRAINDICATIONS

None known

WARNINGS

- Uncured methacrylate resin may cause contact dermatitis and harm the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.
- *37% Phosphoric Acid Etching Gel contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water; get medical attention for the eyes. Do not take internally.
- Bond-1™ One-Component Total Etch Bonding System contains ethanol and is flammable and volatile. Keep away from sources of ignition and excessive heat.

PRECAUTIONS

- Tips and applicator brushes are for single patient use only to prevent cross-contamination between patients.
- If using phosphoric acid to etch the enamel and dentin, never use excessive force on the syringe plunger to dispense etchant as it may lead to over dispensing of the material. If resistance is felt, replace the tip with a new one and check the flow extra-orally.
- Cement-It requires refrigeration and storage at 2-8°C (36-46°F).

ADVERSE REACTIONS

None known

INSTRUCTIONS FOR USE

Cementing a Metal Restoration to Enamel and Dentin

- Apply 37% Phosphoric Acid Etching Gel™ to the enamel and dentin surfaces for 15 seconds. Rinse thoroughly for 15 seconds with water. It is important not to splash etchant into the eyes, onto soft tissue or skin.

- Dry with clean, oil-free air without desiccating the dentin. Do not allow the etched surface to become contaminated with saliva or blood during the bonding procedure.
- Apply two consecutive coats of Bond-1 within 20 seconds. Dry with a gentle stream of air to volatilize any residual ethanol solvent for approximately 10-15 seconds. The tooth should have a shiny appearance. If not, apply an additional coat and gently air dry for 10 seconds. Light cure for 10 seconds. **Optional: Etch the inner surface of the sandblasted metal restoration. Rinse and air dry.**
- Apply one coat of Bond-1 to the inside of the restoration. Remove any excess and residual ethanol with a jet of air for approximately 10-15 seconds.
- Place the Cement-It mixture

inside the restoration. Cement-It has a minimum working time of 45 seconds from mixing and a maximum setting time of 4.5 minutes from mixing.
6. Seat the restoration and apply even pressure with a bite stick, allowing excess cement to vent.
7. Remove any pre-polymerized excess with a small brush or cotton pellet.
After final setting, any remaining excess can be removed with the instruments of your choice.

SHELF LIFE AND DISPOSAL

Refer to the outer packaging for expiration date. Do not use after the expiration date. For proper disposal, always follow Federal or State laws.

The trademarks mentioned in these instructions may be the trademarks of their respective owners.

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