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OsteoBio® Collagenated membranes and feltis - English -

Class III Medical Device, collagenated barrier membranes of heterologous origin, resorbable, engineered to secure the bone grafts in guided tissue regeneration surgery.

Evolution: Origin: equine pericardium / pig parietal peritoneum. Slurdy, stable, soft, flexible, easy to hydrate, can be sutured or glued with fibrin sealant, it provides an effective barrier effect beneficial to underlying tissue regeneration. The resorption time of the material may vary, depending on the grafting site and its thickness, indicated on the label.
ESTIMATED RESORPTION TIME: Sid about 3/4 months; Fine about 3/4 months; X-Fine about 2 months.
CLINICAL INDICATIONS: Evolution is indicated in all surgical operations where the use of a membrane is required to guide the neoformation of vital and vascularized connective tissue and/or the filling or replacement of connective tissues.
DIRECTIONS FOR USE: open the package, take the membrane with sterile presses, cut it and it shape it at will and eventually soak it in sterile solutions; place it in the receiving site and fix it by suture or with screws or pins or other synthetic means, or glue it with fibrin sealant.
PACKAGING: sterile disposable blister.
SIDE EFFECTS: none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded.
CONTRAINDICATIONS: do not use in patients who have pathologies or conditions that contraindicate surgery or who have an acute or chronic infection (ostomyelitis). Not indicated for patients with immunodeficiency or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly or near large vessels or nerves. Do not use on poorly vascularized, necrotic or infected rearing sites. Do not use on sites that do not guarantee the stability and protection of the product after grafting. The product must not be used if the packaging is damaged or opened before use.
WARNING: This product has not been tested on pregnant women. The product is for single use only and, as such, should not be resterilized. The graft site must be properly washed and cleared, and the product must be well stabilized and protected.
The product reserved for surgeons with experience in biomaterial implants: strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on OsteoBio® devices available on the website [www.osteobiol.com](#), as well as the instructions herein. For every single application, the specialist surgeon has the responsibility to ensure proper clinical management of the biomaterials and the patient.

EXPIRY DATE: the product is valid for five years from the date of sterilization and the expiry date is indicated on the label.
STORAGE: The product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label.

DISPOSAL: the product does not contain hazardous substances and it may be disassembled and disposed of as regular municipal waste.

PRODUCT ORIGIN: the last letter on the label of the product code REF indicates the source of the material (S=porcine; E=equine).
Symbols used on the label:

☞ single use ☞ year and month of expiry ☞ oot number ☞ carefully read the instructions
☞ Sterile. Sterilization performed by irradiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed ☞ limitation of the temperature

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