HEAT CURE

- DE** Gebrauchsanweisung
- GB** Instructions for use
- FR** Mode d'emploi
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- PT** Instruções de uso
- NL** Gebruiksaanwijzing
- SE** Bruksanvisning
- DK** Brugervejledning
- GR** Οδηγίες χρήσης
- TR** Kullanma Talimatları
- RU** Инструкция по применению
- UA** Інструкція по застосуванню

66076999/12

KULZER
MITSUI CHEMICALS GROUP

Manufacturer:
Kulzer GmbH
Leipzig Straße 2
63450 Hanau (Germany)

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GB Instructions for use Meliodent® Heat Cure

Colour
See packaging

Composition
Heat curing polymer, powder (polymethylmethacrylate) and liquid (methylmethacrylate, dimethacrylate)

Application
Dental resin for making dentures. For fixed and removable prosthetic devices.

Pack and press technique
Dosage and processing
Hard plaster is recommended for preparing the model and for embedding. The recommended mixing ratio is 35 g : 14 ml (powder : liquid). Introduce the liquid to the mixing beaker, add the appropriate quantity of powder rapidly and mix for approx. 1 minute with a spatula. Meliodent Heat Cure reaches its packable consistency after approx. 10 minutes, the processing time is approx. 25 minutes. The processing time is extended or reduced depending on the ambient temperature.

Polymerisation pack and press technique
Place the flask in boiling water, switch the heat source off and leave for 15 mins. Then boil for 20 mins. Allow the flask to cool down slowly in the water bath after polymerization.

Injection technique
Dosage and processing
The mixing ratio is 35 g : 14 ml (powder : liquid). Place the liquid in the mixing cup, add the specified volume of powder rapidly, stir to form a homogeneous mass for approx. 1 minute and leave to swell for approx. 7 minutes. Then release the mass from the mixing cup, knead briefly by hand (wear protective gloves) and plug into the filling cylinder.

The waiting time before injection from the start of mixing at 23°C (73°F) is 9–10 minutes. The waiting times depend on the room temperature and mixed volume. Close the filling cylinder using the cylinder cover and screw into the Palajet unit. Fasten flask using the centring wheel. After pulling the pressure unit forward as far as the stop, the rocker arm is moved upwards with the filling cylinder inserted and the denture acrylic injected into the flask with the aid of the pressure piston. Keep the flask under pressure in the Palajet unit for 15 minutes (rocker arm up). To check the time, set the timer to 15 minutes.

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Version: 2017-07

Polymerisation injection technique
Polymerisation time in the pressure vessel is 30 minutes, water temperature 100°C (212°F), pressure 2 bar.

Note
MMA/PMMA-based denture materials always contain a small, permissible amount of residual monomer after fabrication. To reduce the risk of intolerance reaction, Kulzer recommends always storing prostheses (both new and repaired) in tepid water for at least 12 hours before placement.

Repairing and relining
Dentures made of Meliodent Heat Cure can be repaired with Meliodent Rapid Repair.

Storage advice
Do not use the material after the expiry date. Do not store above 25°C (77°F). Avoid direct sunlight. Always close containers tightly after use. Keep out of reach of children.

Hazard and safety precautions for liquid (methyl methacrylate)
Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Keep away from heat/sparks/open flames/hot surfaces. No smoking. Do not inhale dust/fumes/gas/mist/vapour/aerosols. Wear protective gloves/protective clothing/eye protection. If skin irritation or rash occurs: Get medical advice/attention. In case of fire: Use the following to extinguish: CO₂, sand, extinguishing powder. Keep container in a well-ventilated place. Keep container tightly closed.

Safety data sheets and more information are available on our website www.kulzer.com

© = registered trademark of Kulzer GmbH

Dated: 2017-07

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© = registrerat varumärke av Kulzer GmbH