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Bone and collagen matrix OsteoBio® - English -

Class III medical device, bone and collagen matrix of heterologous origin, for the replacement of the autologous bone represented by the matrix and the porous form.

Sp-Block: cancellous block-shaped bone matrix, sizes specified on the label.

C-Block: cancellous cylinder-shaped bone matrix, sizes specified on the label.

V-Block: cancellous wedge-shaped bone matrix, sizes specified on the label.

Dual-Block: cortico-cancellous preformed bone matrix, sizes specified on the label.

Cortical-Teck: cortical preformed bone matrix, sizes specified on the label.

CHARACTERISTICS: heterologous bone matrix and collagen, osteoconductive, resorbable, biocompatible, hygroscopic. The product is slightly radiopaque. PACKAGING: sterile disposable blister. CLINICAL INDICATIONS: adjuncts for the reconstruction or partial or total recovery of lost bone portions, such as fillers for non-infected, non-sclerotic and well vascularized bone defects. DIRECTIONS FOR USE: hydrate the blocks with sterile solutions and shape them under sterile conditions; adapt them to the receiving site, carefully freshened in order to ensure proper vascularization of the graft. Always stabilize the preformed in the grafting site using sterile synthetic devices (screws and pins). SIDE EFFECTS: none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded. CONTRAINDICATIONS: do not use in patients who have pathologies or conditions that contraindicate surgery or who have an acute or chronic infection (osteomyelitis). Not indicated for patients with immunodisregulation or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly on near large vessels or nerves. Do not use on poorly vascularized, necrotic or infected recipient sites. Do not use directly on the product after the date of sterilization and the expiry date. The product must not be used if the packaging is damaged or opened before use. WARNING: This product has not been tested on pregnant women. The product is for single use only and, as such, should not be re-sterilized. The graft site must be properly washed and cleaned, and the product must be well stabilized and protected. Product reserved for surgeons with experience in biomaterial implants: strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on OsteoBio® devices available on the website www.osteobiol.com, as well as the instructions herein. For every single application, the specialist surgeon has the responsibility to ensure proper clinical management of the biomaterials and the patient. EXPIRY DATE: the product is valid for five years from the date of sterilization and the expiry date is indicated on the label. STORAGE: The product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label. DISPOSAL: the product does not contain hazardous substances and it may be disassembled and disposed of as regular municipal waste. PRODUCT ORIGIN: the last letter on the label of the product code REF indicates the source of the material (S=porcine; E=equine). Symbols used on the label: single use year and month of expiry lot number carefully read the instructions Sterile. Sterilization performed by radiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed. limitation of the temperature

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