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OsteoBio® Collagenated membranes and felts – English -

Class III Medical Device, collagenated barrier membranes of heterologous origin, resorbable, in guided tissue regeneration surgery. Dermas obtained from derma of porcine origin, its texture and sturdiness promote stabilization and protection of the underlying graft, after careful stabilization and suturing of flaps. Resorbable and biocompatible. The resorption time of the material may vary, depending on the grafting site and its thickness, indicated on the label.

CLINICAL INDICATIONS: membrane engineered to protect reconstruction or recovery of lost tissue portions; collagenated filler; protection of non-infected, non-sclerotic and well vascularized bone defects.

DIRECTIONS FOR USE: open the package and hydrate the membrane in sterile physiological solution until you obtain the desired consistency. After shaping the membrane with sterile instruments, place it on the receiving site, stabilize it and suture the gingival flaps.

PACKAGING: sterile disposable blister.

SIDE EFFECTS: none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded.

CONTRAINDICATIONS: do not use in patients who have pathologies or conditions that contraindicate surgery or who have an acute or chronic infection (osteomyelitis). Not indicated for patients with immunodeficiency or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly or near large vessels or nerves. Do not use on poorly vascularized, necrotic or infected receiving sites. Do not use on sites that do not guarantee the stability and protection of the product after grafting. The product must not be used if the packaging is damaged or opened before use.

WARNING: This product has not been tested on pregnant women. The product is for single use only and, as such, should not be re-sterilized. The graft site must be properly washed and cleared, and the product must be well stabilized and protected. Product reserved for surgeons with experience in biomaterial implants: strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on OsteoBio® devices available on the website www.osteobiol.com, as well as the instructions herein. For every single application, the specialist surgeon has the responsibility to ensure proper clinical management of the biomaterials and the patient.

EXPIRY DATE: the product is valid for five years from the date of sterilization and the expiry date is indicated on the label.

STORAGE: The product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label.

DISPOSAL: the product does not contain hazardous substances and it may be disassembled and disposed of as regular municipal waste.

PRODUCT ORIGIN: the last letter on the label of the product code REF indicates the source of the material (S = porcine; E = equine).

Symbols used on the label:

🕒 single use 📅 year and month of expiry 📖 lot number 📖 carefully read the instructions 🌡️🌡️🌡️ Sterilization performed by irradiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed. 🌡️ limitatiom of the temperature

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