EVICROL™ Chemically-Cured Macrofilled Composite

Characteristics:

Bi-component filling composite for the construction of durable aesthetic fillings. Available in the following shades: class III, and IV. Accidental shades - 25, 27 and 45. The composition substance is a combination of high-molecular materials and a transparent anorganic filler. Firm linking of both components by a chemical bond resulting in good physico-chemical properties allowing multipurpose application possibilities.

Indications: Mechanical properties of EVICROL™ nearly reach those of amalgam which results in substantially broad application possibilities in restorative dentistry. EVICROL™ is most often used for making fillings of class III, and IV, according to Black (carries on contact surface of the incisors and cuspids and tooth neck edges). Apart from that, EVICROL™ can also be used for fractured incisal edges (class IV).

Contraindications: Unsuitable indication in the distal part with increased strain.

Side effects: Not described.

Mixing ratio: In volume: 2 drops of the liquid + 2 scoops of the powder In weight: 0.5 ml of the liquid + 1.5 g of the powder.

Application:

- Prepare the material by mixing the powder of liquid component on a mixing paper slab using a plastic friction spatula. Mixing with a metal spatula is not recommended for it leads to metal erosion resulting in the composite colouring. The correct mixing ratio is achieved by measuring the same number of liquid drops as of leveled measuring spoon of powder. In most cases, 2 drops of liquid and 2 leveled measuring spoon of powder are started with. The weight ratio is 1.5 g of powder to 0.5 g of liquid. The prepared mixture has a paste character, the surface of the mixed material is mildly matte and when friction spatula is drawn off, the material stretches easily. To mix a dose of EVICROL™ prepared from 2 drops, 30 - 40 seconds are needed. The prepared mixture remains workable for about 1 minute at room temperature. Solidifying (ideally under a matrix) follows afterwards which takes about 2 minutes. It is necessary to leave the filling resting during this period and after another 3 minutes it is possible to start further working of the filling. The cavity shall be prepared in the usual manner. It is beneficial to bevel the surface enamel somewhat so as to enlarge the contact surface between the composite and the enamel.
- Common underfilling materials such as zinc-phosphate ADHESOR cement or ADHESOR CARBONITE polycarboxylate are used as cavity liners underneath EVICROL™. Zinc-oxide eugenol cement CARYOSAN cannot be used for this purpose. The liner must protect the whole surface against the exposed dentine! Following liner solidification, the enamel margin shall be etched with EVICROL™ etching liquid (enclosed in package) which shall be applied using a cotton wool stick and left to act for 60 seconds. Effective enamel etching is the substantial condition of microretention connection of the EVICROL™ filling to the enamel. In practice, successful etching is apparent as matte change of the enamel surface. In some cases of very resistant enamel, etching must be performed repeatedly. The etching liquid can be removed by thorough rinsing in water (10 - 20 seconds).
- Following drying, fill the cavity with a slight excess of EVICROL™. The best quality filling result can be obtained by solidifying under a cone strip (cellulose acetate) e.g. in class III, cavities, under a polyethylene foil (in class I, cavities with occlusion) or under a matrix (in class V, and IV, cavities). Foil can be removed after 2 minutes of solidification and the filling can be further shaped using a bar following another 3 minutes and/or polished with „DEPURAL NEO paste“ mixed with vaseline.
- Correct working procedure must result in reliable and permanent microretention junction of Evicol with enamel.

Notice: Do not apply in close pulp proximity.
Waste material disposal: Dispose of unused powder and its packaging with ordinary worksite waste. Residues of unused liquid should be transferred to a company authorized to dispose of waste materials or dispose of it in an incineration plant for dangerous wastes; do not pour into sewerage system.
Warning: EVICROL™ H319 Causes serious eye irritation. H315 Causes skin irritation. H335 May cause respiratory irritation. H317 May cause an allergic skin reaction. H411 Toxic to aquatic life with long lasting effects. P202 Do not get eyes, on skin, or on clothing. P203 Do not eat, drink or smoke when using this product. P273 Avoid release to the environment. Powder: H372 Causes damage to organs (lungs) through prolonged or repeated exposure. H319 Causes serious eye irritation. H317 May cause an allergic skin reaction. H400 Very toxic to aquatic life. P262 Do not get in eyes, on skin, or on clothing. P273 Avoid release to the environment. Etching liquid: H314 Causes severe skin burns and eye damage. P280 Wear protective gloves/protective clothing/eye protection/face protection.

Storage: Store in a dry and dark place at 5 °C to 25 °C, in a well sealed container.
Shelf-life: 2 years
Packing: 1 x 40 g powder - basic shade No 21
3 x 10 g powder - accidental shades No 25, 27, and 45
1 x 26 g liquid
1 x 14 g EVICROL™ etching liquid
30 friction spatulas
1 block of friction papers slab
1 powder measuring spoon

SpofaDental s.r.l.
Markova 238, 506 01 Jičín, Czech Republic
www.pentron.eu
+420 493 583 111
Made in Czech Republic

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