K-ETCHANT Syringe

EN	ENGLISH	English	Dental Phosphoric Acid Etchant
FR	FRANÇAIS	French	Agent de mordantage dentaire à base d'acide phosphorique
ES	ESPAÑOL	Spanish	Agente de grabado de ácido fosfórico para uso dental
IT	ITALIANO	Italian	Mordenzante dentale all'acido fosforico
NL	NEDERLANDS	Dutch	DENTAAL FOSFORZUUR
DE	DEUTSCH	German	Dentales Phosphorsäure-Ätzmittel
SV	SVENSKA	Swedish	Dentalt fosforsyra-etsmedel
NO	NORSK	Norwegian	Dentalt fosforsyre-etsemiddel
FI	SUOMI	Finnish	Hampaan fosforihappoetsausaine
DA	DANSK	Danish	Dentalt fosforsyre ætsningsmiddel
PT	PORTUGUÊS	Portuguese	Condicionador dentário à base de ácido fosfórico
EL	ΕΛΛΗΝΙΚΑ	Greek	Οδοντικό Αδρoπιοητικό Φωσφορικό Οξύ
TR	TÜRKÇE	Turkish	Dental Fosforik Asit
PL	POLSKI	Polish	Wytrawiacz stomatologiczny na bazie kwasu ortofosforowego
RO	ROMÂNĂ	Romanian	Gravant dentar pe bază de acid fosforic
HR	HRVATSKI	Croatian	Dentalno sredstvo za jetkanje na bazi fosforne kiseline
HU	MAGYAR	Hungarian	Fogászati foszforsavas maratóanyag
SL	SLOVENŠČINA	Slovenian	Dentalno jedkalo iz fosforjeve kisline
CS	ČESKY	Czech	Dentální leptadlo – kyselina fosforečná
SK	SLOVENSKY	Slovak	Dentálny leptací prostriedok na báze kyseliny fosforečnej
BG	БЪЛГАРСКИ	Bulgarian	Стоматологичен ецващ агент с фосфорна киселина
UK	УКРАЇНСЬКА	Ukrainian	Стоматологічний протравлювач на основі ортофосфорної кислоти
ET	EESTI	Estonian	Dentaalne fosforhappe söövitusaine
LV	LATVIEŠU	Latvian	Zobārstniecības fosforskābes kodinātājs
LT	LIETUVIŠKAI	Lithuanian	Odontologinė fosforo rūgšties esdinimo medžiaga

Kuraray Noritake Dental Inc.
1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague, The Netherlands

Kuraray Europe GmbH
Philipp-Reis-Str. 4,
65795 Hattersheim am Main, Germany

Phone: +49 (0)69 305 35835 Fax: +49 (0)69 305 98 35835
URL: <https://www.kuraraynoritake.eu>

1562R058R

K-ETCHANT Syringe

I. INTRODUCTION

K-ETCHANT Syringe is an etching gel that consists of 35% phosphoric acid aqueous solution and colloidal silica. It allows for precise and selective placement thanks to its flowable but not runny consistency.

II. INDICATIONS

- Etching of enamel and dentin
- Cleaning of dental restorative materials

III. CONTRAINDICATIONS

Patients with a history of hypersensitivity to this product

IV. POSSIBLE SIDE EFFECTS

K-ETCHANT Syringe may cause inflammation or erosion with its irritancy.

V. PRECAUTIONS

1. Safety precautions

- The product may cause inflammation or erosion with its irritancy. Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
- If the product comes in contact with human body tissues, take the following actions:
 - <If the product gets in the eye>
Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin or the oral mucosa>
Immediately rinse with copious amounts of water.
- This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to any components in the product.
- If the patient demonstrates a hypersensitivity reaction, such as rash, eczema, features of inflammation, ulcer, swelling, itching, or numbness, discontinue use of the product and seek medical attention.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of inflammation, erosion or hypersensitivity that may result from contact with the components in the product.
- Use caution to prevent the patient from accidentally swallowing the product.
- The needle tip is single-use only. Do not reuse it to prevent cross-contamination. Discard it after use.
- If the instruments associated with this product are damaged, use caution and protect yourself; immediately discontinue use.
- If the product adheres to clothing, wash it off with water.
- Dispose of this product as a medical waste to prevent infection. The needle tip must be disposed of after covering the tip of the needle to prevent injury.

2. Handling and manipulation precautions

- The product must not be used for any purposes other than specified in [II.INDICATIONS].
- The use of this product is restricted to licensed dental professionals.
- Use a rubber dam to control moisture, if necessary.
- After each use, remove the needle tip from the syringe and recap the syringe immediately and tightly.
- Do not mix the product with other materials.

3. Storage precautions

- The product must be used by the expiration date indicated on the package.
- The product must be stored at 2 - 25°C / 36 - 77°F when not in use.
- The product must be kept away from extreme heat or direct sunlight.
- The product must be stored in a proper place where only dental practitioners can access.

VI. COMPONENTS

Please see the outside of the package for contents and quantity.

- Etching agent
 - Principal ingredients:
 - Phosphoric acid
 - Water
 - Colloidal Silica
 - Pigment
- Accessory
 - Needle tip (E)

VII. CLINICAL PROCEDURE

1. Syringe preparation

- 1-1. Unscrew the cap of the syringe, and attach the needle tip to the syringe by rotating clockwise. Cover the whole syringe with a disposable barrier (e.g. a poly bag) to prevent contamination from saliva or blood. Disinfect the syringe by wiping it with an absorbent gauze pad moistened with alcohol both before covering and after use.
- 1-2. Dispense a small amount of the gel onto a pad or gauze, and make sure the gel is dispensed adequately.

2. Pretreatment of the adherend area

Pretreat the adherend area according to the instructions for use of the adhesive system being used together.

3. Application of the gel

Dispense the gel slowly and apply it to the adherend surface. Refer to the instructions for use of the adhesive system being used together to confirm the appropriate etching time. The standard etching time of this product is 5-10 seconds.

4. Rinsing and Dry

After the predetermined time, rinse thoroughly and dry in the usual manner.

5. Application of dental adhesive

Proceed with the work using the adhesive system according to the specific instructions for use of its manufacturer.

[CAUTION]

- Regarding the etching conditions, follow the instructions for use of the adhesive system being used together.
- Avoid any contamination to the etched tooth surface with saliva, etc. If contamination occurs, repeat the process of etching, rinsing, and drying.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

If a serious accident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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