Original Operating Instructions

S5

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vhf

Contents

1 Welcome	5	6.3 Switching off the machine	27
1.1 About this document	5	6.4 Mounting & removing blanks	27
1.2 Used symbols	5	6.4.1 Mounting blocks to the block holder	28
1.3 Structure of safety notes	5	6.4.2 Using an optional abutment holder	29
2 General safety instructions	6	6.4.3 Operating the blank changer	29
3 Operating regulations	9	6.4.4 Change position of the blank magazine	30
3.0.1 Intended use	9	6.4.5 Mounting discs & block holders to the blank magazine	30
3.0.2 Controlling the machine through software	9	6.5 Managing tools	33
3.0.3 Maintenance and cleaning	9	6.5.1 Inserting and exchanging tools	33
3.0.4 Spindle	9	6.6 Deactivating and activating the ionizer	35
3.0.5 Unattended operation	9	6.7 Mounting and removing the ionizer cover	35
3.0.6 Transportation and storage	10	6.8 Preparing the machine for wet machining	36
4 Machine overview	11	7 Operation: Executing jobs	38
4.1 Front side of the machine	11	7.1 Starting jobs	38
4.2 Connection panel	11	7.2 Aborting machining	38
4.3 Working chamber door	12	7.3 Job interruptions and job abortions	39
4.4 Working chamber	12	7.3.1 How to proceed in case of a job interruption	39
4.5 Blank magazine	13	7.3.2 How to proceed in case of a machine malfunction	39
4.6 CAM computer	13	7.3.3 How to proceed in case of a tool breakage	39
4.7 Anti-soiling concept	13	7.3.4 How to proceed in case of a power failure	40
4.8 Sound emission	14	7.4 Emergency opening of the working chamber door	41
4.9 Location of the identification plate & serial number	14	7.5 Manually resetting the blank changer to default position	42
4.10 Axes	14	8 Maintenance and do-it-yourself	43
4.11 Technical data	15	8.1 Basic maintenance	43
5 Installing the machine	17	8.2 Maintenance section	43
5.1 Checking the scope of delivery	17	8.3 Preventive maintenance	43
5.2 Choosing the installation site	18	8.4 Where to get service?	43
5.3 Machine installation (schema)	19	8.5 Definition of wear parts	43
5.4 Establishing the electric connection	20	8.6 Cleaning the working chamber	44
5.5 Removing the transport lock	20	8.7 Cleaning the blank magazine	45
5.6 Installing the pneumatics	21	8.8 Cleaning the collet chuck	46
5.6.1 Overview compressed air regulator	22	8.9 Checking the compressed air regulator	48
5.6.2 Installing the pneumatic hose	22	8.9.1 Checking the water separator for condensate	48
5.6.3 Adjusting the air pressure with the compressed air regulator	23	8.9.2 Exchanging / cleaning the contaminated filter cartridge	48
5.7 Installing the air extraction system	24	8.10 Cleaning the housing	49
5.7.1 Requirements for the suction unit	24	8.11 Exchanging the main fuse	50
5.7.2 Installing the suction unit	24	8.12 Calibrating the axes	50
5.7.3 Connecting the suction hose with the optional hose connection	24	8.13 Replacing the tool magazine inserts	51
5.7.4 Installing the switching unit	25	8.14 Updating the software and firmware	52
5.8 Connecting the CAM computer	25	8.14.1 Updating the manufacturing software	52
6 Operation: Preparing jobs	26	8.14.2 Updating the firmware of the machine	52
6.1 Starting the machine	26	8.15 Maintenance table	53
6.2 Starting the machine with a tool in the collet chuck	26	9 Disposal	55
		9.1 Disposing of the cooling liquid	55
		9.2 Disposing of the machine	55

10 Troubleshooting	56
--------------------------	----

Index	59
-------------	----

1 WELCOME

Thank you for your purchase of this S5 dental machine. The machine is delivered to you with pride and confidence. It was produced using the latest techniques and strict quality control.

These operating instructions were prepared to help you understand all functions of your new dental machine. It should also help you keep the machine in good conditions so that you may enjoy many hours of productive work.

You may find updates to this document at:
dentalportal.info – search for S5

1.1 About this document

This document is designed and released for the following groups / individuals:

- End users
- Authorized resellers
- Authorized service technicians

1.2 Used symbols

Instructions

» Single or general instruction

1. Numbered action step

✓ Result

Additional symbols

↗ Cross reference

- List (first level)
 - List (second level)

1. *Numbered image labels*

✓ **Correct** or **Do this**

✗ **Incorrect** or **Do not let this happen** or **Don't do this**

Information to make work more efficient

Important information without any danger for people or objects

Additional information

User interface description

[Buttons]

<KEYS>

User interface text

Text that you need to enter

1.3 Structure of safety notes

SIGNAL WORD

Type and source of hazards

Further explanations and consequences when ignoring the hazard.

» Instructions to avoid the hazard.

The following signal words may be used:

DANGER

DANGER indicates a hazardous situation which will result in death or serious injury.

WARNING

WARNING indicates a hazardous situation which can result in death or serious injury.

CAUTION

CAUTION indicates a hazardous situation which can result in minor or moderate injury.

NOTICE

NOTICE indicates a situation which can lead to physical damage of the product or in the surrounding areas.

2 GENERAL SAFETY INSTRUCTIONS

DANGER

Incorrect operation of the machine

-  » **Before** installing, maintaining and operating the machine, read **all** documents provided with the machine.
- » If it is unclear how to operate the machine in any way, do not use the machine and contact customer service.
- » Ensure that every user has access to the operating instructions.
- » Instruct every user on safe and proper machine handling.

Danger to life due to an electric shock

 If you come into contact with electrically charged parts, you can suffer from an electrical shock. Water increases the risk significantly.

- » Do not remove the housing of the machine.
- » Only have qualified electricians work on any electric equipment.
- » Ensure that an operational Residual Current Device / Ground Fault Circuit Interrupter is installed on the electric circuit of the machine.
- » Run power cables so that they cannot be damaged by sharp edges.
- » **Before** switching on the machine, check power cables for damage.
- » **Before** unplugging the power cable, switch off the machine at the main power switch.

-  » In the following cases, disconnect the machine from the electrical source immediately and prevent it from being restarted:
 - When machine connections or electric cables are damaged
 - When liquid is leaking
 - **Before** you check or run electric cables
- » Replace damaged cables with original manufacturer's spare parts.
- » Never perform any troubleshooting while the machine is operating.
- » Only have authorized service technicians repair the machine.
- » Do not touch the ion jets in the working chamber when the machine is switched on.

- » Do not touch the machine and especially the cables with wet or damp hands.
- » Check the environment of the machine and all accessible internal areas daily for leaked liquid and remove any liquids near or in the machine immediately.
- » Never put any machines or devices which are powered by electricity under the machine.
- » Do not put any objects onto the machine.

WARNING

Respiratory diseases when processing harmful materials

If you inhale harmful materials during their processing, your respiratory tracts may be damaged.

- » Always use a suitable air extraction system during dry machining.
- » Use a suction device with an extra-fine particle filter.
- » Avoid materials which damage your health.

Danger to health if using the ionizer when safety devices are inactive

When the external air extraction system or the compressed air supply are disabled or not working as required, ozone concentrations that can threaten your health may occur in the working chamber.

- » Only use the ionizer when the external air extraction system and the compressed air supply are working properly.
- » Do not inhale the ionized air in the direct vicinity of the ionizer.

Danger to health if using the wrong cooling lubricant

Some cooling liquids can be a serious threat to your health and / or the environment.

- » Only add the cooling lubricant Tec Liquid Pro to the cooling liquid. The mixing ratio is indicated on the bottle label.

Crushing hazard and risk of cutting injuries through moving machine parts

Through the moving axes and the rotating spindle you can suffer bruises and cuts.

- » Only use the machine when the working chamber door is completely closed and undamaged during machining.
- » Store the key for the emergency release in a place where only authorized persons have access.
- » Do not circumvent or deactivate safety devices of the machine.

- » Check the machine regularly for damage, especially the safety devices.
- » Have damaged safety devices repaired by customer service only.
- » Use only original manufacturer's equipment and original spare parts in the machine.
- » Keep children and animals away from the machine.
- » Do not remove the housing of the machine.

Service Mode: Risk of cutting injuries and bruises as well as hazards through ejected debris

If you operate the machine in any "Service Mode" with the working chamber door open, the risk of injury is increased significantly.

- » Operate the machine in "User" mode only unless you have been authorized by the machine's manufacturer to use other modes.
- » Even if you are an authorized user, use the "Service Modes" only when necessary.

- » While in any "Service Mode": Do not reach into the working chamber while the axes are moving or during machining.

- » While in any "Service Mode": Everyone within reach of the machine must wear protective eye wear.

Hearing damage due to loud noise

If you are regularly exposed to loud machining noise, you may suffer from hearing loss and tinnitus.

- » If loud noise cannot be avoided, wear ear protection during machining.

Risk of injuries through loose pneumatic components under air pressure when connections are open

Loose pneumatic components can move extremely fast and unpredictably and may cause injury.

- » **Before** you run the pneumatic hoses, close the compressed air supply valve.
- » **Before** you check the pneumatic hoses and pneumatic connections, set the air pressure to a minimum value.

- » In case of defective machine connections and pneumatic hoses, disconnect the machine from the external compressed air supply and the electrical source.

- » Contact customer service if connections are damaged or defective.

Risk of injuries when opening or closing the working chamber door or the blank magazine door

When you open or close the working chamber door or blank magazine door, the moving working chamber door or blank magazine door may crush your fingers. Objects on the machine may fall and cause injuries or damage.

- » When you open and close the working chamber door or blank magazine door, use one hand and keep the other hand away from the machine.
- » When you close the working chamber door or blank magazine door, ensure your hands do not get caught between the door and the machine housing.
- » Do not place objects on the machine.

Trip, fall and slipping hazards

- » Run cables in such a way that persons cannot trip over them.

- » Keep the working environment and installation site clean.

Risk of cutting injuries and burns

If you touch tools or sharp edges on blanks or the machine, you may suffer from cuts. If you touch the hot spindle body or hot tools, you may suffer from burns.

- » Wear gloves when you perform manual work at the machine or with blanks / tools.

Danger to health if handling the cooling lubricant incorrectly

- » **Before** using the cooling lubricant, read the safety data sheet that was supplied with the cooling lubricant.
- » When handling the cooling lubricant, **always** wear suitable protective clothing.
- » **Always** store the cooling lubricant in the original container.

Reduced ability to act with insufficient lighting

In case of an insufficient lighting your judgment and / or your precision may be reduced.

- » Ensure that the lighting in your working environment is sufficient.

Risk of injury in case of malfunctions caused by insufficient maintenance

If you do not maintain the machine as required, malfunctions may occur which can lead to injuries.

- » Take note of the intervals and conditions mentioned in the maintenance table in the operating

instructions. Carry out the respective maintenance tasks accordingly.

Health risks through constant malpositioning if your working environment is not sufficiently ergonomic

Over the long run, an improper or one-sided positioning can be a risk to your health.

- » Set up an ergonomic work environment.
- » Ensure the seat height and monitor position is ideal and the lighting is sufficient.

3 OPERATING REGULATIONS

If you violate the following regulations, you may lose your entitlement to benefits.

NOTICE

Machine damage if you violate these regulations

If you violate the following regulations, your machine may get damaged and / or cause damage in the surrounding area(s).

- » Thoroughly follow all instructions and information in this section.

3.0.1 Intended use

The machine and the manufacturing software have been designed for the commercial processing of approved dental materials by specially trained persons. The processed objects require additional work before their use on patients.

- » Only process materials that you can select in the manufacturing software.
- » Only use the machine and manufacturing software in a commercial environment.
- » Before creating jobs, verify if the objects being prepared may be utilized at the place of use according to local and / or national regulations or other authorized organizations or entities (e.g. professional associations, health authorities). In particular, verify if the material is approved for the machined object type and if the object type is designed in accordance with applicable regulations. Neither the manufacturing software nor the machine will inform you about possible regulatory infringements, but will execute jobs in accordance with the preferences and materials set by the user.
- » Verify that each object type and each material in your jobs are authorized manufacturing materials. If mandated by local or national regulations, obtain relevant authorization from the responsible organizations or entities (e.g. professional associations, health authorities).
- » Only import objects into the manufacturing software which correspond to the object types that you can select in the manufacturing software. While you can import / manufacture any other objects as well, neither the manufacturing software nor the machine are designed for these other objects and should not be used in this way.
- » Do not manufacture implants or parts of objects that are designed to have contact with implants. These parts include parts of two-part abutments which

contain the connection geometry for the implant. Do not manipulate the connection geometry of pre-fabricated abutments (“prefab abutments”) and you must always check finished objects for accurate connection geometries (i.e. that connection geometries of finished jobs have not been damaged).

- » Only operate the machine with original equipment from vhf. If you use non-authorized equipment, the safety of the machine can not be ensured.
- » Only execute jobs that are processed wet, when you have an operational vhf wet grinding module installed.
- » Do not modify the machine or remove the machine housing.

3.0.2 Controlling the machine through software

You control the machine through specially designed applications which are supplied with the machine.

- » Always use the latest program version that officially supports your machine.
- » Before installing or operating the machine, read the documentation for the applications.
- » Ensure that your CAM computer meets all system requirements.

3.0.3 Maintenance and cleaning

Maintenance and cleaning is part of standard machine usage.

- » Clean and maintain the machine as required. Only then can the machine reach a long service life.

3.0.4 Spindle

The spindle of your machine is a high-precision instrument.

- » Do not use unbalanced tools at high rotational speeds. Such an imbalance puts a great strain on the spindle’s ball bearings, which can cause the bearings to be damaged.
- » When working in the working chamber, do not apply manual pressure against the spindle.

3.0.5 Unattended operation

If the machine runs unattended, the risk of material damage is increased.

- » Only allow unattended operation of the machine to occur if the following conditions are met:
 - The national and local laws allow it.
 - The working chamber of the machine is completely clean.
 - Unauthorized users cannot access the machine.

- The room in which the machine is located has an automatic fire detection system.

3.0.6 Transportation and storage

WARNING

Injuries caused by unsafe transportation

If you transport the machine unsafely, the machine may slip and cause injuries.

- » Always transport **unpacked** machines individually and do not stack them.
- » Ensure that only trained personnel transport the machine to and from the installation site.
- » Ensure that the housing of the machine is completely closed.
- » Always transport the machine in an upright position.
- » Transport and position the machine with as many people as required for the weight of the machine in accordance with local and / or national laws and regulations.
- » Grab unpacked machines only at the left and right grips located at the bottom of the machine. **Do not tilt the machine when carrying it.**

NOTICE

Short-circuit hazard when the machine is too cold

If the machine is transported from a cold environment into a warmer environment, a short circuit may occur caused by condensate.

- » **Before** switching on the machine after transportation, ensure the following:
 - The ambient air has the allowed temperature.
 - The machine has the same temperature as the ambient air. This will take **at least** 48 hours.
 - The machine is completely dry.
- » Ensure that the following conditions are met during the whole transport and / or storage period:
 - Permitted ambient conditions for storage / transport:
 - Ambient temperature (storage / transport): between -20 and 60 °C
 - Relative air moisture: max. 80 %, non-condensing
 - Permitted ambient conditions for operation:
 - Indoor location
 - Dust-free location, pollution degree 2 (IEC 60664-1)
 - Height of the installation site: up to 2000 m (6561 ft) above mean sea level

Preparing transportation or storage

Before transporting or storing your machine, the following preparations are necessary:

1. Remove all blanks from the working chamber.
2. Remove all items from the blank magazine.
3. If the machine was used for wet machining:
 - a. Rinse the cooling liquid system.
 - b. Disassemble the wet grinding module as described in the module's operating instructions.
4. Clean the working chamber. Ensure that the working chamber is completely dry.
5. Clean the blank magazine.
6. Install the transport lock. For this, see the corresponding steps on the supplement.
7. Ensure that the housing of the machine is completely closed.
8. Switch off the machine at the main power switch.
9. Disassemble the machine components by following the installation instructions in reverse order.
10. In case of overseas transport, take proper measures against corrosion.

Repackaging

To repack the machine after preparing its transportation or storage, the following steps are necessary:

1. If possible, use the original packaging. If the original packaging is not available, use a packaging of similar size and quality.
 -  **Original packaging is available from customer service.**
2. Pack the machine and its accessories securely.
3. Protect the packaging against slipping. If machines are properly packed and protected against slipping, they may be stacked.

4 MACHINE OVERVIEW

With your S5 you can process blanks of different materials to create high quality objects for the dental sector. You can find a list of the materials which you can process with the machine in the manufacturing software.

The blank changer of the S5 allows you to process up to 8 jobs in succession without having to start each job manually.

The machine is designed for wet *and* dry machining.

During wet machining, the tools and blanks are constantly being cooled by cooling liquid.

During dry machining, our anti-soiling concept reduces the soiling of sensitive machine parts.

Wet machining requires the optional vhf wet grinding module. The wet grinding module and the required cooling lubricant Tec Liquid Pro are available from customer service.

4.1 Front side of the machine


FIG. 1 FRONT SIDE OF THE MACHINE

1. Blank magazine door
2. View window to the blank magazine
3. Connection panel & main power switch on the rear side
4. Suction opening for the air extraction system
5. Identification plate
6. Recessed grip for opening the blank magazine door
7. Recessed grip for opening the working chamber door
8. View window to the working chamber
9. Working chamber door

4.2 Connection panel


FIG. 2 CONNECTION PANEL OF THE S5

1. Power connection including glass fuse T6,3A L250V
2. Main power switch
3. USB port
4. Suction unit data port (6 pins, IEC 61076-2-106)
5. Switching output for the optional external wet grinding module
6. Cooling liquid connection for the optional external wet grinding module
7. Pneumatic connection (6 mm push-in fitting)

4.3 Working chamber door

The working chamber door locks the working chamber and protects the user from injuries during operation. You can open and close the working chamber door manually.

You *cannot* open the door when the machine is switched off or while the axes are moving.

- » To open or close the working chamber door, pull it upwards or push it downwards with your hand. Use the recessed grip of the door.


FIG. 3 WORKING CHAMBER DOOR

4.4 Working chamber

You can mount blanks and insert tools into the working chamber. This is where the blanks are processed.


FIG. 4 WORKING CHAMBER

1. Outlet for the cooling liquid / air extraction opening
2. Blank changer flap
3. Ionizer
4. Spindle
5. Blank holder; Rotational axis A; Rotational axis B
6. Measuring key with protective cage
7. Tool magazine

Colors of the working chamber lighting

- ! If the working chamber lighting is insufficient, provide additional lighting.

The machine illuminates the working chamber in different colors. The color will change depending on the state of the machine. You will find the colors and respective machine status in the following table:

Color	Status
White	The machine is ready for operation. You can open the working chamber door.
Blue	The machine is operating. The working chamber door is locked.
Red	A machine malfunction has occurred. The working chamber door is locked.

4.5 Blank magazine

The blank magazine is the part of the blank changer where you equip the blanks that you want to machine. The blank changer automatically loads and unloads the blanks into and from the working chamber. The machine can be equipped with up to 8 blanks.

The blank magazine door covers the blank magazine and protects the user from injuries during operation. The blank changer can only operate when the blank magazine door is closed. You can open and close the blank magazine door manually.

- » To open or close the blank magazine door, pull it upwards or push it downwards with your hand. Use the recessed grip of the door.


FIG. 5 BLANK MAGAZINE DOOR

4.6 CAM computer

To operate the machine, you must use a computer running Windows® (“CAM computer”) and specially designed software (“manufacturing software”). The manufacturing software consists of the following components:

- **DENTALCAM** – A CAM application for creating and calculating virtual blanks (“jobs”).
- **DENTALCNC** – A CNC application for machining jobs and maintaining the machine.

To create and design the dental objects you also need a CAD application (sold separately by specialist dealers).

4.7 Anti-soiling concept

The anti-soiling concept decreases the soiling and wear of sensitive machine parts.

The anti-soiling concept does not replace the regular cleaning of the machine. Without regular cleaning, the machine life decreases significantly.

NOTICE

Machine damage when dry machining without an air extraction system

When dry machining, an operational external air extraction system must be installed. Otherwise, the machine will get soiled over time and become defective.

- » Always use a correctly installed and operational external air extraction system during dry machining.

The anti-soiling concept consists of:

Dry machining

- **The external air extraction system** – An external suction unit creates a vacuum in the working chamber and extracts machining debris from it.
- **Air that is emitted into the working chamber** – Blows machining debris away from the spindle and other machine parts.
- **The internal vacuum sensor** – Constantly monitors the vacuum in the working chamber.
- **The ionizer** – The ionizer electrically discharges machining debris of certain materials to reduce their distribution in the working chamber.

Wet machining

The cooling liquid keeps machining debris away from sensitive machine parts.

4.8 Sound emission

The actual sound emission of the machine varies heavily depending on the manufacturing material and the machining conditions.

- » If the machine is exceptionally loud, check the following operating conditions:
 - Cleanliness of the blank holder
 - Condition of the tools
 - Quality of the blanks
- » If loud noise cannot be avoided, wear ear protection during machining.

Sound measurement

Measuring conditions:

- Processed material: CoCr
- Tool status: worn
- Measured value: sound pressure level (distance: 1 m)
- Measurement according to ISO 3746, survey method 3

Established sound emission:

Operating condition	A-weighted sound pressure level
Processing	77 dB(A)
All other operating conditions (tool change, movement of the axes etc.)	<70 dB(A)

4.9 Location of the identification plate & serial number

The identification plate of the machine contains identifying information such as the serial number. You can find the identification plate and machine serial number at the following location: [↗ Connection panel – on page 11](#)

4.10 Axes

This machine has 5 axes: 3 linear axes and 2 rotary axes.


FIG. 6 LINEAR AXES X, Y, Z (THE SPINDLE MOVES ALONG THESE AXES.); ROTARY AXES A, B (THE BLANK HOLDER ROTATES ABOUT THESE AXES.)

4.11 Technical data

Dimensions (W/D/H)

	Unit	Value
Footprint (approx.)	mm	610 x 275
	in	24.0 x 10.8
Housing fully closed (approx.)	mm	700 x 445 x 560
	in	27.6 x 17.5 x 22.1
Minimum required space for operation (approx.)	mm	1050 x 770 x 600
	in	27.6 x 26.4 x 22.6

Base system

	Unit	Value
Weight (approx.)	kg	106
	lbs	234
Axes		5
– Rotational axis A / B		360° / ± 30°
Overvoltage category (IEC 60664-1)		II

External compressed air supply

	Unit	Value
Min / max. air pressure	bar	6 / 8
	psi	90 / 120
Recommended air pressure	bar	7
	psi	100
Air consumption (approx.)		
– Wet machining	l/min	30/35 (at 6/8 bar)
	cfm	1.1/1.2 cfm (at 90/120 psi)
– Dry machining	l/min	60/73 (at 6/8 bar)
	cfm	2.1/2.6 (at 90/120 psi)
Air purity (ISO 8573-1:2010)		Solid particles: class 3 Water content: class 4 Residual oil content: class 3

External air extraction system

	Unit	Value
Minimum extraction capacity	l/min	3,000 (at 220 hPa)
	cfm	105.9 (at 3.2 psi)
Suction unit		Designed for the commercial use in the dental sector Equipped with a filter of the filter class M

Ambient conditions

	Unit	Value
Relative air moisture		80 %, non-condensing
Ambient temperature for storage/transport	°C	-20 – 60
	°F	-4 – 140
Ambient temperature for operation	°C	10 – 35
	°F	50 – 95
Location type		Indoor
Maximum height above mean sea level	m	2000
	ft	6561
Ambient air (IEC 60664-1)		Dust free, pollution degree 2

Spindle

	Unit	Value
Model		SFS 300P (synchronous)
Maximum rotational speed	rpm	60,000
Nominal power under continuous operation (S1)	W	300
Nominal power under uninterruptible periodic operation	W	450

	Unit	Value
(S6)		
Peak power (P_{max})	W	600
Diameter collet chuck	mm	3

Tool changer

	Unit	Value
Maximum tools in magazine		16
Maximum tool length	mm	40

Standard blank holder

	Unit	Value
Min. / max. disc diameter	mm	98.3 / 98.7
Maximum disc height	mm	30
Min. / max. height of the disc rim	mm	10 / 10.5

Cooling liquid system

	Unit	Value
Prepared for wet machining		Requires vhf wet grinding module

Blank changer

	Unit	Value
Number of slots		8
Mounting method		Blank frames required (8 provided)

Connections

	Unit	Value
Pneumatic connection, push-in fitting (Diameter)	mm	6
Power connection	V AC	100 – 240
	Hz	50/60
	W	850, Glass fuse T6,3A L250V
USB port		2.0 A
Data port for supported suction unit or switching unit		Yes
Hose connection for external air extraction system (Diameter)	mm	45
Data port for wet grinding module		Yes
Hose connection for wet grinding module (Diameter)	mm	6

Pin assignment of the suction unit data port

FIG. 7 PIN ASSIGNMENT OF THE SUCTION UNIT DATA PORT

1. + 24 V, switched, output 1
2. + 24 V, switched, output 2 (not used)
3. +24 V, continuous
4. Return signal output 2 (not used)
5. Return signal output 1 (not used)
6. Ground (PE)


5 INSTALLING THE MACHINE

5.1 Checking the scope of delivery

» Unpack the machine and ensure that you have received the following items:

1

2

3

4

5

6

7

8

9

10

11

12

13

14

15

16

1. 1 x Machine S5
2. 1 x Spindle service set
3. 1 x Power cable
4. 1 x USB cable
5. 1 x Pneumatic hose
6. 1 x Compressed air regulator
7. 1 x Torx screwdriver (TX 10)
8. 1 x Hexagon screwdriver (2.5 mm)
9. 1 x QuickFrame for tool-less blank fixation
10. 1 x Key for the emergency release of the working chamber door
11. 1 x Calibration set: 1 micrometer, 3 blanks for manufacturing test and calibration specimen, 1 radius cutter with 2 teeth (P200-R2-40)
12. 2 x Tool magazine inserts
13. 1 x Drill bit (2.8 mm) for tool positions
14. 1 x Measuring pin
15. 1 x Hose connection for the external air extraction unit
16. 1 x Crevice nozzle (for cleaning the working chamber)

Not depicted:

- This document
 - 1 x Transport lock in the working chamber
 - 1 x Supplement about removing the transport lock
 - 1 x Measuring pin for customer service
 - 1 x Set of replacement screws for the blank holder and the tool magazine cover
- » Keep the packaging of the machine, the carrying aid and the transport lock for future transports.

5.2 Choosing the installation site

The installation site must meet the following criteria:

- Firm and even surface, must be able to carry the weight of the machine.
- Alternating current source.
- An operational Residual Current Device / Ground Fault Circuit Interrupter on the electric circuit of the machine.
- Machine requires an external air extraction system.
- Machine requires an external compressed air supply.
- Access to the internet.

You can find specific values and additional requirements in the chapter on technical data. [Technical data – on page 15](#)

Distances to maintain

NOTICE

Damaging of the machine if safety distances are not maintained

If you do not maintain the safety distances, the movable parts of the housing can collide with obstacles when being opened and get damaged. If the ventilation openings are covered, the machine may overheat and get severely damaged.

- » Ensure that the following safety distances are always maintained.


FIG. 8 DISTANCES TO MAINTAIN – TOP VIEW


FIG. 9 DISTANCES TO MAINTAIN – SIDE VIEW

5.3 Machine installation (schema)

! You can *either* use the switching unit including the control cable *or* the data cable of supported suction units. The data cable must be provided by the manufacturer of the suction unit.

i How to install the vhf wet grinding module is described in the operating instructions for the module.

The operating instructions are delivered with the vhf wet grinding module. They are also available for download at:
<https://www.dentalportal.info/NO1>


FIG. 10 MACHINE INSTALLATION (SCHEMA)

1. Power connection
2. USB cable
3. Pneumatic hose
4. Suction unit
5. External compressed air supply
6. Data cable of supported suction units (optional)
7. Control cable of the switching unit (optional)
8. Suction hose
9. Switching unit (optional)

5.4 Establishing the electric connection

NOTICE

Damaging of the machine through heavy voltage fluctuations and power surges

Heavy voltage fluctuations and power surges can disrupt the control unit and can cause system failures.

- » Connect the machine to a dedicated current circuit or ensure that no devices are connected that can cause heavy voltage fluctuation when switched on.
- » If heavy voltage fluctuations cannot be avoided, install a surge protector that protects the machine from heavy voltage fluctuations.

NOTICE

Short-circuit hazard when the machine is too cold

If the machine is transported from a cold environment into a warmer environment, a short circuit may occur caused by condensate.

- » **Before** switching on the machine after transportation, ensure the following:
 - The ambient air has the allowed temperature.
 - The machine has the same temperature as the ambient air. This will take **at least** 48 hours.
 - The machine is completely dry.

The machine requires a continuous power supply for proper operation.

1. Plug the provided power cable into the power connection at the connection panel of the machine.
2. If power failures occur regularly at the installation location or if there are frequent voltage fluctuations, install an Uninterruptible Power Supply (UPS), type online / VFI (IEC 62040-3, class 1).

 If a power failure occurs during job execution, the tool may break and the blank may be destroyed.

3. Insert the plug of the cable into a socket that is protected by a Residual Current Device / Ground Fault Circuit Interrupter.

5.5 Removing the transport lock

Before operating the machine for the first time, you must remove the transport lock. The transport lock prevents the spindle from getting damaged during transport.

1. Ensure the following:
 - The machine is connected to the electrical source.
 - The CAM computer *is not* connected to the machine.
2. Switch on the machine at the main power switch.
3. Open the working chamber door.
4. Switch off the machine at the main power switch.
5. Remove the transport lock as shown in the supplement.

5.6 Installing the pneumatics

WARNING

Risk of injuries through leaking compressed air and lashing pneumatic hoses

Open or loose pneumatic connections can cause severe injuries.

- » Ensure that **during installation and maintenance** of the pneumatic hoses and of the compressed air regulator compressed air is not conducted through the hoses and connections.
- » **Before** conducting compressed air through the hoses and connectors, verify that the hoses are securely inserted into the correct connectors and that they are not damaged. This also applies to the compressed air regulator.
- » Do not conduct compressed air through damaged hoses and connectors.

NOTICE

The spindle may suffer bearing damage and electrical damage if the compressed air is contaminated

The incoming compressed air must be dry and oil-free according to ISO 8573-1:2010 because the compressed air regulator only serves as an **indicator** for contaminated air.

Air purity according to ISO 8573-1:2010

Solid particles	class 3	Filtration degree better than 5 µm for solid particles
Water content	class 4	Maximum pressure dew point +3 °C
Residual oil content	class 3	Maximum oil content: 1 mg/m ³

- » Ensure that the compressed air meets the above requirements.
- » Only connect the machine to the compressed air supply if the compressed air regulator is properly installed.
- » Connect the machine to the compressed air supply only via the provided compressed air regulator.

You can find specific values and additional requirements in the chapter on technical data. [↗ Technical data – on page 15](#)

The machine requires the compressed air for the following tasks:

- For the opening and closing of the collet chuck during tool changes.

- For the spindle sealing air which prevents foreign bodies from entering the spindle.
- For the sealing air in the working chamber that keeps machining debris away from sensitive machine parts.
- For the ionizer.

5.6.1 Overview compressed air regulator

The machine is connected to the external air supply via a compressed air regulator. You can use this regulator to monitor and regulate the pressure of the incoming air.

The compressed air regulator is supplied with the machine and needs to be mounted on the side of the machine housing when installing the machine. The regulator has the following connections:

- 1/8" internal thread, fitted with male compressed air connector to connect the external compressed air supply
- 6 mm push-in fitting to connect the machine.

NOTICE

Failure of the water separator caused by a wrong alignment of the compressed air regulator

The compressed air regulator must **always** be mounted **in an upright position** because otherwise the water separator will not work.

- » Mount the compressed air regulator in an upright position.

On the left side of the machine are two bores which you can use to mount the compressed air regulator on the machine.

- » Mount the compressed air regulator in an upright position, using the oval-head screws provided in the bores.


FIG. 11 MOUNTING THE COMPRESSED AIR REGULATOR

5.6.2 Installing the pneumatic hose


FIG. 12 INSTALLING THE PNEUMATIC HOSE

- ! The blue compressed air socket must not be confused with the stainless steel socket for cooling liquid. Otherwise the machine will not work properly.

1. Close the external compressed air supply valve.
2. Use the provided pneumatic hose to connect the *right* pneumatic connection of the compressed air regulator to the pneumatic connection of the machine.
3. Connect the external compressed air supply to the *left* pneumatic connection of the compressed air regulator.
4. Thoroughly verify that all external pneumatic hoses are properly seated in their corresponding connections and that the hoses and connectors are undamaged.
5. If all hoses and connectors are properly installed and undamaged, open the external compressed air supply valve.

5.6.3 Adjusting the air pressure with the compressed air regulator

Setting the air pressure is only necessary if the air pressure shown by the pressure gauge does not lie between the minimum and maximum air pressure. You can find specific values and additional requirements in the chapter on technical data.


FIG. 13 COMPRESSED AIR REGULATOR: REGULATING AND MONITORING THE AIR PRESSURE

1. Rotary knob for pressure regulation
2. Pressure gauge for monitoring the outgoing air pressure
3. Bowl of the water separator
4. Discharge screw

1. Pull the rotary knob on top of the compressed air regulator slightly upwards.
 2. Turn the rotary knob in the desired direction:
 - Turn it towards "+" to increase the pressure
 - Turn it towards "-" to decrease the pressure
 3. Push the rotary knob down again.
- ✓ The knob is locked and cannot be changed inadvertently.


FIG. 14 SETTING THE AIR PRESSURE

5.7 Installing the air extraction system

Components of the air extraction system:

Component	Source	Required?	Prerequisite
Suction unit incl. suction hose	Customer service, specialist dealers	Yes	–
Switching unit	Customer service*	No	Data cable <i>not used</i>
Data cable of supported suction units	Manufacturer of the suction unit	No	Supported suction unit, switching unit <i>not used</i>
Hose connection	Customer service	If the suction hose does not fit	–

*The switching unit is not available worldwide.

5.7.1 Requirements for the suction unit

You can find specific values and additional requirements in the chapter on technical data. [↗ Technical data – on page 15](#)

» Use a suction device with the following properties only:

- Designed for the commercial use in the dental sector
- Equipped with a filter of the filter class M
- Equipped with safety devices which protect you from static discharges (e. g. through an anti-static suction hose)

5.7.2 Installing the suction unit

You can install the suction unit as follows:

1. Read the documentation for the suction unit. Follow the operating and safety instructions at any time.
2. Check if the connection of the suction hose has an outer diameter of 45 mm. If the diameter is different, either adjust the hose or use the optional hose connection.
3. Insert the suction hose of the suction unit into the opening for the air extraction of the machine. Ensure that the suction hose is properly seated.


FIG. 15 INSERTING THE SUCTION HOSE INTO THE MACHINE

4. If you want the machine to automatically switch the suction unit on and off, choose 1 of the following options:
 - Install the switching unit (extra equipment). The switching unit is not available worldwide.
 - Connect a data cable provided by the manufacturer of the suction unit to the suction unit data port of the machine. The operating manual of the suction unit should contain additional details.
5. Continue with the installation of the suction unit as described in the documentation of the unit.

5.7.3 Connecting the suction hose with the optional hose connection

If you cannot connect the suction hose of the suction unit directly to the machine, install the hose connection as follows:

1. Turn the thread of the hose connection counter-clockwise until the connection is completely open. If the thread gets detached from the hose connection, place it onto the connection again and turn it clockwise once so that it is screwed to the connection again.
2. Insert the suction hose of the suction device completely into the hose connection on the side of the thread.


FIG. 16 INSERTING THE SUCTION HOSE INTO THE HOSE CONNECTION

3. Turn the thread of the hose connection clockwise as far as you can.
- ✓ The suction hose is firmly attached to the hose connection.
4. Insert the hose connection into the opening for the air extraction system of the machine. Ensure it is firmly connected.
- ✓ The installation of the suction hose with the optional hose connection is complete.


FIG. 17 INSERTING THE HOSE CONNECTION INTO THE OPENING FOR THE AIR EXTRACTION SYSTEM

5.7.4 Installing the switching unit

If you want the machine to automatically switch the suction unit on and off, but a data cable is not available, you can use the optional switching unit. The switching unit is not available worldwide.

1. Connect the power cable of the suction device to the switching unit.
2. Connect the control cable of the switching unit to the suction device data port at the connection panel of the machine.
3. Plug the switching unit into a power socket.


FIG. 18 CONNECTING THE SWITCHING UNIT TO THE SUCTION DEVICE AND THE MACHINE

5.8 Connecting the CAM computer

! The transport lock *must not* be installed when you connect the CAM computer to the machine. [Removing the transport lock](#) – on page 20

» If you want to control multiple machines with 1 CAM computer, use our multi-machine control feature. See the documentation for the manufacturing software.

1. Switch on the machine at the main power switch.
2. Close the working chamber door.
3. Start the CAM computer.
4. Use the provided USB cable to connect a USB port of the CAM computer to the USB port at the connection panel of the machine.

i To avoid connection errors, we highly recommend that you use the USB or Ethernet cable which was provided with the machine. Original spare parts are available from customer service.

5. Install the newest version of DENTALCAM & DENTALCNC that is released for the machine. If you want to control multiple machines with 1 CAM computer, use our multi-machine feature. See the documentation for the manufacturing software.
6. To determine the USB port number and connect to the machine, select the depicted icon in the **Application settings** view in DENTALCNC. 

✓ DENTALCNC will try to establish a connection to the machine. If this is successful, the application will display the port number beside the icon. The machine references.

! The machine will not reference if the working chamber door is open.

7. If the working chamber door was open during the previous step, close it. If the machine does not reference as a result, quit DENTALCNC and restart the application.

6 OPERATION: PREPARING JOBS

Before you can machine blanks, you need to prepare the machine. The corresponding jobs must have been transferred to DENTALCNC where they display in the job list.

6.1 Starting the machine

NOTICE

Short-circuit hazard when the machine is too cold

If the machine is transported from a cold environment into a warmer environment, a short circuit may occur caused by condensate.

- » **Before** switching on the machine after transportation, ensure the following:
 - The ambient air has the allowed temperature.
 - The machine has the same temperature as the ambient air. This will take **at least** 48 hours.
 - The machine is completely dry.

You usually start the machine as follows:

1. Ensure that the machine is properly installed.
 2. Switch on the machine at the main power switch.
 3. If the machine controls the suction unit, switch on the unit and select the necessary extraction level.
- ✓ The suction unit is not running. If it is, either the switching unit or the data cable of supported suction units is not properly installed.

i If you control the suction unit manually, you will switch it on immediately before job execution.

4. Close the working chamber door.

! The machine will not reference if the working chamber door is open.

5. Start the CAM computer.
 6. Start DENTALCNC.
- ✓ The following happens:
- a. The machine references.
 - b. The working chamber is illuminated in white.
7. If the machine doesn't reference because the working chamber door is open, close the door. Wait until the machine has referenced.
- ✓ After the machine has referenced, it is operational.

6.2 Starting the machine with a tool in the collet chuck

Under some circumstances like a power outage, there can be a tool in the collet chuck of the spindle when you start the machine. You need to remove the tool from the collet chuck before you can use the machine.

CAUTION

Danger of cuts and burns when touching tools with

your bare hands

If you handle tools on the cutting surface, you may be injured. As the tool may be very hot, you may also suffer from skin burns.

- » Only touch tools at their shank.
- » When handling tools, wear protective gloves.

NOTICE

Machine damage if you do not remove the tool

If the tool remains in the spindle after you have confirmed the message, it will collide with machine parts such as the measuring key and severely damage them.

- » **Always** follow the instructions below when you start the machine with a tool in the collet chuck.

1. Start the machine.
- ✓ DENTALCNC displays that there is a tool in the collet chuck.
2. Open the working chamber door.
3. **CAUTION!** Wear gloves.
4. Hold the tool in the collet chuck in place.
5. Confirm the current message.
- ✓ The following happens:
 - a. The collet chuck opens.
 - b. The current dialog window closes.
 - c. A dialog window opens.
6. Remove the tool from the collet chuck.


FIG. 19 REMOVING THE TOOL FROM THE COLLET CHUCK

7. Confirm the current message.
- ✓ The machine is ready for operation.

6.3 Switching off the machine

**DANGER**

Danger of an electric shock if unplugging the power cable before switching off the machine

If you unplug the power cable while the main power switch is still in “ON” position, residual voltage in the power cable may cause you to suffer from an electric shock.

- » **Before** unplugging the power cable, switch off the machine at the main power switch.

To switch off the machine, do the following:

1. Ensure that the working chamber is clean.
2. Switch off the machine at the main power switch.
3. (Optional) Unplug the power cable.
4. (Optional) Switch off the main power switch of your workstation / facility.

6.4 Mounting & removing blanks

The machine can process the following blanks:

- Discs with a diameter of 98.5 mm – 98.8 mm
- Blocks, max. size: 45 x 20 x 20 mm (L/D/H)*
- Prefabricated abutments*

*requires extra equipment

You can obtain extra equipment from customer service.

Tools used for cobalt-chrome blanks are shorter than tools for other materials. Only cobalt-chrome blanks with a **maximum thickness of 18 mm** can be machined.

How to mount the different blank types

Blank type	Holder required?	How to mount
Discs	No	Discs >> Blank frame* >> Blank magazine
Blocks	Yes	Blocks >> Block holder >> Blank magazine
Prefabricated abutments	Yes	Prefabricated abutments >> Prefabricated abutment holder >> Blank

Blank type	Holder required?	How to mount
		frame* >> Blank magazine

* 8 blank frames provided with the machine

6.4.1 Mounting blocks to the block holder

You can mount up to 3 blocks to the block holder.

- !** The block type, the number of blocks and the block positions must correspond to the corresponding job in DENTALCAM.


FIG. 20 THE S5 BLOCK HOLDER

1. Upper positioning groove
2. Clamping bolt
3. Lower positioning groove
4. Positioning pin (x 3)
5. Hole for block shafts (x 3)
6. Fixing screw (x 3)

1. Use the provided torque wrench to loosen the fixing screw at the desired position.

FIG. 21 UNSCREWING THE SCREW AT THE DESIRED POSITION

2. If the blank is a multilayer blank, orientate the blank so that the top color layer is on top.

3. Position the block so that the positioning pin at the block holder lies in the groove in the block shaft.


FIG. 22 POSITIONING PINS (LEFT MARKINGS) MUST LIE IN THE GROOVE IN THE SHAFT (RIGHT MARKING)

4. Insert the shaft of the block into the desired position in the block holder until it is firmly seated.


FIG. 23 INSERTING A BLOCK INTO THE BLOCK HOLDER

FIG. 24 INSERTING A BLOCK INTO THE BLOCK HOLDER

5. Use the provided torque wrench to fixate the block with the corresponding fixing screw. Tighten the screw firmly.


FIG. 25 SCREWING DOWN THE SCREW TO IMMOBILIZE THE BLOCK

- ✓ You can now mount the block holder.

6.4.2 Using an optional abutment holder

With an optional abutment holder, your machine can process prefabricated abutments with prefabricated connection geometries.

- i** Abutment holders for common abutment systems can be obtained from customer service.

You can find more information about abutment holders and abutment blanks at the following internet address:

dentalportal.info/abutments

6.4.3 Operating the blank changer

You can mount up to 8 blanks to the magazine slots of the blank changer. The machine automatically loads blanks into the working chamber when they are required for the next job. This allows you to process up to 8 jobs in succession without having to start each job manually.

CAUTION Crushing hazard caused by the moving blank magazine and blank gripper

If you reach into the blank magazine while the blank magazine or blank gripper is moving, you may suffer from bruises.

- » **Before** moving the blank magazine via DENTALCNC, remove your hands from the blank magazine and close the blank magazine door.
- » While the machine is operating, always expect the blank magazine to suddenly start moving.
- » When you are working in the blank magazine while the machine is operating, always keep your hands at the change position and **never** reach into the blank magazine.
- » If you manually move the blank magazine, ensure that you do not pinch your fingers.

To the left of each blank magazine slot is a laser engraved number labeling the corresponding slot.


FIG. 26 THE BLANK MAGAZINE (SLOT NUMBERS MARKED ORANGE, ARROWS INDICATE THE SLOTS THEY CORRESPOND TO)

- The blank types in the physical slots *must* correspond to the blank types which are assigned to the virtual slots in DENTALCNC.
- To move the blank magazine or to manually load blanks into the working chamber, use the **Machining** view in DENTALCNC.

- i** How to manually operate the blank changer and how to assign jobs to blank magazine slots is described in the documentation for DENTALCNC.

6.4.4 Change position of the blank magazine

When you mount blanks to a blank magazine slot, you *must* use the change position. Otherwise, the blanks in your machine will not correspond to the blank positions assigned in DENTALCNC.

- !** Never mount blanks to a different position, even if the position is easily accessible. The machine may load the wrong blanks, which can result in unusable machining results and broken tools.

The change position is the position where the blank frame is oriented vertically in the blank magazine. Every blank magazine slot can be moved into the change position. The change position is shown in the following figure:


FIG. 27 THE CHANGE POSITION (COLORED ORANGE)

6.4.5 Mounting discs & block holders to the blank magazine

You can mount blanks to the blank magazine while the machine is operating.

- !** To relocate blanks to different blank magazine positions, *do not simply swap them around*. The blanks in the blank magazine must correspond to the blanks in the virtual blank magazine in DENTALCNC.

1. Move the desired blank magazine slot into the change position via DENTALCNC. Make sure the blank magazine slot is available in DENTALCNC.
2. Open the blank magazine door.

3. Remove the blank frame or block holder from the change position:
 - a. With one hand keep pressed the blank frame release (marked orange).
 - b. With the other hand pull the blank frame or holder (marked blue) to the right and remove it.


FIG. 28 BLANK FRAME RELEASE (MARKED ORANGE) AND BLANK FRAME IN CHANGE POSITION (MARKED BLUE)

4. To mount a disc to a blank frame, do the following:
 - a. Loosen the 4 fixing screws (marked orange) of the blank frame and open the fixing brackets (marked blue). Do not remove the screws.


FIG. 29 OPENING THE FIXING BRACKETS OF THE BLANK FRAME

- b. Remove any object from the blank frame.
- c. Insert the desired disc into the blank frame.
- d. If the disc has been marked for remounting, rotate the disc until the marking is aligned with the corresponding recesses (marked orange) in the blank frame ([↗ Fig. 30 on the facing page](#)).


FIG. 30 RECESSES (MARKED ORANGE) TO BE ALIGNED WITH MARKINGS FOR SAFER REMOUNTING

- e. Close the fixing brackets of the blank frame (marked blue) and screw in the fixing screws (marked orange).


FIG. 31 CLOSING THE FIXING BRACKETS OF THE BLANK FRAME

- ✓ The disc is mounted to the blank frame.
5. (Optional) To mount a disc to a *QuickFrame*, do the following:

i The *QuickFrame* is a blank frame that uses magnets to fixate the blank. It is suitable for PMMA and zirconia discs with a disc rim height of 10 to 10.3 mm. The *QuickFrame* is available from customer service.

- a. Place the blank frame on an even surface.
- b. With one hand hold the blank frame and the mounted disc (if any).

! The blank frame cover is held in place by magnets. When you remove the cover, the blank is no longer secured and may fall out when you tilt the blank frame.

- c. To remove the blank frame cover, tilt it side-wards with your other hand. Use the recess

(marked blue) on the left and right side of the cover.


FIG. 32 REMOVING THE BLANK FRAME COVER

- d. Remove any object from the blank frame.

NOTICE

Damages to the machine or bad machining results due to improper blank fixation

If you do not properly clean the blank frame and especially the magnets or if you do not correctly place the blank frame cover, the blank may not be properly fixated. As a result, the machine and tools may get damaged and the machining results may not be satisfactory.

- » **Before** placing a disc, thoroughly clean the blank frame, especially the magnets, positioning pins and recesses.
- » **Always** place the blank frame cover so that it fits the positioning pins and is flush with the bottom part of the blank frame.

- e. Clean the bottom part and the cover of the blank frame with a cloth. Especially clean the magnets (marked green), positioning pins (marked orange) and recesses (marked blue). Use a mild cleaning agent if necessary.

FIG. 33 CLEANING THE BLANK FRAME; MAGNETS MARKED GREEN, POSITIONING PINS MARKED ORANGE, RECESSES MARKED BLUE

! The rim height of the disc that you are mounting has to be in the range of 10 to 10.3 mm. Otherwise the blank change process might fail and machining may be aborted.

- f. Insert the desired disc into the blank frame.
- g. If the disc has been marked for remounting, rotate the disc until the marking is aligned with the corresponding recesses (marked orange) in the blank frame ([↗ Fig. 34 see below](#)).


FIG. 34 RECESSES (MARKED ORANGE) TO BE ALIGNED WITH MARKINGS FOR SAFER REMOUNTING

- h. To place the blank frame cover, position the recesses (marked blue) of the cover on the positioning pins (marked orange) of the blank frame. Make sure that the cover is flush with the bottom part of the blank frame.


FIG. 35 PLACING THE BLANK FRAME COVER (POSITIONING PINS MARKED ORANGE, RECESSES MARKED BLUE)

✓ The cover is held in place by magnets. The disc is mounted to the blank frame.

6. To mount a blank frame or block holder to the blank magazine, do the following:

- a. With one hand keep pressed the blank frame release, with the other hand pick up the blank frame or block holder.
- b. Orientate it so that the fixing screws face you and the clamping bolt (marked orange) is in the lower left corner.
- c. Push the blank frame or block holder into the blank magazine slot at the change position and hold it in place. The clamping bolt must lie completely in the recess.


FIG. 36 INSERTING A BLANK FRAME INTO THE BLANK MAGAZINE

- d. Remove the other hand from the blank frame release.
7. Check if the blank frame or block holder is immobilized in the blank magazine slot. Otherwise remove it and repeat the previous step.
- ✓ The blank or designated holder can be loaded into the working chamber.

6.5 Managing tools

NOTICE

Damaging of the spindle or the tool positions if you use improper tools

Improper tools can damage the collet chuck of the spindle and / or the tool positions.

- » Only use tools with a sufficiently large chamfer at the tool shank.
- » Install a retaining ring as a stop ring according to DIN 471-A3.
- » Only insert tools with a maximum diameter of 3 mm at the thickest part into the collet chuck.
- » Only insert tools with a maximum cutting edge diameter of 2.6 mm into the tool magazine.

We recommend original tools as they are designed especially for the designated jobs.

- » If you use tools from other manufacturers, ensure that the following specifications are met:
 - Maximum tool length: 40 mm
 - Maximum shank diameter: 3 mm
 - Minimum shank length: 7 mm (recommended clamping length: 13-14 mm)

You can insert up to 16 tools into the tool magazine.


FIG. 37 TOOL MAGAZINE IN THE WORKING CHAMBER (MARKED ORANGE)

Normally, the machine loads tools from the tool magazine of the machine. In the following cases, this is not possible:

- A tool magazine position is worn so that it can no longer properly hold a tool
- You use tools which do not fit into the tool magazine

In these cases, you can force a manual tool change in DENTALCNC (see the corresponding documentation).

- i** If tool magazine inserts are worn, you should replace them. [↗ Replacing the tool magazine inserts – on page 51](#)

6.5.1 Inserting and exchanging tools

You need to insert or exchange tools in the following cases:

- Upon first use of the machine
- After exchanging tool magazine inserts due to wear
- When the tool life of a tool has expired
- When tools are damaged or worn
- If the next jobs require additional / different tools than those in the tool magazine

- i** You can check the remaining tool life of all tools in the **Tool management** view in DENTALCNC.

You can equip the tool magazine of your machine in two ways:

- By manually inserting the tools into the tool magazine. This is the fastest way. It is described below.
- Via the spindle – you insert a tool into the collet chuck and the spindle will deposit the tool in the tool magazine. This option takes more time, but might be more convenient for some. It is described in the documentation for the manufacturing software.

You insert tools manually as follows:

1. Remove worn or damaged tools from the tool magazine if any.
2. Insert the new tools:
 - a. Ensure that the positions of the tools in the tool magazine match the tool positions in DENTALCNC.
 - b. Insert the tools straight into the tool positions with the cutting edge pointing downwards. Push them in until the ring touches the rubber.

! If the positions of the tools in the tool magazine do not match the tool positions in DENTALCNC, the machine will use the wrong tool(s) during job execution and the job result will become unusable.

3. After replacing a tool, reset the tool life value. For this, use the Tool management view in DENTALCNC.



FIG. 38 INSERTING TOOLS INTO THE TOOL POSITIONS

1	U030-R2-40	9	P250-F1-40
2	U050-F2-40	10	Z100-R2-40
3	U060-R2-40	11	Z100-R2D-40
4	U120-F2-40	12	Z120-F2D-40
5	P100-R1-40	13	Z200-R3-40
6	P100-R2-40	14	Z200-R3D-40
7	P200-R1-40	15	Z060-R2D-40
8	P200-R2-40	16	G260-T-35

FIG. 39 TOP: TOOL POSITIONS 1 – 16 IN THE TOOL MAGAZINE
 BOTTOM: TOOL POSITIONS 1 – 16 IN DENTALCNC

6.6 Deactivating and activating the ionizer

The ionizer considerably reduces the required cleaning effort after dry machining. However, using the ionizer also contributes to the overall compressed air consumption of the machine.

You can decide whether you want to activate or deactivate the ionizer via DENTALCNC. This is a general setting that will affect all jobs that you machine.

- » To deactivate or activate the ionizer, do the following:
 - a. Start DENTALCNC.
 - a. Open the DENTALCNC **Application settings** with the following icon in the main icon bar: 
 - a. Open the **General settings** with the following icon in the local icon bar: 
 - b. Deactivate or activate the **Ionizer activated** option.

6.7 Mounting and removing the ionizer cover

The ionizer can only operate during dry machining. A removable cover provides protection from moisture during wet machining. You need to remove the ionizer cover for dry machining and mount it for wet machining.

NOTICE

Machine damage if you do not mount or remove the ionizer cover as required

If the ionizer cover is not mounted for wet machining, the cooling liquid can cause a short circuit and damage the machine. If the ionizer cover is not removed for dry machining with ionizer, the pressurized air from the ion jets can create an overpressure and damage the machine.

- » **Before** you start wet machining, mount the ionizer cover.
- » **Before** you start dry machining with ionizer, remove the ionizer cover.
- » Before you start dry machining with ionizer, remove the ionizer cover:
 - a. Open the working chamber door.
 - b. Switch off the machine at the main power switch.
 - c. Pull down the ionizer cover (marked orange) at the rear side, and remove it from the ionizer.


FIG. 40 REMOVING THE IONIZER COVER

- d. Place the cover in the storage area above the ionizer: position the recesses (marked blue) of the cover on the pins (marked orange) of the storage area.


FIG. 41 PLACING THE IONIZER COVER IN THE STORAGE AREA (RECESSES MARKED BLUE, PINS MARKED ORANGE)

- ✓ The ionizer cover is held in place by magnets.
- » Before you start wet machining, mount the ionizer cover:
 - a. Open the working chamber door.
 - b. Switch off the machine at the main power switch.
 - c. Remove the ionizer cover from the storage area, by pulling the upper edge towards yourself.
 - d. Place the cover on the ionizer: position the recesses (marked blue) of the cover on the pins (marked orange) next to the ionizer.


FIG. 42 INSTALLING THE IONIZER COVER (RECESSES MARKED BLUE, PINS MARKED ORANGE)

- ✓ The ionizer cover is held in place by magnets.

6.8 Preparing the machine for wet machining

For wet machining, your S5 requires the vhf wet grinding module. The wet grinding module and the required cooling lubricant Tec Liquid Pro are available from customer service.

- ! » **Before you operate the machine with the vhf wet grinding module, read all documents provided with the module.**

The operating instructions are delivered with the vhf wet grinding module. They are also available for download at: <https://www.dentalportal.info/N01>

You prepare the machine for wet machining as follows:

1. Install the wet grinding module and connect it to the machine.
2. Clean the wet grinding module and fill it with sufficient cooling liquid.
3. Switch on the machine at the main power switch.
4. Open the working chamber door.
5. Mount the ionizer cover.
6. Clean the working chamber.
7. Close the working chamber door.
8. Start DENTALCNC on the CAM computer.
9. If DENTALCNC asks you to unlock the wet grinding module, enter the unlock code that you received with the module. If you did not receive an unlock code, contact customer service.
10. If a tool is already inserted in the collet chuck, go to step 18 ([page 36](#)).
11. Select the following icon:

- ✓ The following happens:
 - a. The spindle moves to the tool change position.
 - b. A dialog window opens. DENTALCNC prompts you to confirm the opening of the collet chuck.

12. **CAUTION!** Wear gloves.
13. Open the working chamber door.
14. Confirm the current message.

- ✓ The following happens:
 - a. The collet chuck opens.
 - b. The current dialog window closes.
 - c. A dialog window opens.

15. Insert the tool into the collet chuck and hold it in place. Be sure to insert the tool until the tool ring reaches the tip of the collet chuck.


FIG. 43 INSERTING THE TOOL INTO THE COLLET CHUCK

16. Confirm the current message.

- ✓ The following happens:
 - a. The collet chuck closes. You can release the tool.
 - b. The current dialog window closes.
 - c. A dialog window opens.

17. Close the working chamber door.

18. Confirm the current message.

19. Activate the cooling liquid pump by selecting the depicted icon in the **Machin-** **ing** view in DENTALCNC.

20. Wait until the cooling liquid spray is free of air and all jets hit the tool tip ([Fig. 44](#) on the facing page). This can take several minutes.


FIG. 44 CORRECT COOLING LIQUID SPRAY (LEFT) AND SPRAY WITH INSUFFICIENT PRESSURE (MIDDLE, RIGHT)

21. Deactivate the cooling liquid pump by selecting the depicted icon in the **Machining** view in DENTALCNC.

7 OPERATION: EXECUTING JOBS

You control and start manufacturing with DENTALCNC. In this chapter, we will give you a brief overview. For the complete instructions, see the documentation for DENTALCNC.

NOTICE

Damaging of the machine when using damaged tools or blanks

If tools or blanks are damaged, parts can break off and damage the machine during job execution.

» Check the blanks and tools **thoroughly** for damage before every job execution.

When you have prepared your jobs and your machine, you can start machining. Machining is a fully automated process and only requires your attention in case of unexpected events.

! Do not move the machine during job execution, otherwise the results may become imprecise.

7.1 Starting jobs

- Ensure the following:
 - You created a job on your CAM computer. It has been transferred to DENTALCNC.
 - All required tools are in the right positions in the tool magazine and are neither worn nor damaged. They have also been added to the virtual tool magazine in DENTALCNC.
 - The required blanks are mounted in the blank magazine and are assigned in DENTALCNC.
 - For wet machining: The wet grinding module is properly installed and there is enough cooling liquid in the tank.
 - For wet machining: The ionizer cover is installed.
 - For dry machining: The ionizer cover is removed.
 - The compressed air supply is set correctly.
- Close the working chamber door.
- If you manually control the suction device, switch it on and set it to the required level.
- Start machining via the depicted icon in DENTALCNC.

7.2 Aborting machining

You can abort machining as follows:

- Select one of the following icons:

Abort machining immediately

Abort machining after the current job (only when using the blank changer)

- Confirm the current message.

✓ The following happens:

- Machining stops according to your selection.
 - You are prompted to select whether the tool in the spindle can be inserted safely into the tool magazine.
- If you want the machine to automatically put the tool into the tool magazine, confirm the dialog. If you want to remove the tool manually from the collet chuck, answer the question in the negative.

Removing the tool manually is necessary in case that putting it into the tool magazine would damage the tool magazine inserts.

- ✓ Depending on your choice, the spindle either puts the tool into the tool magazine or moves into the tool change position. In the latter case, continue with the next step.

CAUTION

Danger of cuts and burns when touching tools with your bare hands

If you handle tools on the cutting surface, you may be injured. As the tool may be very hot, you may also suffer from skin burns.

- » Only touch tools at their shank.
 » When handling tools, wear protective gloves.

- CAUTION!** Wear gloves.
- Open the working chamber door.
- Hold the tool in the collet chuck in place.
- Confirm the current message.

✓ The following happens:

- The collet chuck opens.
 - The current dialog window closes.
 - A dialog window opens.
- Remove the tool from the collet chuck.


FIG. 45 REMOVING THE TOOL FROM THE COLLET CHUCK

9. Close the working chamber door.
10. Confirm the current message.
 - ✓ The following happens:
 - a. The collet chuck closes.
 - b. If your machine is equipped with an automatic working chamber door, the door automatically closes.
 - c. The spindle moves to its default position.

7.3 Job interruptions and job abortions

A job will be *interrupted* in the following cases:

- The compressed air supply is not sufficient
- The vacuum in the working chamber is not sufficient
- For wet machining: The cooling liquid flow rate is not sufficient

An *interrupted* job will normally be continued automatically after the error is corrected.

A job will be *aborted* in the following cases:

- In case of a machine malfunction
- In case of a tool breakage
- In case of a power failure

If a job was *aborted*, you have to restart it.

i DENTALCNC allows you to resume the job from the last machining step. For more information, see the corresponding documentation.

7.3.1 How to proceed in case of a job interruption

If the job was interrupted, DENTALCNC displays a corresponding message.

If the compressed air is insufficient

» Check the following:

- The manometer of the compressed air regulator
- The installation of the pneumatic hoses
- Your compressor

If the vacuum is insufficient

» Check the suction hose and your suction unit.

If the cooling liquid flow rate is insufficient

1. Check if the wet grinding module is properly installed and operational.
2. Check if there is enough cooling liquid in the tank
3. Check if the cooling liquid, the tank and the filter are clean.

7.3.2 How to proceed in case of a machine malfunction

A machine malfunction is recognized by the internal control unit in case of a critical event. The working chamber will be illuminated in red. DENTALCNC displays the error message and error code that was sent by the control unit.

1. Record the error message and error code that is displayed.
2. Restart the machine and the CAM computer. If the problem persists, continue with the next step.
3. Disconnect the machine from the electrical source and prevent it from being restarted.
4. Contact customer service. Have the error message and error code readily available.
5. If you need to remove a blank from the working chamber, perform an emergency opening of the working chamber door.

7.3.3 How to proceed in case of a tool breakage

If a tool breaks during machining, the machine will not recognize this immediately. Instead, the spindle will continue to move with the broken tool. The tool breakage will be recognized upon the following events:

- The next regular tool change

A tool breakage can be caused by the following:

- The tool was damaged or worn
- The tool was put into the wrong tool position or was manually inserted into the spindle at the wrong time. As a consequence, it was not suitable for the processing step.
- The distribution of the objects in the blank (“nesting”) was not correct.

If a tool breaks, do the following:

1. Open the working chamber door.
2. Remove all parts of the broken tool from the working chamber and the collet chuck.
3. Re-add the tool to the virtual ATB *and* virtual tool magazine in DENTALCNC.
4. If the spindle picked up the tool from the tool magazine, check if the tool was inserted into the correct position. Insert a spare tool into the correct position in the tool magazine.
5. If you manually inserted the tool into the collet chuck, check if the broken tool corresponds to the tool type which you were prompted to insert. Have a correct spare tool ready.
6. Close the working chamber door.
7. Resume the job.

 DENTALCNC allows you to resume the job from the last tool change. For more information, see the corresponding documentation.

 If tools break regularly, you can find additional information in the troubleshooting section:
[🔗 Troubleshooting – on page 56](#)

7.3.4 How to proceed in case of a power failure

As long as the machine is not powered, you do not have access to the working chamber.

- » After a short power failure, restart the machine and the CAM computer.
- » If you need to access the working chamber in case of a longer power failure, perform an emergency opening of the working chamber door.

If the blank changer is moving when a power failure occurs, it may not move back to the default position. This will produce an emergency stop signal when the machine is started and the machine will not reference.

- » If the blank changer does not move and the machine does not reference, perform a manual reset of the blank changer to default position. [🔗 Manually resetting the blank changer to default position – on page 42](#)

7.4 Emergency opening of the working chamber door

WARNING

Crushing hazard and cutting injuries if the working chamber door is open

If the working chamber door is open during machining, it will not protect users from bruises and cuts.

- » **Do not** open or close the working chamber door during machining.
- » **Never** operate the machine with the working chamber door open.
- » Perform the emergency opening only if you are authorized to do so and if you have received training on it.
- » Activate the safety interlock **immediately** after completing your work in the working chamber.

CAUTION

Cutting injuries when touching a rotating tool

If a power failure or a machine malfunction occurs during machining, the spindle including the inserted tool keeps rotating. If you touch the rotating tool, you will suffer from cutting injuries.

- » Wait until the spindle has stopped rotating before performing an emergency opening.

The working chamber door is equipped with a safety interlock that only allows you to open the door when the machine is powered. You can use the emergency release of the working chamber door if the working chamber door is permanently locked.

Location of the opening for the emergency release:


FIG. 46 OPENING FOR THE EMERGENCY RELEASE

You can perform an emergency opening as follows:

1. Switch off the machine at the main power switch. Disconnect the machine from the electrical source.

2. Open the blank magazine door.
3. Deactivate the safety interlock of the working chamber door, as shown in the figure (↗ Fig. 47 see below):
 - a. Insert the key for the emergency release of the working chamber door through the opening.
 - b. Deactivate the safety interlock of the working chamber door by turning the key counter-clockwise by 90°


FIG. 47 USING THE KEY FOR THE EMERGENCY RELEASE

*Top: Remove the cap
Middle: Insert the key
Bottom: Turn the key*

4. Open the working chamber door.
5. Carry out your work in the working chamber.
6. Close the working chamber door.
7. Reactivate the safety interlock of the working chamber door by turning the key for the emergency release clockwise by 90°.
8. Check if the working chamber door actually is locked.
9. If you can still open the working chamber door, repeat the 3 previous steps.
10. Close the opening for the emergency release again with the protective cap.
11. Close the blank magazine door.

7.5 Manually resetting the blank changer to default position

Cutting injuries when touching a rotating tool

If a power failure or a machine malfunction occurs during machining, the spindle including the inserted tool keeps rotating. If you touch the rotating tool, you will suffer from cutting injuries.

» Wait until the spindle has stopped rotating before performing an emergency opening.

In some cases (e. g. when a power failure occurs) the blank changer may not be able to move back into default position. This will produce an emergency stop signal when the machine is started and the machine will not reference.

In this case a manual reset of the blank changer to default position is required.

1. Open the working chamber door. If necessary, perform an emergency opening of the working chamber door. [Emergency opening of the working chamber door – on the previous page](#)
2. Switch off the machine at the main power switch. Disconnect the machine from the electrical source.
3. Disconnect the machine from the external compressed air supply.
4. Open the blank magazine door.
5. With one hand open the blank changer flap (marked orange) in the working chamber.
6. With the other hand push the blank changer carriage (marked blue) towards the blank gripper (marked green).


FIG. 48 BLANK CHANGER FLAP (MARKED ORANGE), BLANK CHANGER CARRIAGE (MARKED BLUE) AND BLANK GRIPPER (MARKED GREEN)

7. If there is a blank frame in the blank gripper, open the gripper and remove it.
8. Close the blank gripper.

9. Push the blank changer carriage completely towards the blank magazine.
10. Push the blank gripper into the end position at the left side of the machine.
11. Open the blank gripper.
 - ✓ The blank changer is in the default position.
12. Close the working chamber door.
13. If you performed an emergency opening, reactivate the safety interlock of the working chamber door. [Emergency opening of the working chamber door – on the previous page](#)
14. Close the blank magazine door.
15. Connect the machine to the external compressed air supply.
16. Connect the machine to the electrical source
17. Switch on the machine at the main power switch.
18. Start the CAM computer.
 - ✓ The machine references.

8 MAINTENANCE AND DO-IT-YOURSELF

Some day-to-day basic maintenance and preventive maintenance is essential to keep the machine mechanics and electric components in good condition for proper machining results.

It is your responsibility to make sure that preventive maintenance, as well as basic maintenance, is performed.

You are the only one who can ensure that your machine receives the proper maintenance care. You are a vital link in the maintenance chain.

8.1 Basic maintenance

Basic maintenance includes tasks which are part of everyday operation. You are responsible to ensure that these tasks are carried out according to the maintenance table. You only need minimal manual skills for these tasks and most required tools are provided with the machine.

 A maintenance kit is available for your machine. It contains the required items for the maintenance tasks that you need to perform after 500 and 1,000 operating hours and for calibrating the machine. The maintenance kit is available from customer service.

8.2 Maintenance section

For your convenience, DENTALCNC lists all basic maintenance tasks in the **Maintenance** section. In the **Maintenance** section, you can see when the individual tasks are due.

When you have marked a maintenance task as complete, its time interval is reset and the list is updated accordingly.

» When you have performed a maintenance task, mark it as complete in the **Maintenance** section. This will keep the list up-to-date.

8.3 Preventive maintenance

Preventive maintenance for this machine has to be scheduled every 2 years, or after 2,000 operating hours at the latest.

» To schedule preventive maintenance, contact customer service.

8.4 Where to get service?

Customer service is your main contact for all service related questions. They will provide you with spare

parts, maintenance tips and they will perform preventive maintenance for you on request.

» When your machine is delivered or installed, ask the service technician for contact details for the customer service team. We also recommend scheduling the first preventive maintenance appointment at this point to ensure your machine gets proper maintenance.

8.5 Definition of wear parts

The machine and the extra equipment are warranted for a period of 24 months or 2,000 operating hours, whatever comes first. The warranty covers defects in materials or fabrication as long as the regulations for using the machine in all documents are followed.

Of course, the warranty also covers wear parts as long as their failure cannot be attributed to the function-related abrasion. The wear parts that are mentioned in the maintenance table can already wear down within the warranty period due to their normal function. The average useful life of the wear parts can be seen in the maintenance table.

Use these values to determine operating costs, to plan your spare part stock as well as to create individual maintenance and service plans.

Maintenance table

[page 53](#)

8.6 Cleaning the working chamber

Cleaning the working chamber includes the following components:

- Measuring key
- Blank holder
- View window

NOTICE

Damaging of the linear guides or the spindle when cleaning with compressed air

If you clean the working chamber with compressed air, material chips can reach the linear guides or the spindle bearings.

» **Never** clean the working chamber with compressed air.

CAUTION

Breathing difficulties caused by processing dust

Processing dust that gets into your lungs can cause breathing difficulties.

» Clean the machine only if the air extraction system is properly installed and activated.

» Wear a face mask of class FFP2 during the entire cleaning

We recommend cleaning the working chamber after all other required maintenance tasks.

1. Have ready:
 - A wet cloth
 - A mild cleaning agent (optional)
 - A dry cloth (if the machine was used for dry machining only)
 - A vacuum cleaner (if the machine was used for dry machining only)
 - A wet brush for the measuring key
2. Close the working chamber door.
3. Move the axes into the cleaning position by selecting the depicted icon in the **Machining** view in DENTALCNC. 
4. Open the working chamber door.
5. Mount the ionizer cover.
6. If the machine was used for dry machining only:
 - a. Suck up the coarse dirt in the working chamber with the vacuum cleaner.
 - b. Clean all surfaces and cracks in the working chamber thoroughly with the dry cloth. Use a mild cleaning agent if necessary.

7. If the machine was used for wet machining: Clean all surfaces and cracks in the working chamber thoroughly with a wet cloth. Use a mild cleaning agent if necessary.
8. Clean the measuring key with the wet brush:
 - a. Clean every opening of the protective cage (marked orange) with the wet brush.
 - b. Clean the measuring key from all sides with the wet brush, reaching into the openings of the protective cage.
 - c. Clean the protective cage with a cloth.


FIG. 49 MEASURING KEY (MARKED BLUE) WITH PROTECTIVE CAGE (MARKED ORANGE)

9. Thoroughly clean the blank holder from all sides with a brush. Especially clean all openings and movable parts of the blank holder.
10. If necessary, clean the view window with a dry cloth. Use a wet cloth if the machine was used for wet machining. Use a mild cleaning agent if necessary.
11. Allow the working chamber to dry.
12. Close the working chamber door.
13. Move the axes into the default position by selecting the depicted icon in the **Machining** view in DENTALCNC. 

8.7 Cleaning the blank magazine

Cleaning the blank magazine is important because otherwise machining debris will damage sensitive machine parts.

1. Open the blank magazine door.
2. Remove all blank frames and block holders from the blank magazine.
3. Remove all blanks.
4. Clean the blank magazine and the interior of the blank changer with a vacuum cleaner and a wet cloth. Use a mild cleaning agent if necessary. Lift the blank changer flap in the working chamber to reach the blank changer carriage.
5. Thoroughly clean all blank magazine slots, blank frames and block holders with a brush.
6. Equip the blank frames and block holders and insert them into the blank magazine. Ensure that the blank magazine slots correspond to the slots in DENTALCNC.

8.8 Cleaning the collet chuck

You need to clean the collet chuck with the spindle service set which was provided with your machine.

NOTICE

Damaging of the spindle when cleaning with compressed air

compressed air

If you clean the collet chuck with compressed air, the spindle bearings can be damaged.

» Clean the collet chuck **only** with the appropriate service set.


FIG. 50 SPINDLE SERVICE SET

1. Cleaning brush
2. Knurled nut
3. Tube of collet chuck grease
4. Cleaning cone

To clean the collet chuck do the following:

1. Have the spindle service set ready.
2. Close the working chamber door.
3. Open DENTALCNC's general application settings with the depicted icon in the main icon bar. 
4. Open the **Machine parameters** view with the depicted icon in the local icon bar. 
5. Move the axes into the tool change position by selecting the depicted icon in the **Machining** view in DENTALCNC. 

- ✓ The spindle moves into the position for collet chuck cleaning.
6. Open the working chamber door.
 7. Untighten the collet chuck with the knurled nut:
 - a. Insert the pin of the knurled nut into the collet chuck.
 - b. Ensure that the collet chuck fully sits in the recess of the knurled nut.


FIG. 51 RECESS IN KNURLED NUT (MARKED ORANGE)

- c. Twist the knurled nut counter-clockwise.


FIG. 52 UNTIGHTENING THE COLLET CHUCK WITH THE KNURLED NUT

8. Unscrew and remove the collet chuck with your hand.


FIG. 53 UNSCREWING THE COLLET CHUCK (MARKED ORANGE IN LEFT ILLUSTRATION) WITH YOUR HAND

9. Put the knurled nut and the collet chuck aside within reach.
10. Clean the inner cone of the spindle with the cleaning cone of the service set.


FIG. 54 CLEANING THE INNER CONE OF THE SPINDLE

11. Clean the collet chuck with the brush of the service set.
 - a. Insert the brush into the collet chuck until the end of the brush aligns with the end of the collet chuck.

! Do not push any further because otherwise dirt may get into the thread of the collet chuck.

FIG. 55 THE END OF THE BRUSH MARKED WITH A VERTICAL ORANGE LINE

- b. **Rapidly** pull the brush out of the collet chuck and twist it at the same time.


FIG. 56 CLEANING THE COLLET CHUCK

- c. Repeat the last 2 steps multiple times.

NOTICE Damaging of the spindle when using the wrong grease or applying the grease incorrectly

If you use unsuitable grease or if grease gets into the longitudinal slots of the collet chuck, the machine may get damaged.

- » Before applying the grease, ensure that the collet chuck is perfectly clean.
- » Ensure that no grease gets into the longitudinal slots of the collet chuck.
- » Only use a very small, less than pinhead-size amount of the grease.
- » Only use the provided grease of the service set.

12. Thoroughly clean the complete outer surface of the collet chuck with a clean and dry cloth.


FIG. 57 CLEANING THE COMPLETE OUTER SURFACE OF THE COLLET CHUCK WITH A CLOTH

13. Put a small amount of the collet grease on the index finger and smear it with the thumb.
14. Apply the smeared collet grease to the flanks of the collet chuck.


FIG. 58 GREASING THE COLLET CHUCK; SURFACE IN WHICH GREASE IS TO BE APPLIED IS MARKED IN BLUE; SLOT THAT MAY NOT BE GREASED IS MARKED IN RED

15. Insert the collet chuck into the spindle. Screw the collet chuck in by twisting it counter-clockwise with your hand.


FIG. 59 SCREWING IN THE COLLET CHUCK (MARKED ORANGE IN THE RIGHT ILLUSTRATION) INTO THE SPINDLE WITH YOUR HAND

16. Insert the pin of the knurled nut into the collet chuck. Ensure that the collet chuck is properly seated in the recess of the knurled nut.
17. Tighten the collet chuck with the knurled nut:
 - a. Insert the pin of the knurled nut into the collet chuck.
 - b. Ensure that the collet chuck fully sits in the recess of the knurled nut.


FIG. 60 RECESS IN KNURLED NUT (MARKED ORANGE)

- c. Twist the knurled nut clockwise.

! The collet chuck must be screwed in tightly. Otherwise rotational imperfections can occur during operation which will worsen your processing results


FIG. 61 TIGHTENING THE COLLET CHUCK WITH THE KNURLED NUT; RECESS IN KNURLED NUT MARKED ORANGE

18. Store the spindle service set in a safe place.
19. Close the working chamber door.
20. Clean the nozzle plate.

8.9 Checking the compressed air regulator

NOTICE

Damaging of the machine when compressed air is contaminated

Compressed air that does not fulfill the guidelines for purity according to ISO 8573-1 can damage the machine.

- » Check the water separator of the compressed air regulator daily for contamination.
- » **Never** use the machine if there is water, oil or solid particles in the water separator.

8.9.1 Checking the water separator for condensate

Condensate in the separator usually points to compressed air not being dry enough.

1. Check if water, oil or solid particles piled up in the water separator.
2. If this is the case, switch the machine off immediately and proceed as follows.
3. Disconnect the machine from the compressed air supply.
4. Check the compressed air supply and make sure that the compressed air fulfills the requirements for air purity according to ISO 8573-1. Do not use the machine until the compressed air fulfills this requirement.
5. Drain the water separator by turning the discharge screw counter-clockwise.
- ✓ The condensate is blown out downwards under pressure.
6. Close the discharge screw again by turning it clockwise.


FIG. 62 OPENING / CLOSING THE DISCHARGE SCREW OF THE COMPRESSED AIR REGULATOR

8.9.2 Exchanging / cleaning the contaminated filter cartridge

You have to clean or exchange the filter cartridge in the water separator in case of strong contamination.

! A strongly contaminated cartridge can lead to a pressure loss.

If the compressed air fulfills the requirements for air purity according to ISO 8573-1, the filter cartridge usually does not have to be changed.

» If the filter cartridge is contaminated, check the purity of your compressed air.

You exchange or clean the filter cartridge as follows:

1. Disconnect the machine from the compressed air supply.
2. To remove the bowl of the water separator, do the following:
 - a. Press the release and turn the bowl counterclockwise as far as it will go.
 - b. Pull the bowl off the housing.


FIG. 63

3. Unscrew the filter screw below the filter cartridge.
4. Pull out the filter cartridge and clean it if necessary.

i A new filter cartridge is available as spare part from customer service.

5. Insert the new or cleaned filter cartridge and reassemble the water separator. While doing so, insert the release into the recess in the housing and turn the bowl until it locks into place.

8.10 Cleaning the housing

NOTICE

Damaging of the housing when using an unsuitable cleaning agent

If you use an unsuitable cleaning agent and/or cleaning tool for cleaning the machine housing, the surface or the adhesive foil can get damaged.

- » To avoid scratches, only use a microfiber cloth to clean the housing.
- » Take care of the adhesive symbols so that they do not peel off. The adhesive foil is especially sensitive to rubbing and strong cleaning agents.
- » If the usage of a special cleaning agent is necessary to eliminate certain dirt, we recommend checking the suitability of the cleaning agent at a hidden place of the part first.

1. Clean the surface with a dry microfiber cloth.
2. If some dirt cannot be removed this way, moisten the cloth. Use a pH neutral cleaning agent if necessary.

8.11 Exchanging the main fuse

The internal power supply of the machine has a main fuse that is accessible from the outside and can be replaced if necessary.

» As a replacement fuse, only use a fuse of the following type: T6,3A L250V

i A new main fuse is available as spare part from customer service.

1. Switch off the machine at the main power switch.
2. Remove the power cord from the connection panel.
3. Remove the cover of the fuse.


FIG. 64 COVER OF THE FUSE (MARKED ORANGE)

4. Remove the defective fuse and replace it with a new fuse.
5. If you do not have a replacement fuse ready, take the replacement fuse from the right side of the fuse cover and put it into the left side of the fuse cover.
6. Remount the fuse cover.

8.12 Calibrating the axes

NOTICE

Deterioration of machining results caused by an incorrect calibration

calibration

At delivery, your machine is already calibrated. As long as your machining results are accurate, a new calibration is not necessary. A calibration takes much time and will deteriorate the machining results if it is improperly executed.

- » In case of inaccurate machining results, try adjusting the machining conditions first: Check the fixation and quality of the blank and the state of the tool.
- » **Before** calibrating the machine, contact customer service.
- » **Be very careful** when measuring and entering data during calibration. When in doubt, abort the calibration.

By calibrating the machine with test and calibration specimens, it may be possible to improve the machining results.

i The documentation for the manufacturing software contains all information on calibrating the machine. Therefore, you will only find information specific to this machine in this document.

Your machine is delivered with a calibration set. It contains the following parts:

- Calibration blanks which are used to mill calibration and test specimens
- A tool for milling the calibration or test specimens
- A micrometer for measuring the machining precision

You can calibrate the machine as follows:

1. Have the calibration set readily available.
2. Mount the calibration blank into the blank holder.
3. Follow the instructions on calibrating the machine in the documentation for DENTALCNC.
4. Keep all parts of the calibration set except used calibration blanks for further use.

8.13 Replacing the tool magazine inserts

When tool magazine inserts are worn, they should be replaced. New inserts are delivered without holes for the tools. The holes must be drilled into the inserts with the machine.

- Your machine comes with tool magazine inserts as spare parts and with the drill tool.
- Additional inserts and drill tools are available via customer service.

i The documentation for the manufacturing software contains step-by-step instructions for drilling the holes into the inserts. Below you find a description of how to replace the tool magazine insert in the machine.

You can replace tool magazine inserts as follows:

1. Have the spare tool magazine inserts readily available.
2. Open the working chamber door.
3. Remove all tools from the tool magazine.
4. Unscrew the 3 screws on the upper side of the tool magazine and lift the cover.


FIG. 65 REMOVING THE COVER OF THE TOOL MAGAZINE

5. Remove the 2 existing tool magazine inserts and replace them with 2 new ones.


FIG. 66 NEW TOOL MAGAZINE INSERTS (RIGHT, WITHOUT TOOL POSITIONS) REPLACE OLD INSERTS (LEFT, WITH TOOL POSITIONS)

6. Put the cover back onto the tool magazine and screw it down.
7. Follow the instructions in the documentation for DENTALCNC and drill the tool positions using the provided drill bit.

8.14 Updating the software and firmware

Updating the manufacturing software and the firmware is an interdependent procedure.

» Be sure to read all information in this chapter.

NOTICE

Damaging of the control unit when a firmware update is interrupted

If the firmware update is interrupted, the control unit of the machine may become permanently damaged.

- » Only update the firmware if a permanent electrical supply of the machine and the computer is guaranteed.
- » Only update the firmware if the computer's operating system is stable and is free of malware.
- » Only update the firmware if the connection between the computer and the machine is stable. Always use a wired connection during the update.
- » Do not disconnect the machine or the computer from the electrical source nor switch off the machine or computer during a firmware update.
- » Do not close DENTALCNC during a firmware update.

8.14.1 Updating the manufacturing software

DENTALCAM and DENTALCNC are updated regularly.

There are 2 ways to update your existing installation:

- Automatically via the update function in DENTALCNC (recommended)
 - Manually via an installation package. Ask customer service for the latest installation package.
- » To update DENTALCAM & DENTALCNC, follow the instructions in the documentation for the manufacturing software.

8.14.2 Updating the firmware of the machine

The firmware is the internal control software of your machine. New versions may introduce new functions and improve existing ones. New firmware versions come as part of new versions of DENTALCNC.

New firmware versions will be automatically installed after a new DENTALCAM & DENTALCNC version has been downloaded via the automatic update function. In the following cases, you need to update the firmware manually:

- You do not use the automatic update function.

- No connected machine was found during the automatic download of the new DENTALCAM & DENTALCNC version.

You can manually update the firmware of your machine as follows:

If DENTALCNC contains a new firmware, you will be prompted to update the firmware after DENTALCNC connected to the machine or when you start machining. The firmware update should take 5 – 15 minutes to complete.

1. To update the firmware, confirm the message. You can cancel the message to postpone the update, but we recommend updating the firmware immediately.

✓ A window opens.

2. In the current window, select the following icon:

✓ The following happens:

- a. DENTALCNC updates the firmware of the connected machine.
- b. DENTALCNC prompts you to restart the machine.

3. Switch off the machine at the main power switch.

4. Wait 10 seconds.

5. Switch on the machine at the main power switch.

✓ The machine references. The machine is ready for operation.

8.15 Maintenance table

- !** If your machine is used for wet machining, additional maintenance is required that is not listed in this table. The maintenance tasks are described in the operating instructions for the vhf wet grinding module. The operating instructions are delivered with the vhf wet grinding module. They are also available for download at: <https://www.dentalportal.info/NO1>

Several times per day

Task	Recommended interval	Procedure / Parts	Tool illustration
Checking the external compressed air lines for damage	Before switching on the machine	Visual check	
Checking the suction hose for damage	Before switching on the machine	Visual check	

Once per day

Task	Recommended interval	Procedure / Parts	Tool illustration
Checking the compressed air regulator (page 48)	Before work Clean or exchange cartridge in case of visible contamination or every 2 years	If necessary, new filter cartridge	
Daily cleaning of the working chamber (page 44)	After work If soiled	Vacuum cleaner, Brush, Dry cloth	

Once per week

Task	Recommended interval	Procedure / Parts	Tool illustration
Cleaning the collet chuck (page 46)	Once per week In case of rotational imperfections	Spindle service set	
Cleaning the blank magazine (page 45)	Once per week If soiled	Wet cloth, Dry cloth	

When necessary

Task	Recommended interval	Procedure / Parts	Spare part illustration
Updating the software and firmware (page 52)	When an update is available		
Cleaning the housing (page 49)		Microfiber cloth, Water, Mild cleaning agent (optional)	
Exchanging the main fuse (page 50)		Replacement fuse T6,3A L250V	

Wear parts that you can exchange yourself

Wear part	Recommended interval	Procedure / Parts	Spare part illustration
Collet chuck (page 46)	1,000 operating hours* Every year*	Removing & inserting the collet chucks (old/new) as during cleaning	
Tool magazine inserts (page 49)	500 operating hours* Every year*		

Wear parts that customer service exchanges for you

Wear part	Recommended interval	Procedure / Parts	Spare part illustration
Spindle bearings (requires spindle exchange by customer service)	2,000 operating hours*		

*These are recommendation guidelines. Depending on the processing material and how well the machine is cleaned, these values may differ.

9 DISPOSAL

9.1 Disposing of the cooling liquid

The cooling liquid must not be disposed of via the sewers. When disposing of the cooling liquid / machining residues, obey the following regulations.

- » Avoid entry of hazardous machining residues into the soil, water or into sewers.
- » Dispose of the machining residues as described by the manufacturer of the material.
- » Obey the national and local laws of the disposal location in any case.
- » If necessary, have the machining residues disposed of by an approved disposal company.
- » Have the cooling liquid and machining residues disposed of by an approved disposal company.
- » Keep a reference sample of the disposal product for at least 6 months.

9.2 Disposing of the machine

The machine must not be disposed of with the residual waste. This is indicated by the icon  which depicts a crossed out trashcan. In the European Union (EU), this is in accordance with Directive 2012/19/EU.

We will dispose of the machine at no cost. The owner will bear the costs for disassembly, packaging and transport.

- » Before sending in the machine for disposal, contact customer service.
- » If you dispose the machine yourself, obey the national and local laws of the disposal location.
- » If necessary, have the machine disposed of by an approved disposal company.

Dismantling, transport and packaging

[page 10](#)

Copyright

Distribution or duplication of all content is only allowed by written consent of vhf camfacture AG. This includes reproduction by presentation and broadcasting.

This document is published by:

vhf camfacture AG

Lettenstraße 10

72119 Ammerbuch, Germany

dentalportal.info

10 TROUBLESHOOTING

In case something is not working as intended, take a look at the following troubleshooting guide.

NOTICE

Machine damage due to improper troubleshooting

In case of improper troubleshooting, your machine may get damaged.

- » If you are unsure of how to perform certain steps during troubleshooting or cannot solve the problems, abort the troubleshooting and contact customer service.

Additional symbols in this chapter

- ❓ Question to narrow down the problem
- 💬 Suggested solution

I cannot open the working chamber door

❓ Is the machine operating?

While the axes are moving, you cannot open the working chamber door.

💬 If applicable:

- » Wait until the machine has finished.

❓ Has a power failure occurred at the installation site of the machine?

💬 If applicable:

- » Depending on the duration of the power failure, restart the machine or perform an emergency opening.

❓ Is electricity available at the installation site of the machine?

💬 If applicable:

1. Connect the machine to the electrical source.
2. Switch on the machine at the main power switch.
3. If the working chamber lighting does not illuminate, check if the power cable is properly seated in the machine and is connected to the electrical source.
4. Try connecting the machine to a different socket.

❓ Is the door blocked?

💬 If applicable:

1. Switch on the machine at the main power switch.
2. If the working chamber lighting does not illuminate, check if the power cable is properly seated in the machine and is connected to the electrical source.
3. If the machine is not powered, perform an emergency opening when necessary.

I have installed all components, started the software but the machine does not reference

❓ Is the working chamber door open?

The machine does not reference with the working chamber door open.

💬 If applicable:

- » Close the working chamber door.

❓ Is an Ethernet / USB cable connected to the machine?

💬 If applicable:

- » Check if the Ethernet / USB cable is properly seated in the connector and is undamaged. If possible, use the provided cable.

❓ Is the working chamber illuminated in red?

In this case a machine malfunction occurred.

💬 If applicable:

1. Restart the machine.
2. If the working chamber continues to be illuminated in red, contact customer service.

The machining results are not satisfactory and / or tools keep breaking

❓ Do the tool positions in the Tools section correspond to the tools in the respective tool magazine?

If not, the machine uses the wrong tools during job execution.

💬 How to check this:

1. In DENTALCNC, compare the tool positions in the with the tools in the corresponding tool magazine.
2. Replace the wrong tools in the tool magazine with the correct ones.

❓ Is the blank properly mounted?

💬 How to check this:

- » Remove the blank and remount it. If using blocks: The groove of blocks has to properly sit on the corresponding positioning pin.

❓ Are bolts, fixing mechanisms, gaps and openings of the blank holder contaminated by processing dust?

💬 If applicable:

- » Clean the mentioned components thoroughly.

❓ Is the measuring key contaminated?

💬 If applicable:

- » Clean the measuring key with a brush.

❓ Are the tools worn?

💬 How to check this:

1. Visually inspect all tools.
2. Check the tool life values in DENTALCNC.
3. Replace worn tools with new ones.

? **Do rings on the tools sit in the groove on the tool shank?**

How to check this:

- » Visually inspect all tools and push rings that have moved into the groove again.

? **Are the tool magazine inserts worn?**

If applicable:

- » Replace the tool magazine inserts with new ones.

? **Do the parameters of the job in the software correspond to the parameters of the blank?**

How to check this:

- » Make sure that the following parameters of the job and of the blank in the machine correspond to each other. Also make sure that they are suitable for the objects that you want to machine.

- Material type
- Blank dimensions
- Indications (types) of the individual objects

? **Do you use the latest version of the manufacturing software that is released for the machine?**

- » To update DENTALCAM & DENTALCNC, follow the instructions in the documentation for the manufacturing software.

? **Are the object files of sufficient quality?**

How to check this:

1. Check the quality of object files (STL files) in your CAD application or an STL viewer. Especially, consider the manufacturer information on wall thickness and border thickness.
2. If necessary, adjust the settings of your scanner and scan application.

? **Is the collet chuck of the spindle contaminated or is it loosely seated in the spindle?**

If applicable:

1. Clean the collet chuck with the provided spindle service set.
2. When you insert the collet chuck into the spindle, make sure that it is properly seated.

? **Did you exchange the collet chuck at the recommended interval?**

How to check this:

- » Verify the recommended interval for exchanging the collet chuck in the maintenance table. If necessary,

exchange the collet chuck.

The computer displays that the flow rate is too low

i You can resolve this problem while the job execution is interrupted. DENTALCNC continues the job as soon as the problem is resolved.

? **Is cooling liquid lacking in the wet grinding module? Is the cooling liquid soiled?**

If applicable:

- » Refill the tank with clean cooling liquid.

If necessary:

- » Clean the tank. Fill in clean cooling liquid.

? **Is the coarse filter in the wet grinding module clogged?**

If applicable:

- » Clean the filter and the tank. Fill in clean cooling liquid.

? **Is the wet grinding module not properly installed or defective?**

If applicable:

- » Ensure that the wet grinding module is properly installed and operational.

? **Are the liquid nozzles in the working chamber clogged?**

If applicable:

- » Clean the nozzles with the inter-dental brush.

The computer interrupts the job and displays that the air pressure is too low

i You can resolve this problem while the job execution is interrupted. DENTALCNC continues the job as soon as the problem is resolved.

? **Is the compressed air regulator properly set?**

- » Set the air pressure at the compressed air regulator to a correct value.

[↗ Adjusting the air pressure with the compressed air regulator – on page 23](#)

? **Is the error caused by the external compressed air supply system?**

How to check this:

1. Close the main external compressed air supply valve.
2. Check if all pneumatic hoses are properly seated in their connections and are undamaged.
3. Check if your compressor is switched on and set correctly.
4. Open all required valves of your compressed air supply system.

? Does the air pressure fluctuate significantly so that jobs are frequently interrupted?

 If applicable:

1. Verify that your compressor can *permanently* generate the required air pressure and volume flow.

 Not every compressor is designed for commercial use with dental machines.

2. If necessary, replace your compressor with one that meets the recommendations.

 [Installing the pneumatics – on page 21](#)

The computer interrupts the job and displays that the vacuum is too low

 You can resolve this problem while the job execution is interrupted. DENTALCNC continues the job as soon as the problem is resolved.

? Is the suction unit switched on and operating?

 How to check this:

1. Check if the suction hose is properly seated in the corresponding opening and is undamaged.
2. If the machine controls the suction unit:
 - a. Check if the switching unit or data cable is properly installed.
 - b. Try operating the suction unit without the switching unit or data cable.
3. Switch on the suction unit.
4. Set the extraction level higher until the machine continues the job.

? Is the filter or container of the suction unit full?

 If applicable:

» Insert an empty filter into the suction unit or empty the container.

? Does the suction unit have an automatic tapping function?

 If applicable:

» Choose a shorter tapping interval.

I have exchanged the tool magazine inserts and now there are no holes for the tools anymore

Tool magazine inserts are delivered without drilled holes for tools. You drill them with the machine.

» Use DENTALCNC to drill new holes into the new inserts.

Index

A

Anti-soiling concept 13
 Aspiration 24
 Axes 14

B

Blank changer 29
 Change position of the blank magazine 30
 Manually resetting the blank changer to default position
 42

C

CAM computer 13
 Change position of the blank magazine 30
 Collet chuck 46
 Compressed air regulator 23, 48
 Connection panel 11

D

Disposal 55
 Machine 55

E

Electric connection 20
 Emergency opening of the working chamber door 41

F

Firmware update 52
 Front side of the machine 11

H

Hose connection 24

I

Identification plate 11, 14
 Implants
 Operating regulations 9

Installation 17

 Air extraction system 24
 CAM computer 25
 Compressed air regulator 23
 Pneumatic hose 22
 Pneumatics 21
 Switching unit 25

Installation site 18

Ionizer

 Deactivating and activating the ionizer 35
 Mounting and removing the ionizer cover 35

J

Job interruption 39

M

Machine

 Malfunction 39

Machine location 18

Main fuse 50

Maintenance 9

 Cleaning the collet chuck 46
 Cleaning the housing 49-50
 Cleaning the working chamber 44
 Compressed air regulator 48
 Exchanging the main fuse 50
 General information 43
 Preventive maintenance 43
 Replacing the tool magazine inserts 51
 Wear part 43

Maintenance table 53

Malfunction 39

Manufacturing software 9

Measuring key 44

Mounting & removing blanks 27

Multi-machine control 25

O

Operating temperatures 18

P

- Pneumatic connection 21-22
- Pneumatic hose 21-22
- Power connection 20
- Power failure 40
- Prefabricated abutments
 - Operating regulations 9
- Preventive maintenance 43

R

- Repackaging 10

S

- Scope of delivery 17
- Sound emission 14
- Spindle 9
- Starting jobs 38
- Starting the machine with a tool in the collet chuck 26
- Storage 10
- Suction unit 24
 - Requirements 24
- Switching unit 25

T

- Tool breakage 39
- Tool magazine inserts 51
- Transport 10
- Transport lock 20
- Troubleshooting 56

U

- Unattended operation 9

W

- Wear part 43
- Wet machining
 - Preparing the machine for wet machining 36
- Working chamber 12
 - Clean 44
 - Emergency opening 41
- Working chamber door 12

Original EC Declaration of Conformity

according to EC directive for machinery 2006/42/EC Annex II A

We,

vhf camfacture AG
Lettenstraße 10
72119 Ammerbuch
Germany

hereby declare expressly that the

Machine: CNC milling machine
Type: S5
Serial No.: S5ID100000000 – S5ID299999999

fulfills all the relevant provisions of the following directives:

- 2006/42/EC Machinery directive
- 2014/30/EU EMC directive

References of the applied harmonized standards according to article 7 passage 2:

- EN 614-1:2006 + A1:2009
- EN ISO 12100:2010
- EN ISO 16090-1:2018
- EN ISO 13849-1:2015
- EN ISO 13849-2:2012
- EN 60204-1:2018
- EN IEC 61000-3-2:2019
- EN 61000-3-3:2013
- EN 61326-1:2013
- EN 61326-2-1:2013

References of other standards:

- IEC 61010-1:2010 + COR:2011 + A1:2016, modified + A1:2016/COR1:2019

The manufacturer undertakes to electronically transmit relevant information on the machinery in response to a reasoned request by the national authorities. Person established in the community who is authorized to compile the technical file:

Dipl.-Ing. (FH) Frank Benzinger
Vorstandsvorsitzender / Chief Executive Officer (CEO)
vhf camfacture AG
Lettenstraße 10
D-72119 Ammerbuch

Ammerbuch, 2021/06/30

(Frank Benzinger, CEO)

Original Declaration of Conformity

according to Supply of Machinery (Safety) Regulations 2008

We,

vhf camufacture AG
Lettenstraße 10
72119 Ammerbuch
Germany

hereby declare expressly that the

Machine: CNC milling machine
Type: S5
Serial No.: S5ID100000000 – S5ID299999999

fulfills all the relevant provisions of the following directives:

Supply of Machinery (Safety) Regulations 2008
Electromagnetic Compatibility Regulations 2016

References of the applied harmonized standards:

- BS EN 614-1:2006 + A1:2009	- BS EN ISO 13849-2:2012	- BS EN 61326-1:2013
- BS EN ISO 12100:2010	- BS EN 60204-1:2018	- BS EN 61326-2-1:2013
- BS EN ISO 16090-1:2018	- BS EN IEC 61000-3-2:2019	
- BS EN ISO 13849-1:2015	- BS EN 61000-3-3:2013	

References of other standards:

- IEC 61010-1:2010 + COR:2011 + A1:2016, modified + A1:2016/COR1:2019

The manufacturer undertakes to electronically transmit relevant information on the machinery in response to a reasoned request by the national authorities. Person established in the community who is authorized to compile the technical file:

Dipl.-Ing. (FH) Frank Benzinger
Vorstandsvorsitzender / Chief Executive Officer (CEO)
vhf camufacture AG
Lettenstraße 10
D-72119 Ammerbuch

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