RAYPEX® 6

RAYPEX® 6

CE
0123

 Manufacturer
VDW GmbH • Bayerwaldstr. 15 • 81737 Munich • Germany
Phone +49 89 62734-0 • Fax +49 89 62734-304
info@vdw-dental.com • www.vdw-dental.com

M 5L RAY6 DFU PRT
Rev.3/31.05.2018

RAYPEX® 6 de | en | fr | it | es

Gebrauchsanweisung
Operating Manual
Mode d'emploi
Manuale d'uso
Manual de instrucciones

 VDW®

Endo Easy Efficient®

 VDW®

Endo Easy Efficient®

Inhalt | Contents | Contenu | Indice | Índice

Gebrauchsanweisung	4
Operating Manual	30
Mode d'emploi	56
Manuale d'uso	82
Manual de instrucciones	108
Annex Annex Annexe Allegato Anexo	134

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Electromagnetic compatibility (EMC)

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Congratulations on your purchase of the
RAYPEX® 6 touch apex locator.

en

Please do not hesitate to contact VDW GmbH for help with any doubt or problem related to this manual. Kindly keep this manual for further reference.

VDW GmbH reserves the right to change the information and data contained in these Directions For Use at anytime and without prior notice.

These Directions For Use are available in other languages upon request.

This booklet has been compiled with the greatest care, but even with all our efforts, mistakes can never be entirely excluded.

Suggestions for improvements are welcome at any time. In this case, kindly contact VDW GmbH directly.

VDW GmbH

Bayerwaldstr. 15
81737 Munich
Germany

Phone +49 89 62734-0
Fax +49 89 62734-304
info@vdw-dental.com
www.vdw-dental.com

Contents

1	Identification of Symbols	32	8	Maintenance, cleaning and sterilization.....	49
1.1	Symbols used in these Directions for Use.....	32	8.1	General information.....	49
1.2	Symbols used on Packaging, Device and Parts	32	8.2	Cleaning, Disinfection and Sterilisation ...	50
2	Indications for Use	33	8.2.1	Pre-treatment	50
3	Contraindications	33	8.2.2	Manual cleaning and disinfection	51
4	Warnings	33	8.2.3	Inspection/Maintenance	52
5	General Precautions.....	34	8.2.4	Packaging	52
6	Adverse Reactions	34	8.2.5	Sterilization.....	52
7	Step-by-Step Instructions	34	8.2.6	Storage	53
7.1	Standard Components.....	34	8.2.7	Material resistance.....	53
7.2	Installation	35	9	Conformity.....	53
7.2.1	Charger	35	10	Manufacturer	53
7.2.2	Rechargeable Battery.....	35	11	Warranty.....	54
7.3	Description of User Interface	36	12	Disclaimer	54
7.4	Operation.....	37	13	Technical Data	55
7.4.1	Connecting the Device	37	ANNEX		
7.4.2	Starting Length Determination.....	38	Electromagnetic Compatibility.....		
7.4.3	Apex Localisation	39	134		
7.4.4	Sound Level Selection	41			
7.4.5	Automatic Shutdown	41			
7.5	Tips for Successful Length Determination	42			
7.6	Device Setup	43			
7.6.1	Settings Menu Features.....	44			
7.6.2	Optional DR'S CHOICE Apical Line	44			
7.6.3	CHECK Mode.....	46			
7.6.4	DEMO Mode.....	48			


1 Identification of Symbols

en

1.1 Symbols used in these Directions for Use

 WARNING	If the instructions are not being followed properly, operation may result in hazards for the product or the user/patient.	 NOTE	Additional information, explanation on operation and performance.
---	---	--	---

1.2 Symbols used on Packaging, Device and Parts

 SN	Serial number		Refer to directions manual/booklet
	Manufacturer		Manufacturing date
	Class II product		Special disposal of waste electrical and electronic equipment (Directive 2002/96/EEC)
	Type BF applied part		Direct current (connection for power supply)
	Gost symbol, product is in conformity with Russian safety standards (GOST-R)		Caution
 REF	Catalogue number (re-order number)		CE marking
	Keep dry!		Fragile
	Temperature limit		Content (quantity)
	Atmospheric pressure limitation		Humidity limitation (Non-Condensing)

2 Indications for Use

FOR DENTAL USE ONLY!

RAYPEX® 6 is a microprocessor controlled device used for root canal length determination.

3 Contraindications

Use of the RAYPEX® 6 is contraindicated on patients or by operators having implanted electronic devices such as pacemaker, etc.

4 Warnings

In this chapter a description of serious adverse reactions and potential safety hazards for the product or the user/patient is included.

Please read the following warnings before use.

WARNINGS

- *The device may be used only by qualified dentists acting in accordance with the national regulations.*
- *The device may be used only in suitable locations and not outdoor.*
- *Make sure that the cords do not hinder the free passage of people.*
- *Do not expose the device to direct or indirect sources of heat. Store and use the device in a safe environment.*

- *Do not use the device in the presence of inflammable anaesthetic mixtures.*
- *Do not immerse in liquid.*
- *Use only original accessories.*
- *The external charger to which the device is connected should comply with standards in force.*
- *Do not use the device if it appears to be damaged or defective.*
- *Do not carry out repairs or modifications to the device without prior authorization by VDW GmbH. If any fault occurs, contact your local supplier rather than have it repaired by an unauthorized person.*
- *Do not connect the device to or use it in combination with any other apparatus or system.*
- *Do not use the device as an integral component of any other apparatus or system. VDW GmbH will not be responsible for accidents, equipment damage, bodily injury or any other trouble which results from ignoring this prohibition.*
- *In case of any doubt, contact your local supplier or the After Sales Service of VDW GmbH in Munich.*
- *The device complies with the Electromagnetic compatibility standard (IEC 60601-1-2). Nevertheless, make sure that any possible electromagnetic interference does not pose any additional risk.*

5 General Precautions

Read these safety precautions thoroughly prior to use. These precautions allow you to use the product safely, preventing harm to you and to others.

It is of the utmost importance that this manual is preserved for future consultation. The manual must accompany the system in all cases of sale or other transfer in order that the new owner can observe the precautions and warnings.

The manufacturer declines any responsibility in the case of:

- Use of the device for applications other than those specified in the directions for use.
- Modifications or repairs performed by persons not authorized by the manufacturer.
- Use of non-original components or components other than those specified in the STANDARD COMPONENTS chapter (7.1.).

6 Adverse Reactions

There are no known adverse reactions.

7 Step-by-Step Instructions

! WARNING

The device may be used only by qualified dentists acting in accordance with the national regulations.

Refer to the WARNINGS chapter (4) to verify any special care to exercise before starting to use the device.

When opening the package and prior to installation, check the device for completeness. Report any damage sustained during shipping or any missing parts to your retailer within 24 hours of receipt of the device.

Ambient Conditions for Operation

- Use: indoor
- Ambient temperature: +10°C to + 40°C (50°F to 104°F)
- Relative humidity: 10% to 90%
- Atmospheric pressure: 70 kPa to 106 kPa

! WARNING

Do not install the device in damp places or in places where it will come into permanent contact with liquids of any kind.


Fig. 1 Incorrect placement of device

7.1 Standard Components

- 1 RAYPEX® 6 apex locator
- 1 Charger
- 1 Tester for operation check
- 1 Accessory set including: 1 measurement cable, 2 lip clips, 2 file clips, 1 touch probe
- Directions For Use


🔑 NOTE

None of the accessories is delivered disinfected or sterilized!

7.2 Installation

7.2.1 Charger

- Select the charger plug adapter which matches your power outlet (Fig. 2).

Europe

USA

Australia

UK

Fig. 2 Plug adapters for charger

- Slide the plug adapter downwards into the slots until it locks in place with a click.

WARNINGS

- *Prior to first use, the battery must be charged for 6 hours!*
- *Use the original charger only.*

7.2.2 Rechargeable Battery

RAYPEX® 6 is powered by a Nickel-Metal Hydride (NiMH) rechargeable battery. Battery status during the operation is shown on the main screen:

Battery Symbol

White: indicates battery level from fully charged to approximately 20% of capacity.

Flashing red: When battery level drops under 20%, the indicator will turn red and start flashing. The battery needs to be recharged (Fig. 3).


Fig. 3 Battery symbol is flashing red

NOTE

The device with flashing red battery symbol is still functional and may be used for several treatments before the device shuts down.

Charging Screen

Battery status during battery charging is shown on the charging screen:

Red: charging, low battery level

Yellow: charging, medium battery level

Green: charging completed, battery fully charged

To charge the battery, follow the next steps:

- Make sure that the measurement cable is not connected to the patient.
- Remove the measurement cable from RAYPEX® 6.
- Connect the external charger to the device jack and plug it into the power outlet.
- If the battery is fully discharged, e.g. when it has not been used for a long period of time, it must be recharged for 6 hours.
- If the battery is low, 4 hours recharging is sufficient.

WARNINGS

For removing the charger, follow the sequence below:

- Before removing the external charger from the power outlet, disconnect the charger from the device.
- The device cannot be used while charging.

Battery Replacement

The battery compartment is located at the rear of RAYPEX® 6 and its cover is secured by screws.

- Release the screws and remove the battery compartment cover.
- Remove the battery from battery compartment and disconnect the battery cable plug from the battery jack of RAYPEX® 6.
- Insert the cable plug of the new battery into the battery jack.
- Insert the battery into the battery compartment.
- Close the battery compartment and secure the cover with the screws.

WARNINGS

- Before replacing the battery, disconnect the charger from the device!
- The used battery must be disposed of in accordance with the local regulations

7.3 Description of User Interface

RAYPEX® 6 has a tilt adjustable front panel with a large graphic TFT touch screen.

The main screen shows the following icons and symbols:

	Settings icon
	Sound icon
	Measuring cable connector and lip clip/ file clip symbols
	Full canal image
	Apical zoom image

RAYPEX® 6 Overview


Fig. 4 RAYPEX® 6 front view


Fig. 5 RAYPEX® 6 rear view

7.4 Operation

- Important Advice for successful determination**
- Always use a rubber dam to isolate the working area
 - Always wear gloves during the measuring procedure.
 - Prior to measurement the access cavity must be dried with a cotton pellet in order to prevent incorrect measurement.
 - Select the size of the measuring file matching the root canal diameter.

7.4.1 Connecting the Device

- Press the On/Off button to switch the device on. After welcome melody and welcome screen, the main screen is displayed.
- In the framed window “plug-and-start” symbols are shown to indicate how to connect the device properly.

 indicates that the measurement cable is disconnected (Fig. 6).


Fig. 6 Measurement cable is disconnected

- Make sure that the measurement cable is not connected to the patient.
- Plug the measurement cable into the receptacle on the right side of RAYPEX® 6.

The connector symbol will turn from red to grey and the two electrode symbols turn yellow.



Fig. 7 Measurement circuit is open

To close the measurement circuit and starting the length determination respectively, see next chapter 7.4.2.

Optional cable connection test

From time to time it is recommended to check the cables.

- Connect the file clip contact to the lip clip

The connector symbol and the lip clip/file clip symbols will turn green indicating proper connection (Fig. 8).



Fig. 8 Cable connection test

NOTE

Measurement cable with attached lip clip and file clip constitute Applied Parts of the device.

7.4.2 Starting Length Determination

- Connect the file clip and the lip clip to the measurement cable.
- Before connecting the measurement cable with attached lip clip and file clip to the patient, make sure that the measurement cable is plugged into device receptacle.
- Place the lip clip on the patient's lip on the opposite side of the tooth to be treated.
- Insert the file into the root canal and clamp the file clip to the file (Attach the file clip on the metal part – directly underneath the plastic handle).

NOTE

Touch probe (included in the accessories set) for convenient measuring in the molar region. It is rather easy to handle as it is not necessary to fix the file clip to the measuring file, simply touch the metal part of the measuring file with the fork shaped end of the probe.

Two initial beeps indicate closed measurement circuit and beginning of length determination. The file movement in the canal is shown on the full canal image in the left part of the screen (Fig. 9).


Fig. 9 Beginning of length determination

NOTE

Absence of two initial beeps and no progress of the file indicate faulty connection:

- *Disconnect the measurement cable from the patient.*
- *Check proper connection between the measurement cable and the accessories.*
- *Clean the file clip/touch probe contact.*
- *Irrigate the canal, if necessary, and start again.*

WARNINGS

- *We do not recommend to continue the measurement if the initial 2 beeps are not emitted.*
- *Go to the Settings Menu and select the CHECK Mode to check the performance of the device.*

7.4.3 Apex Localisation

Coronal and Medial Section

Slowly insert the measuring file into the canal.

File movement along the coronal and medial section towards the apical region is represented on the full canal image by the ellipse continuously moving down (Fig. 10).


Fig. 10 Coronal / medial section

Apical Section

Zoomed view of the file progression is shown on the enlarged image of the apical part of the canal – the Apical Zoom (Fig. 11).


Fig. 11 Apical section - blue

In the apical section the indication line indicates the exact position and changes accordingly from blue to green and then to yellow (Fig. 12, Fig. 13).


Fig. 12 Apical section - green


Fig. 13 Apical section - yellow

File movement in the Apical Zoom is accompanied by audio signals, which serve as additional indication of the file tip position. The interval between the beeps becomes shorter the more the file approaches the apex.

When the file tip reaches the apical foramen the indication line is marked red and a constant sound is emitted (Fig. 14).


Fig. 14 Apical foramen – red bar

NOTE

The apical indication line shows the file tip position inside the canal:

- blue section: warning zone very close to the apical region
- green to yellow sections: apical region
- red bar: apical foramen

Over instrumentation

Once the file tip has passed the apical foramen, the warning red dot appears beneath the apical zoom image and short warning beeps will sound (Fig. 15).


Fig. 15 Apical foramen is passed - red dot

⚠ WARNING

As with all electronic length determination devices, the bars shown in the apical zoom do not represent the distance in millimeters.

Interruption of measurement

During the length determination the file clip may be disconnected from the file and reconnected at any time (e.g. when the file is changed to a larger size or when the length of another canal should be determined). The device detects automatically that a new length determination cycle is initiated and indicates it with two short beeps.

Completion of measurement

- Before unplugging the measurement cable from the device receptacle, disconnect the lip clip and the file clip from the patient.
- Move the file stopper to the selected reference point on the tooth.
- Gently remove the file from the canal and measure the apical length between the stopper and the file tip.

7.4.4 Sound Level Selection

To adjust the sound volume of RAYPEX®6, tap the sound icon on the screen (Fig. 16).


Fig. 16 Tap the sound icon

Adjust the preferred sound volume (Fig.17).

The sound volume refers to the warning beeps and to the ON/OFF button melody.


Fig. 17

Select the sound volume

NOTE

When RAYPEX®6 is turned off, the selected sound volume is stored in the device memory and is activated automatically when the device is switched on.

7.4.5 Automatic Shutdown

RAYPEX®6 automatically turns off after 5 minutes when not in use.

NOTE

To prolong the lifetime of the battery, it is recommended to switch the device off after each measurement.

7.5 Tips for Successful Length Determination

en

Please review the below checklist to better understand any implausible measuring results and to take appropriate actions:

Too fast movement or even jumping to the apex is indicated for the following reasons:

Symptom	Solution
Excess liquid in the pulp chamber or root canal (rinsing solution, blood or saliva), creating wrong conductive path and incorrect measurements	Dry the access cavity with a cotton pellet/air-blower. Wait until excess bleeding can be stopped.
Gingival proliferation can lead to direct contact with the measuring file causing a short circuit and incorrect measurements	Isolate the access cavity by: <ul style="list-style-type: none"> • adequate preparation filling • placing a rubber dam • electrocauterizing
The measuring file contacting metallic restorations (crown, parapulpal post, amalgam filling) may cause a short circuit and incorrect measurements	Carefully enlarge the access cavity and isolate with flow composite. Widen carefully the opening at the top of the crown

Too lazy or extremely delayed movement is indicated for the following reasons:

Symptom	Solution
Obliterated root canal Impeding the conductive path and preventing normal device functioning.	<ul style="list-style-type: none"> • Check the comparative x-ray for hints. • Catheterize with ISO 06/ 08 file till the working length.
Retreatment: Blockage by old canal filling material residue, impeding the conductive path and preventing normal device functioning	<ul style="list-style-type: none"> • Take an x-ray to re-check and try to completely remove the old canal filling material prior to measuring
Blockage by remnants of a medicated substance (e.g. calcium hydroxide) impeding the conductive path and preventing normal device functioning	Completely remove the remnants prior to measuring.
Extremely dry root canal impeding the conductive path and preventing normal device functioning	Rinse root canal with irrigating solution such as NaCl or NaOCl and dry the access cavity with a cotton pellet/air-blower.


WARNING

In some cases precise determination of file position cannot be obtained.

Special condition symptom: Exceptionally large apical foramen due to lesion or incomplete formation	May lead to shorter measurement than the actual length.
Root fracture or perforation	May lead to incorrect measurements.

Comparison of electronic length determination versus radiography:

X-ray radiography represents two-dimensional projection of three-dimensional root canal system. There are some cases when the radiographic length and the electronic length do not match.

In the case of a lateral canal curvature the x-ray may show a shorter working length than with RAYPEX® 6.

Electronic length determined using RAYPEX® 6 is usually more accurate than the radiographic length.

7.6 Device Setup

To enter the Settings Menu, tap the settings icon (Fig. 18).

Fig. 18 Tap the settings icon

To select the settings features, use the scrolling bar (Fig. 19). (see 6 features in chapter 7.6.1)


Fig. 19 Scroll between settings

To exit the Settings Menu, tap the back arrow icon or the RAYPEX® 6 “home” icon (Fig. 20)





Fig. 20 Exit settings menu

7.6.1 Settings Menu Features

	DR'S CHOICE: set optional Dr's Choice apical line in the Apical Zoom
	CHECK mode: check the operation of both apex locator/cables
	DEMO mode: activate for demo purpose
	Display brightness: adjust the brightness
	Display background: select between dark or inverse white colour
	Sound type: select between 2 different sounds

7.6.2 Optional DR'S CHOICE Apical Line

This feature enables to mark an individual pre-determined reference position at the required distance from the apex. This variable apical line can be set between the first green bar and the last yellow bar.

When DR'S CHOICE apical line is set, clear visual and audio indication is given that the file tip has reached this pre-selected position.

To set the DR'S CHOICE apical line or to modify the apical line position, follow the next steps:

- Enter the Settings Menu and select DR'S CHOICE (Fig. 21).


Fig. 21 Select DR'S CHOICE

- Place your finger on the blue dot and drag it to the preferred position in the apical section (Fig. 22).


Fig. 22 Slide up the blue dot

- Having set this variable apical line at the preferred position, e.g. at the last green bar (Fig. 23).


Fig. 23 DR'S CHOICE is set

- To cancel the DR'S CHOICE apical line return to the Settings Menu, select DR'S CHOICE and tap OFF (Fig. 25).


Fig. 25 Tap OFF

If DR'S CHOICE is set, the DR'S CHOICE apical line is additionally shown on the main screen during measurements (Fig. 24).


Fig. 24 Main screen with DR'S CHOICE set

When the file tip reaches the DR'S CHOICE apical line and during further advance of the file special beeps sound, clearly distinguished from the regular beeps.

When the apical foramen is reached, a solid tone sounds as usual. If overinstrumentation occurs, an audio warning signal sounds.

- To exit DR'S CHOICE, tap the back arrow icon or the RAYPEX® 6 "home" icon (Fig. 26).


Fig. 26 Exit DR'S CHOICE

7.6.3 CHECK Mode

This built-in CHECK feature enables to test automatically the basic function of the device in a first step, followed by the test of its accessories in a second step using the dedicated tester (Fig. 27).

Fig. 27 Tester

To profit from the CHECK feature, follow the next steps:

- Make sure that the measurement cable is not connected to the patient.
- Disconnect the measuring cable/charger from the device.
- Enter the Settings Menu and select the CHECK feature (Fig. 28).

Fig. 28 Select CHECK

Device Functional CHECK

- Plug the tester in the device as shown on the screen (Fig. 29).

Fig. 29 Insert the tester

- Device CHECK will start automatically and the results will be shown on the screen
- either OK (Fig. 30) to indicate that the device fully operates or
- ERROR (Fig. 31) to indicate any malfunction.

Fig. 30 Device OK message

⚠ WARNING

ERROR message indicates that the device is not functioning properly.

Fig. 31 Device ERROR message

For assistance call your local distributor or VDW GmbH directly.

- Disconnect the tester from the device and prepare the test of the cable with accessories.

Cables Functional CHECK

🔑 NOTE

If the device functional CHECK is OK, then you have to proceed to the cables functional CHECK.

- Connect the measuring cable to the device (Fig. 32).


Fig. 32 Connect the measuring cable

- Insert the file clip and lip clip (or use instead of the lip clip a second file clip) in the measuring cable.
- Connect file clip and lip clip (or a second file clip) to the contact strips on the tester as shown on the screen (Fig. 33).


Fig. 33 Connect file clip and lip clip to tester

- Cables CHECK will start automatically and the results – either OK (Fig. 34) or ERROR (Fig. 35) - will be shown on the screen.

Fig. 34 Accessories are OK

! WARNING

ERROR message indicates that the accessories are not functioning properly (cable breakage) or the contact area is dirty.


Fig. 35 Accessories ERROR message

For assistance call your local distributor or VDW GmbH directly.

- To exit the CHECK mode, tap the back arrow icon.
- Disconnect the measuring cable from the device.

7.6.4 DEMO Mode

The DEMO mode allows you to get acquainted with the device and to demonstrate its operation to your patients.

To activate the DEMO mode, follow the next steps:

- Disconnect the measuring cable or charger (if connected) from the device.

NOTE

The DEMO mode will not start if the measuring cable is connected.

- Enter the Settings Menu and select the DEMO feature (Fig. 36).


Fig. 36 Select DEMO

- During the DEMO sequence “DEMO” is written on the screen, indicating a measuring simulation.
- Tap the screen to pause/resume the demo simulation.

- To exit the DEMO mode, tap the RAYPEX® 6 “home” icon (Fig. 37).


Fig. 37 Exit DEMO mode

NOTE

If the measuring cable is connected during the DEMO mode running, RAYPEX® 6 will automatically exit the DEMO mode.

NOTE

Servicing and repair work should only be carried out by service engineers trained at the factory.

The measuring cable and the surface of the equipment should be cleaned with a paper towel or soft cloth or wipe that has been wetted with aldehyde-free disinfection and cleaning solution (bactericidal and fungicidal), e.g

- DENTIRO® wipes, Oro Clean Chemie AG, Switzerland
- TopActiv wipes, ad-Arztbedarf GmbH, Germany
- SprayActiv, ad-Arztbedarf GmbH, Germany
- DY wipes, DENTSPLY, France

! WARNINGS

- *To disinfect the equipment and the measuring cable, wipe the surface with a clean cloth or wipe that has been wetted only slightly with a non-aggressive disinfectant.*
- *Do not use any liquids or sprays directly on the equipment, particularly the display.*
- *Do not use strong alcohol to disinfect surfaces.*

8 Maintenance, cleaning and sterilization

8.1 General information

The RAYPEX® 6 is maintenance-free and does not have any components that need to be maintained by the user.

Before preparing them for reuse, remove the lip clip, file clip or touch probe from the measuring cable.

The lip clip, file clip and touch probe must always be cleaned, disinfected and sterilized before they are used. This also applies before they are used for the first time. The only way to ensure that these accessories are effectively sterilized is to thoroughly clean and disinfect them. When doing this, follow the special instructions set out in section 8.2, "Cleaning, disinfection and sterilization".

In addition, when you use any of the other equipment that you have in your practice, make sure you follow the directions for use for this equipment. As part of your responsibility for ensuring that the accessories are sterile, always make sure that only validated cleaning, disinfection and sterilization methods are used, that the equipment (sterilizer) is regularly maintained and inspected, and that the validated parameters are adhered to each time.

In addition, ensure that you always comply with the applicable statutory rules and regulations with regard to hygiene in your practice or clinic. In particular, this applies to the guidelines regarding effective prion inactivation.

For your own safety, when handling contaminated accessories, always wear gloves, protective eyewear and a surgical mask.

WARNINGS

- *The measuring cable cannot be autoclaved.*
- *The use of substances and materials other than those referred to above may damage the equipment and accessories.*
- *Do not employ sterilization methods that use heat, radiation, formaldehyde, ethylene oxide or plasma.*

8.2 Cleaning, Disinfection and Sterilisation

(in accordance with DIN EN ISO 17664)

The cleaning, disinfection and sterilization process applies only to the lip clip, file clip and touch probe. As long as they are treated with due care and are not damaged or contaminated, these accessories can be reused multiple times.

We accept no liability in the event that these instructions are disregarded or processes that have not been validated are used to prepare the accessories for reuse.

Make sure undamaged sterile packaging is used.

8.2.1 Pre-treatment

- Pulp and dentine remnants must be removed immediately from the accessories; if necessary, use a brush to carefully clean the products (no later than 2 hours after use). Do not allow remnants to dry on. After using the accessories on the patient, place them directly into a dish filled with a suitable cleaning and disinfection solution (e.g. CIDEZYME®, ENZOL® Enzymatic Detergent Solutions, Johnson & Johnson Medical, 0.8% for between 1 minute and 2 hours) for cleaning, pre-disinfection and interim storage. Make sure that the products are fully immersed. Then wash the accessories under running sterile, deionized water or in a disinfection solution at least three times for one minute each time, in order to remove all visible traces of contamination and remnants. The disinfectant should be aldehyde-free (aldehyde causes blood to stain permanently), tested for efficacy (e.g. VAH/DGHM or FDA certification or CE mark), suitable for disinfecting the accessories, and compatible with the accessories (see section 8.2.7, "Material resistance").

- Only use clean, soft brushes to manually remove contamination and remnants, or a clean, soft cloth or wipe that is only used for this purpose. Do not use metal brushes or wire wool. Check that no visible contamination or remnants remain, and repeat the pre-cleaning process if necessary.

- The file clip must be squeezed together and then released five times during cleaning to allow the inner parts to be cleaned more effectively. Please note that disinfectant used for pre-treatment is only for personal protection and does not obviate the need for disinfection once cleaning has been completed. Pre-treatment should never be omitted.

WARNINGS

Do not use automated processes for cleaning or disinfecting the accessories.

8.2.2 Manual cleaning and disinfection

When choosing cleaning agents and disinfectants, make sure that

- They are suitable for cleaning or disinfecting instruments
- You use a disinfectant with proven efficacy (e.g. with VAH/DGHM or FDA certification or CE mark) and that the disinfectant is compatible with the cleaning agent
- The chemicals used are compatible with the accessories (see section 8.2.7, "Material resistance").

Combined cleaning/disinfection products should only be used if the instruments are only slightly contaminated (no visible contamination/remnants).

Comply with the concentrations and contact times specified by the manufacturers of the cleaning agents and disinfectants, as well as their instructions regarding the intensity of subsequent rinsing.

Use only freshly prepared solutions, water that is sterile or has a low microbe content (< 10 cfu/ml) and a low endotoxin content (< 0.25 EU/ml, e.g. purified water (PW/HPW)), and filtered, oil-free air for drying.

Make sure that the accessories are not in direct contact with one another.

Step-by-step procedure

Cleaning

- Place the pre-cleaned accessories into the cleaning bath for the prescribed contact time (e.g. CIDEZYME®, ENZOL® Enzymatic Detergent Solutions, Johnson & Johnson Medical, 0.8% for 1 minute); make sure that the products are fully immersed (if necessary, use a soft brush to carefully brush them down). The file clip must be squeezed together and then released five times during cleaning to allow the inner parts to be cleaned more effectively.

- Then, remove the accessories from the cleaning bath and rinse them off thoroughly at least three times for one minute each time with sterile, deionized water; when you rinse it, squeeze the file clip together and then release it five times.

- Next, place the accessories in an ultrasonic bath with a cleaning agent (e.g. CIDEZYME®, ENZOL® Enzymatic Detergent Solutions, Johnson & Johnson Medical, 0.8% for 20 minutes); make sure that the products are fully immersed (if necessary, use a soft brush to carefully brush them down). The file clip must be squeezed together and then released five times during cleaning to allow the inner parts to be cleaned more effectively.

- Then, remove the accessories from the ultrasonic bath and rinse them off thoroughly at least three times for one minute each time with sterile, deionized water; when you rinse it, squeeze the file clip together and then release it five times.

Disinfection

- Once the accessories have been cleaned and inspected, place them into the disinfection bath for the prescribed contact time (e.g. Cidex OPA, Johnson & Johnson Medical, 100% for 20 minutes); the accessories must be sufficiently immersed in the solution. The file clip must be squeezed together and then released five times during disinfection to allow the inner parts to be disinfected more effectively.

- Then, remove the accessories from the disinfection bath and rinse them off thoroughly at least five times for one minute each time with sterile, deionized water; when you rinse it, squeeze the file clip together and then release it five times.
- Dry the accessories by blowing them down fully with oil-free, filtered compressed air and then leaving them to dry further in a clean place for at least 20 minutes.

Once the accessories are dry, inspect and pack them as soon as possible (see the section entitled "Inspection/maintenance" and "Packaging").

8.2.3 Inspection/Maintenance

Inspect all accessories once they have been cleaned or cleaned/disinfected. Dispose of any defective accessories immediately. These defects include:

- Deformation of the plastic
- Corrosion
- Discoloration of the plastic

Accessories that are still contaminated must be cleaned and disinfected again. Maintenance is not required. Do not use instrument lubricant.

8.2.4 Packaging

Pack accessories into single-use sterilization packaging (single-layer packaging) that fulfils the following requirements:

- Complies with DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (withstands temperatures of up to 142 °C (288 °F) or more, sufficient vapour permeability)

8.2.5 Sterilization

Use only the sterilization methods listed below; other sterilization methods are not permitted.

- Steam sterilization
- Fractionated vacuum/pre-vacuum method (at least three vacuum cycles) or gravity method with sufficient product drying. The drying time that is actually required depends directly on parameters that are the sole responsibility of the user (loading configuration, how many items are loaded and how closely together they are loaded, condition of the sterilizer, etc.) and must therefore be established by the user. However, the drying time must never be less than 20 minutes.

- The gravity method is less effective and must only be used if it is not possible to use the fractionated vacuum method.

- Steam sterilizer in accordance with DIN EN 13060 or DIN EN 285, ANSI/AAMI ST 79

- The sterilization process must be validated in accordance with DIN EN ISO 17665 (valid installation and operational qualification (IQ and OQ) and product-specific performance qualification (PQ))

- Maximum sterilization temperature 135 °C (275 °F) plus tolerance in accordance with ISO DIN EN ISO 17665

- Sterilization time (exposure time at sterilization temperature) at least 3 minutes at 134 °C (273 °F)

Rapid sterilization and sterilization with unpackaged accessories are not permitted.

Dry heat sterilization, radiation sterilization and sterilization using formaldehyde, ethylene oxide or plasma are also not permitted.

WARNINGS

- *Regional requirements have to be followed additionally.*

8.2.6 Storage

Once they have been sterilized, the accessories must be kept dry and free of dust in their sterilization packaging.

- Sterility cannot be guaranteed if packaging is open or damaged (check the packaging before using the accessories).

8.2.7 Material resistance

When choosing cleaning agents and disinfectants, make sure that they do not contain any of the following substances:

- Phenol
- Strong acids (pH < 6) or strong alkalis (pH > 8)
- Aldehydes
- Anti-corrosive substances (especially di- or triethanolamine)
- Oxidants (hydrogen peroxide, sodium hypochlorite over 5% strength)
- Solvents

The material can withstand temperatures of up to 137 °C/279 °F (maximum exposure temperature).

9 Conformity

The device is CE-marked according to MDD 93/42/EEC, revised by the Directive 2007/47/EC.

For details see the Declaration of Conformity under www.vdw-dental.com

Notified Body:
TÜV Produkt Service GmbH
Ridlerstraße 65
80339 Munich, Germany

10 Manufacturer

VDW GmbH
Bayerwaldstr. 15
81737 Munich
Germany

Phone: +49- (0) 89-6 27 34-0
Fax: +49- (0) 89-6 27 34-190
Service: +49- (0) 89-6 27 34-555

Website: www.vdw-dental.com
Email: info@vdw-dental.com

11 Warranty

In addition to the warranty applicable under the sales contract with the dental dealer, VDW directly offers customers the following service warranty:

1. VDW guarantees proper manufacturing of the product, the use of top quality materials, performance of all required tests and adherence to all applicable pertinent laws and regulations relating to the product.

The full functionality of RAYPEX®6 apex locator is covered by a 12 month warranty, beginning on the date of delivery to the customer (according to the shipping papers with the respective product serial number, issued by the seller at the time of purchase).

Cables and battery are covered by a warranty period of 6 months.

The customer is only entitled to warranty services within the warranty period and under the condition that VDW has been informed of the defect in writing within a two month period from the date of discovery of the defect.

2. In the event of a warranty claim, the VDW Service Center in Munich, Germany will handle any repairs/replacement within 3 business days from receipt at the VDW facility in Munich plus transportation time required to return the goods to the customer.

3. This warranty covers only the exchange or repair of individual components or parts affected by manufacturing defects. The cost for personnel to be dispatched for purposes of technical assistance from the dealer to the customer, and/or packaging expenses of the customer will not be covered by VDW.

Any customer claims beyond the realm of repair, such as claims for damages, will not be covered.

This warranty does not include any compensation for direct or indirect personal injuries or material damage of any kind.

The customer will not be entitled to damages covering for downtime of the equipment.

4. The warranty does not extend to damage for which evidence can be provided by VDW that such damage has been caused by negligence on the part of the user as related to Directions for Use, particularly during battery charging or battery replacement.

The warranty explicitly excludes any defects arising on the basis of:

- damage occurring during transportation to VDW for the purposes of repair,
- damage occurring through environmental events such as lightning strikes, fire and/or humidity.

This warranty shall automatically become cancelled and void if the product has been improperly repaired, modified or manipulated in any manner by the user, by non-authorized persons or by third party personnel.

5. This warranty is valid only if the device sent for repair includes the invoice with confirmation of the product's shipping date.

6. Legal claims, such as under product liability law or claims against suppliers from whom the customer has purchased the product, in particular the dental dealer, remain unaffected.

12 Disclaimer

VDW, its representatives and dealers shall have no liability or responsibility for damages caused by the clinical use of their products, whether the use is accidentally associated with other electro-medical devices (e.g. pacemakers) or not.

VDW, its representatives and dealers shall bear no liability or responsibility for damages caused by improper use of the device, not according to the Directions For Use.

13 Technical Data

RAYPEX® 6 electronic apex locator is a programmable electrical medical device.

RAYPEX® 6 belongs to the following category of medical devices:

INTERNALLY POWERED DEVICE	NiMH rechargeable battery
PROTECTION AGAINST ELECTRIC SHOCK	Type BF
LEVEL OF SAFETY IN PRESENCE OF INFLAMMABLE ANAESTHETIC MIXTURES OR OXYGEN	Not suitable for use in proximity of inflammable anaesthetic mixtures containing air, oxygen or nitrogen oxide
OPERATING MODE	Continuous
PROTECTION AGAINST LIQUID PENETRATION	No protection against ingress of liquids
TRANSPORTATION AND STORAGE CONDITIONS	Temperature range: -20 °C to +60 °C (4 °F to 140 °F) Relative humidity: 10 % to 90 %, non-condensing Atmospheric pressure: 50 kPa to 106 kPa
OPERATING CONDITIONS	Temperature range: +10 °C to + 40 °C (50 °F to 104 °F) Relative humidity: 10 % to 90 %, non-condensing Atmospheric pressure: 70 kPa to 106 kPa
DIMENSIONS	Folded: Width: 101 mm, Length: 110 mm, Thickness: 27 mm Open: Width: 101 mm, Length: 110 mm, Height: 97 mm
WEIGHT	350 gr.
DISPLAY TYPE	3.5" TFT colour display with touch panel
DISPLAY / ACTIVE AREA	70 mm x 53 mm
EXTERNAL CHARGER	Input: AC 100-240V, 50/60 Hz
EXTERNAL CHARGER	Output: 5V DC, 1000 mA

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en

NOTES

- The RAYPEX® 6 requires special precautions with regard to electromagnetic compatibility.
- It must be installed and prepared for use as described in section 7.4 "Operation".
- Certain types of RF wireless communication equipment such as mobile telephones are likely to interfere with the RAYPEX® 6.
- The recommended radiation levels of RF wireless communication equipment specified in this paragraph must therefore be complied with.
- The RAYPEX® 6 must not be used near or on top of another device. If this cannot be avoided, it is necessary – before clinical use – to check the equipment for correct operation under the conditions of use.
- The use of accessories other than those specified or sold by VDW as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the unit.

Electromagnetic Emissions

NOTES

- RAYPEX® 6 is intended for use in the professional healthcare facility or home healthcare electromagnetic environment specified in the tables below.
- The user and/or installer of the unit must ensure that it is used in such an environment.

Guidance and Manufacturer's Declaration – Electromagnetic Emissions – RAYPEX® 6

The RAYPEX® 6 is intended for use in professional healthcare facility or home healthcare electromagnetic environment specified below;
The customer or the user of the RAYPEX® 6 should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The RAYPEX® 6 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The RAYPEX® 6 is suitable for use in professional healthcare facilities or domestic establishments.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity – RAYPEX® 6			
The RAYPEX® 6 is intended for use in the electromagnetic environment specified below; The customer or the user of the RAYPEX® 6 should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact	± 8 kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
	± 15 kV air	± 15 kV air	
Electrical fast transients/bursts, IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines	Mains power quality should be that of a typical public low-voltage power supply network that supplies buildings used for domestic purposes, commercial or hospital, clinic environment.
	± 1 kV for input/output lines	Not Applicable	
Surges IEC 61000-4-5	± 1 kV Line-to-line	± 1 kV Line-to-line	Mains power quality should be that of a typical public low-voltage power supply network that supplies buildings used for domestic purposes, commercial or hospital, clinic environment.
	± 2 kV Line-to-ground	± 2 kV Line-to-ground	
Voltage dips, Voltage interruptions IEC 61000-4-11	0 % UT; 0.5 cycle	0 % UT; 0.5 cycle	Mains power quality should be that of a typical public low-voltage power supply network that supplies buildings used for domestic purposes, commercial or hospital, clinical environment. If the user of the RAYPEX® 6 requires battery charging during power mains interruptions, it is recommended that RAYPEX® 6 charger be powered from a separate power supply (UPS, etc.).
	0 % UT; 1 cycle and 70 % UT; 25/30 cycles	0 % UT; 1 cycle and 70 % UT; 25/30 cycles	
Rated power frequency magnetic fields IEC 61000-4-8	0 % UT; 250/300 cycles	0 % UT; 250/300 cycles	Power frequency magnetic fields should be at levels characteristic of a typical public low-voltage power supply network that supplies buildings used for domestic purposes, commercial or hospital, clinic environment.
	30 A/m 50 or 60 Hz	30 A/m 50 or 60 Hz	
<p> NOTE <i>UT is the a.c. mains voltage prior to application of the test level.</i></p>			

Guidance and Manufacturer's Declaration – Electromagnetic Immunity – RAYPEX® 6

The RAYPEX® 6 is intended for use in the electromagnetic environment specified below;
The customer or the user of the RAYPEX® 6 should assure that it is used in such an environment.

Immunity test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic environment - guidance
Conducted disturbances induced by RF fields	3Vrms 150 kHz to 80 MHz	3Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the RAYPEX® 6, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
IEC 61000-4-6	6 Vrms in ISM bands 150 kHz to 80 MHz 80% AM at 1 kHz	6 Vrms in ISM bands 150 kHz to 80 MHz 80% AM at 1 kHz	Recommended separation distance $d = 1.17\sqrt{P}$ $d = 1.17\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.7 GHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manu- facturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b

Interference may occur in the
vicinity of equipment marked with the
following symbol:

NOTES

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures objects and people.

- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the RAYPEX® 6 is used exceeds the applicable RF compliance level above, the RAYPEX® 6 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the RAYPEX® 6.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Specifications for enclosure port immunity to RF wireless communications equipment

The RAYPEX® 6 is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

The user and/or installer of the unit can help prevent electromagnetic interference by maintaining radiation levels of RF wireless communications equipment (emitters) within the compliance limits specified in the table below.

Recommended radiation levels of RF wireless Communications Equipment			
Frequency band	EC 60601-1-2 Test level	Compliance level	Minimum separation distance
380 - 390 MHz	27 V/m	27 V/m	0.3 m
430 - 470 MHz	28 V/m	28 V/m	0.3 m
704 - 787 MHz	9 V/m	9 V/m	0.3 m
800 - 960 MHz	28 V/m	28 V/m	0.3 m
1700 - 1990 MHz	28 V/m	28 V/m	0.3 m
2400 - 2570 MHz	28 V/m	28 V/m	0.3 m
5100 - 5800 MHz	9 V/m	9 V/m	0.3 m

NOTE
These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.