VOCO

Futurabond® U bottle

EN Instructions for use
MD EU Medical device

Product description:

Futurabond U bottle is a dental adhesive for use with methacrylate-based restorative, core build-up or luting materials. It can be applied with phosphoric acid etching (selective enamel etching or total-etch) of the dental hard tissue as well as without (self-etch).

Indications:

- [A1] Adhesive for direct restorations of all classes of cavities using light-curing composite / compomer / ORMOCER® materials with a methacrylate base.
- [A2] Adhesive for core build-ups made of light-curing composite materials.
- [A3] Adhesive for indirect restorations with light-curing luting composites.
- [B] Intraoral repair of composite, compomer or ORMOCER® restorations, ceramic veneers and all-ceramic restorations without an additional primer.
- [C] Treatment of hypersensitive tooth necks.
- [D] Protective varnish for glass ionomer cement restorations.
- [E1] Sealing of cavities prior to amalgam restorations.
- [E2] Sealing of cavities and core preparations prior to the temporary luting of indirect restorations.

In combination with **Futurabond U bottle DCA (Dual Cure Activator)**:

- [F1] Adhesive for direct self-curing or dual-curing composite restorations of all classes of cavities and for core build-ups.
- [F2] Adhesive for indirect restorations using dual-curing or self-curing luting composites.
- [G] Luting of root posts with dual-curing or self-curing luting composites.

Contraindications:

Futurabond U bottle contains organic acids, (di)methacrylates and BHT. **DCA** contains amines. **Futurabond U bottle / DCA** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Futurabond U bottle is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Futurabond U bottle should only be used by a professionally trained dental practitioner.

Use:

Futurabond U bottle not only ensures a secure bond to natural dental hard tissue, but also to composites, metals and ceramics. It is necessary to isolate the operating site from moisture. Any contamination of the cavity or material surface with blood or saliva must be prevented. The use of a rubber dam is recommended.

The subsequent stages of application depend on the indication (see sections [A] - [G]).

When using light-curing materials, **Futurabond U bottle** can be applied without the addition of **Futurabond U bottle Dual Cure Activator** (see sections [A] - [E]). The **Futurabond U bottle Dual Cure Activator** needs to be added only when working with self-curing or dual-curing materials (see sections [F] - [G]).

[A] Direct or indirect restorations with light-curing restorative, core build-up or luting composites

1. Prepare the cavity and core according to the principles of the adhesive restoration technique, and remove any plaque or tartar in the case of non-carious class V cavities. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
2. Etching of the dental hard tissue (optional)
In self-etch mode, **Futurabond U bottle** achieves high levels of adhesive strength on enamel and dentine. Prior etching of the dental hard tissue may optionally be carried out. Only enamel which is **unprepared** must always be etched at the beginning of the bonding process.
 - Etch the surfaces to be etched using a dental acid-etch agent based on phosphoric acid (see applicable instructions for use).
 - Aspirate the acid-etch agent, rinse with water for approx. 20 s. Dry off excess moisture with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance.
3. Place one drop of **Futurabond U bottle** on a mixing palette.
4. Apply the adhesive evenly to the surfaces of the cavity and **rub it in** for 20 s with a disposable applicator.
5. Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents.
6. Cure the adhesive layer for 10 s using a commercially available polymerisation device (LED or halogen light with an output of > 500 mW/cm²).
7. Apply the restorative, core build-up or luting material in accordance with the respective instructions for use.

[B] Intraoral repair of fillings or restorations

- (made from composite, glass ceramic, zirconium oxide ceramic, aluminium oxide ceramic, metal (incl. amalgam))
1. Thoroughly roughen the surface of the restoration to which the adhesive is to be applied and sandblast it, preferably with aluminium oxide < 50 µm. In the case of ceramics, remove the varnish 1 mm beyond the margin.
 2. Clean thoroughly with water spray and dry with moisture- and oil-free air.
 3. Continue with the application of **Futurabond U bottle** immediately afterwards (see [A] 3. - 6.). The use of an additional primer is not necessary.
 4. Carry out the repair using a light-curing composite according to the manufacturer's instructions for use.

[C] Treatment of hypersensitive tooth necks

1. Clean the tooth neck in the usual way (e.g. with polishing paste or using a powder jet).
2. Clean thoroughly with water spray and dry with moisture- and oil-free air. Do not overdry. The objective is a surface with a silky matt appearance.
3. Afterwards, continue with the application of **Futurabond U bottle** (see [A] 3. - 6.).
4. The oxygen inhibition layer should be removed using a cloth or a foam pellet.

[D] Protective varnish for glass ionomer cement restorations

1. If a matrix was used to facilitate shaping, remove the matrix first.
2. Place one drop of **Futurabond U bottle** on a mixing palette.
3. Apply **Futurabond U bottle** to the modelled filling with a disposable applicator.
4. Dry off **Futurabond U bottle** with dry, oil-free air for at least 5 s in order to remove any solvents.

5. Cure **Futurabond U bottle** on every surface of the restoration using a commercially available polymerisation device for 10 s (LED/halogen light with an output of > 500 mW/cm²).
6. Once the restoration has completely cured, grind the occlusion and seal the ground areas once more using **Futurabond U bottle** in order to prevent the restoration from drying out.
7. The oxygen inhibition layer should be removed using a cloth or a foam pellet.

[E] Sealing of cavities prior to amalgam restorations or prior to temporary luting

Note: Do not seal cavities with **Futurabond U bottle** if temporary luting using a methacrylate-based material is scheduled afterwards.

1. Prepare the cavity according to the further treatment planned. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
2. Place one drop of **Futurabond U bottle** on a mixing palette.
3. Apply **Futurabond U bottle** to the cavity surface to be treated with a disposable applicator and **rub it in** for 20 s.
4. Dry off **Futurabond U bottle** with dry, oil-free air for at least 5 s in order to remove any solvents.
5. Cure **Futurabond U bottle** using a commercially available polymerisation device for 10 s (LED or halogen light with an output of > 500 mW/cm²).
6. Depending on the indication, continue to work with the desired material afterwards in accordance with the corresponding instructions for use.

[F] Direct or indirect restorations with self- or dual-curing restorative, core build-up or luting composites

In combination with **Futurabond U bottle DCA**, **Futurabond U bottle** creates a strong bond to dual-curing and self-curing methacrylate-based composites and luting materials.

1. Prepare the cavity or core according to the principles of the adhesive restoration technique. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
2. Etching of the dental hard tissue (optional)
In the self-etch mode, **Futurabond U bottle** achieves high levels of adhesive strength on enamel and dentine. Prior etching of the dental hard tissue may optionally be carried out. Only enamel which is **unprepared** must always be etched at the beginning of the bonding process (e.g. in the luting of Maryland bridges).
 - Etch the surfaces to be etched using a dental acid-etch agent based on phosphoric acid (see applicable instructions for use).
 - Aspirate the acid-etch agent, rinse with water for approx. 20 s. Dry off excess moisture with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance.
3. Mix one drop of **Futurabond U bottle** and one drop of **Futurabond U bottle DCA** in a mixing palette with a disposable applicator for 3 s.
4. Apply the mixed adhesive evenly on the surfaces of the cavity and **rub it in** for 20 s with a disposable applicator.
5. Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents.
6. Cure the adhesive layer for 10 s using a commercially available polymerisation device (LED or halogen light with an output of > 500 mW/cm²).
7. Depending on the indication, continue to work with the desired material in accordance with the corresponding instructions for use.

[G] Luting of root posts

1. Thoroughly mix one drop of **Futurabond U bottle** and one drop of **Futurabond U bottle DCA** in a mixing palette with a suitable disposable applicator (e.g. **Endo Tim**) for 3 s.
2. Use the applicator to apply the mixed adhesive into the prepared root canal and to the occlusal surface of the core and **rub it in** for 20 s.
3. Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents. Remove any excess material by dabbing with an absorbent paper strip. Do not light-cure the adhesive layer.
4. Introduce the root post, which has been prepared in accordance with the manufacturer's instructions, into the root canal using a suitable self- or dualcuring composite for post luting (see applicable instructions for use) in such a way that a slight excess of material is achieved.
5. To fix the root post and cure the occlusally-applied **Futurabond U bottle**, the material should be light-cured for 40 s.

Warnings, precautionary measures:

- Phenolic substances, especially preparations containing eugenol or thymol, inhibit the curing of **Futurabond U bottle**. Therefore, do not use zinc oxide eugenol cements or other eugenol-containing substances in combination with **Futurabond U bottle**.
- Avoid contact with the oral mucosa. Contact with **Futurabond U bottle** may cause whitening of the oral mucosa as a result of the coagulation of proteins. This is a reversible effect which subsides within a few days.
- In case of eye contact, rinse well with plenty of water and consult an ophthalmologist.
- **Futurabond U bottle** is light-curing. Hence, please avoid exposure to strong lighting during application.
- The mixed adhesive is self-curing and light-curing, further processing must therefore take place immediately after mixing.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Futurabond U bottle: BisGMA, HEMA, ethanol, water, HEDMA, methacrylate phosphoric acid ester, methacrylate functionalised polyacrylate, UDMA, initiators, stabilisers
Futurabond U bottle and DCA: Ethanol, BisGMA, HEMA, water, HEDMA, methacrylate phosphoric acid ester, methacrylate functionalised polyalkenoic acid, UDMA, initiators, stabilisers

Storage:

Tightly close **Futurabond U bottle** and **Futurabond U bottle DCA** immediately after use.
Storage at 4°C-23°C. Do not expose to direct sunlight or the operating light. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Futurabond U bottle / Futurabond U bottle DCA** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Futurabond U bottle / Futurabond U bottle DCA** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 AA1577 E1 0724 99 © by VOCCO