Mojo™

Veneer Cement

English

INSTRUCTIONS FOR USE:

Mojo Veneer Cement is a light cure, esthetic cementation system that is designed for use with porcelain, ceramic, and composite veneer restorations. The simple, highly versatile shade selections allow this material to be used with a wide range of veneer cases, while offering the least amount of shade shift available of all materials tested. The corresponding water soluble try-in gels allow for a perfect match of shades each time, producing beautiful, esthetic restorations. The Mojo kit includes Bond-1™ Primer / Adhesive and a silane primer for an easy to use, comprehensive veneer bonding system.

INDICATIONS:

- Bonding of porcelain, ceramic, and composite veneers.

CONTRAINDICATIONS:

- Not to be used under opaque metals or crowns and bridges where light curing of the material is not possible.
- Not to be used on patients with known allergies to methacrylate resins
- Not for use on exposed pulp

CEMENTATION

1. Prepare the tooth according to the principles of adhesive dentistry. Make sure to clean and remove surface debris thoroughly from the preparation site.
2. Apply a small amount of Mojo Veneer Cement try-in gel to the restoration. Seat it on the tooth using light pressure. Assess shade and remove the restoration. Once the desired shade has been selected, remove the veneer and rinse thoroughly with water. Remove any organic debris by cleaning with alcohol or acetone, followed by cleaning in water in an ultrasonic cleaner. Clean the tooth thoroughly with both water and pumice.
Note: Try-in gels have unlimited working times and will not set.
3. Ensure all restorations have been properly treated and etched according to the manufacturer's instructions.
4. Apply 1-2 coats of silane to the restoration(s) and allow the silane to stand for 60 seconds. Air dry.
Note: Pentron recommends using Bond-1 Primer / Adhesive (following steps 5-8). If you choose to use another dentin/ enamel bonding system,

English

INSTRUCTIONS FOR USE

Français

French

MODE D'EMPLOI

Deutsch

German

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Italiano

Italian

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Español

Spanish

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Čeština

Czech

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please follow the manufacturer's instructions for use.

5. Apply 1 coat of Bond-1 Primer / Adhesive to the inside of the prepared restoration(s) and air dry for 10 seconds using a gentle stream of air to remove solvent.
6. Apply 37% Phosphoric Acid Etching Gel* to enamel and dentin surfaces for 15 seconds. Rinse thoroughly for 15 seconds with water. It is important not to splash etch into eyes, onto soft tissues or skin.
7. Dry with clean, oil-free, air without desiccating the dentin. Do not allow etched surface to become contaminated with saliva or blood during the bonding procedure. NOTE: Never use excessive force on syringe plunger to dispense etch which may lead to over dispensing of the material. If resistance is felt replace current tip with a new tip and check flow extraorally.
8. Apply two consecutive coats of Bond-1 Primer / Adhesive to tooth within 10 seconds to ensure a shiny surface.
9. Air dry for 10 seconds using a gentle stream of air to remove solvent. Hold the air syringe 1 inch from preparation or restoration; avoid disturbing the resin surface. If surface shine disappears after drying, re-apply with Bond-1 Primer/Adhesive and re-dry.
10. Place the selected shade of Mojo Veneer Cement onto the silanated veneer. Carefully place veneer. Allow excess cement to exude.
Note: Mojo Veneer Cement shades can be mixed for customization to further "warm" or "brighten" the shade. If shades are mixed, be careful not to incorporate air into the mixture.
11. Light cure for 1-2 seconds to tack the veneers into place. Remove excess cement as needed.
12. Light cure each surface of the restoration for 20 seconds.
13. Finish with diamond/carbide burs, if necessary.

NOTE: Tips and applicator brushes are for single patient use only, to prevent cross-contamination between patients.

NOTE: Refer to outer package for expiration date. Do not use after expiration date.

***CAUTION:** 37% Phosphoric Acid Etching Gel contains phosphoric acid. Avoid contact with skin, eyes, and soft tissue. In case of contact with skin or eyes, flush immediately with water; get medical attention for eyes. Do not take internally.

CAUTION: Bond-1 Primer/ Adhesive contains acetone and is flammable and volatile. Keep away from sources of ignition and excessive heat.

CAUTION: Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

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Veneer Cement

Rx Only

Refrigeration required.
Store at 2-8°C. (36-46°F)

LIMITED WARRANTY –

Pentron's technical advice, whether verbal or in writing, is designed to assist dentists in using Pentron's product. The dentist assumes all risk and liability for damages arising out of the improper use of Pentron's product.

In the event of a defect in material or workmanship, Pentron will replace the defective product or part thereof, or reimburse the actual cost of the defective product. In order to take advantage of this limited warranty, the defective product must be returned to Pentron.

PENTRON

 Manufactured by Pentron Clinical

1717 West Collins Ave. Orange, CA 92867 USA

800-551-0283 / 714-516-7557 / www.pentron.com

 European Union Representative

SpofaDental a.s.

Markova 238

CZ-506 01 Jičín

Phone: +420 493 583 111