VOCO

Ionobond

EN **Instructions for use**
MD EU Medical device
In accordance with DIN EN ISO 9917-1

Product description:

Ionobond is a glass ionomer cement for lining and core build-up.
Ionobond sets rapidly and is radiopaque.

Indications:

For cavity lining and core build-up.

Contraindications:

Ionobond contains polyacrylic acid and parabens. **Ionobond** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Ionobond is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Ionobond should only be used by a professionally trained dental practitioner.

Use:

After preparing the cavity, clean and dry enamel and dentine thoroughly.

Shake powder before use.

Use 1 drop of liquid to 1 level measuring scoop of powder (hold bottle vertically). This corresponds to a mixing ratio of 3.6 g powder to 1 g liquid.

Only mix it in the specified mixing ratio.

Place powder and liquid onto a glass plate or suitable mixing pad and mix the powder into the liquid, a little at a time to begin with and then in larger quantities.

Use a suitable spatula for mixing.

Mixing time	Working time	Setting time
30 s	1.5 - 2 min	3 - 4 min

The mixing time is approx. 30 s. **Ionobond** can be worked for approx. 1.5 - 2 min.

Apply **Ionobond** in a thin, uniform layer to all dentine surfaces. **Ionobond** will set in the mouth within 3 - 4 min.

In case the definite treatment is not carried out immediately, fill the whole cavity with **Ionobond** and during the second treatment session, remove the upper layer to insert the definite filling material. Fill the cavity with a suitable filling material.

Warnings, precautionary measures:

- Rinse metal instruments with cold water immediately after use and before the cement residues can set.
- For very deep cavities (remaining dentine layer less than 1.5 mm) pulpal protection with an appropriate product is required.
- Certain substances (e.g. varnishes, insulating agents, eugenol residues etc.) can impair adhesion and should therefore not be used if possible, with the exception of pulp protection.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Calcium aluminium fluorosilicate glass, polyacrylic acid, water, barium sulfate, tartaric acid, fumed silica, parabens

Storage instructions and application method:

Tightly close powder and liquid container immediately after use. Store in dry conditions at 4-23°C. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Ionobond** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Ionobond** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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