VOCO

x-tra base®

EN Instructions for use
MD EU Medical device
 In accordance with DIN EN ISO 4049

Product description:

x-tra base is a radiopaque, light-curing flowable composite. In contrast to conventional flowables, it is characterised by an increased depth of cure while at the same time showing reduced shrinkage stress, which means that it can be cured in layers of 4 mm thickness.

When used as a base/liner **x-tra base** has to be overlaid occlusally with an at least 2 mm thick surface layer of a methacrylate-based universal / posterior composite. **x-tra base** is used with a dentine enamel bond. **x-tra base** contains 75 % w / w inorganic fillers in a methacrylate matrix.

Shades:
 Universal, A2

- Indications:**
- Base in class I and II cavities
 - Cavity lining under direct restorative materials in class I and II cavities
 - Small, non occlusal stress-bearing class I restorations according to minimally invasive filling therapy
 - Class III and V restorations
 - Extended fissure sealing
 - Undercut blockout
 - Repair of small enamel defects
 - Repair of small defects in aesthetic indirect restorations
 - Repair of temporary C&B-materials
 - Core build-up

Contraindications:
x-tra base contains methacrylate and BHT. **x-tra base** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:
x-tra base is suitable for use in all patients without any age or gender restrictions.

Performance features:
 The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:
x-tra base should only be used by a professionally trained dental practitioner.

Use:
Preparation:
 Clean the teeth being treated. If applicable, mark occlusal contact points. Bring material to room temperature before use.

Shade selection:
 Select a shade of the final composite chosen for the occlusal surface by comparing the shade guide with the cleaned tooth before anaesthesia. Moistened the shade sample and the tooth and select the shade, preferably in daylight.

Cavity preparation:
 In principle, the cavities should be prepared according to the rules for adhesive restorative treatment; preparation should be minimally invasive in order to protect healthy dental hard tissue. Then clean and dry the cavity.

Creating a dry working field:
 Ensure that the work area is sufficiently dry. Use of a rubber dam is recommended.

Matrix placement:
 Take appropriate measures for areas where the applied **x-tra base** can flow away. Translucent matrices that are wedged in the approximal area are advantageous. Minimal separation facilitates the shaping of the approximal contact and placement of the matrix.

Pulp protection:
 A suitable pulp-protective liner that can be covered with a stable cement, if necessary, should be applied in close proximity to the pulp.

Bonding materials:
x-tra base is applied with a dentine enamel bond, using the adhesive technique. Any light-curing bonding materials can be used. Observe the relevant instructions for use during preparation (etching technique) and processing.

Application:
Syringe: Lock the enclosed intraoral tip on the syringe with a clockwise turn and directly apply **x-tra base**. In order to guarantee the functionality of the non-directing syringe do not pull back the plunger during or after use. The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. The use of a protective barrier is recommended to protect the syringe against contamination with bodily fluids, dirt from hands or oral tissue.

Caps: Remove protective cap. Insert Caps into the opening of the dispenser (following the instructions for use). Bend the metal cannula into the desired direction. Begin applying the material directly into the prepared cavity at the deepest point and fill the cavity from bottom to top. Keep the tip of the metal cannula immersed all times. To prevent the formation of air bubbles ensure a slow, even flow. If necessary, air bubbles or excess material can be removed before polymerisation. When used as a base/liner the distance to the occlusal margin of the cavity must be at least 2 mm. **x-tra base** can be polymerised in increments of up to 4 mm.

Light-curing:
 Conventional polymerisation devices are suitable for light-curing the material. Ensure that all areas of the filling are sufficiently light-cured. Additional light-curing of the accessible approximal areas and from the oral as well as the vestibular side is recommended. If metal matrices are used, these should be removed prior to additional light-curing. Hold the light emission tip of the device as close as possible to the surface of the filling. Otherwise, the curing depth may be reduced. Incomplete curing may lead to discolouration of the material, loss of stability of the composite filling and discomfort.

Depending on the light output and the chosen shade the curing times are:

Shade	universal	A2
Licht output		
500 - 800 mW/cm ²	10 s	40 s
> 800 mW/cm ²	10 s	20 s

Final occlusal layer:

When used as a base / liner **x-tra base** has to be overlaid with an at least 2 mm thick surface layer of a methacrylate-based universal / posterior composite. For working and finishing directions follow the manufacturer's instructions for use.

Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, lead to curing failures in restorative composites. Consequently, the use of zinc oxide/eugenol cements or other materials containing eugenol should be avoided in combination with restorative composites.
- Each **x-tra base** Caps unit should be used for one patient only and should be applied exerting uniform but not excessive pressure.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, BisEMA, TCDDMA, ytterbium trifluorid, fumed silica, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C-23°C. Seal the syringes immediately after removing the material in order to avoid exposure to light and subsequent polymerisation. Do not use after the expiry date.

Disposal:
 Dispose of the product in accordance with local regulations.

Reporting obligation:
 Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **x-tra base** must be reported to VOCO GmbH and the responsible authority.

Note:
 The Summary of Safety and Clinical Performance of **x-tra base** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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VC 60 BB1794 E1 0822 99 © by VOCO