

Your practice is *our* inspiration.™

Premise | Strength and beauty you can trust



Premise

Evolving beyond nano with trimodal technology




Building on an excellent, clinically proven reputation, the top-rated Premise™ universal trimodal composite is newly enhanced with an optimized filler blend, with all the strength and beauty you have come to trust. This improved combination of fillers in Kerr’s proprietary filler technology has increased strength and radiopacity and made handling even better.

*Inspired by:
The clinical need to
easily create lifelike
restorations that last
was the inspiration
behind Premise.*

Research

Kerr’s patented, proprietary trimodal technology is comprised of three different fillers—prepolymerized filler (PPF), patented Point 4 filler and 0.04 micron nanofiller—that work in concert to meet all the demands of today’s practitioner: higher compressive strength along with increased durability, easy placement, and superior polishability for beautiful esthetic restorations.




Source: Kerr Dental Materials Center (Watts method).
Data available upon request.

Features

Improved trimodal filler technology. Same great strength and performance with smoother, more consistent handling.

Ultra-low shrinkage. 84% filler load minimizes shrinkage to 1.6%—up to 45% less than other brands.

Easy polish. Clinically proven, quick, sustained polish and luster that lasts over time.




24-hour testing conducted by Kerr R&D. Data available upon request.



Trimodal Filler Technology

Three distinct fillers, including a nanofiller, work in harmony to meet the demands of today's discerning clinicians: low shrinkage, easy handling, and quick, sustained polish.




Prepolymerized Filler (PPF) – Large enough to increase loading and already preshrunk, our PPF is a proprietary blend of low-shrinkage resin, barium glass nanoparticles, so plucking becomes a nonissue.

- Decreased shrinkage
- Optimal handling

Point 4 Filler Technology

(barium glass filler of 0.4 μm average size)

- Clinically proven durability
- Clinically proven polish



Silica Nanofiller (.04 μm nanoparticles)

- Enhanced polish
- Decreased shrinkage
- Optimal handling

Premise Clinical Images

Anterior case before



Completed anterior case



© 2003 Montage Media Corporation. Case images courtesy of Dr. Douglas Terry.

Posterior case before amalgam removal



Completed posterior case



Case images courtesy of Dr. Stephen Poss.

Your practice is *our* inspiration.™



Ordering

32611 Master Syringe Kit*

Contains: 1 Syringe (4 g) of each shade: A1, A2, A3, A3.5, B1, C2, D2, XL1, XL2, Clear, Amber, Gray, A3 Opaque, A3.5 Opaque, A3 Packable and A3.5 Packable

32612 Master Unidose Kit*

Contains: 10 Unidose Tips (0.2 g each) of each shade: A1, A2, A3, A3.5, B1, C2, D2, XL1, XL2, Clear, Amber, Gray, A3 Opaque, A3.5 Opaque, A3 Packable and A3.5 Packable

*Accessory Kit (included with all Master Kits)

Contains: 1 Bottle OptiBond® Solo™ Plus (5 ml), 1 Bottle OptiGuard® (5 ml), 1 Syringe Kerr Gel Etchant (3 g), 10 Disposable Gel Etchant Syringe Tips, 25 Disposable Mixing Wells, 50 Kerr Applicators, 1 Premise Shade Guide, 1 Kerr Hawe Sample Polishing Kit

Refill Part Numbers

Body Shades – Syringe

32617 A1
32618 A2
32619 A3
32620 A3.5
32621 A4
32622 B1
32623 B2
32624 B3
32625 B4
32626 C1
32627 C2
32628 C3
32629 C4
32630 D2
32631 D3
32632 D4
32633 XL1
32634 XL2

Opaque Shades – Syringe

32635 A2
32636 A3
32637 A3.5
32638 A4
32639 B1
32640 B2
32641 C2
32642 D2

Translucent Shades – Syringe

32643 Amber
32644 Gray
32645 Clear
32646 Super Clear

32870 A2
32647 A3
32648 A3.5
32649 A4

Body Shades – Unidose

32650 A1
32651 A2
32652 A3
32653 A3.5
32654 A4
32655 B1
32656 B2
32657 B3
32658 B4
32659 C1
32660 C2
32661 C3
32662 C4
32663 D2
32664 D3
32665 D4
32666 XL1
32667 XL2

Opaque Shades – Unidose

32668 A2
32669 A3
32670 A3.5
32671 A4
32672 B1
32673 B2
32674 C2
32675 D2

Translucent Shades – Unidose

32676 Amber
32677 Gray
32678 Clear
32679 Super Clear

Packable Shades – Unidose

32871 A2
32680 A3
32681 A3.5
32682 A4

Accessories

21782 Unidose Gun/Dispenser Large, Gray
31569 Unidose Gun/Dispenser Small, White
32555 Premise Shade Guide

33371 Premise Flowable Assortment Kit

Contains: 4 Syringes (1.7 g of each shade: A1, A2, A3 and B1) and 40 Tips

Premise Flowable 4-pack Shade Kits

Contains: 4 Syringes (1.7 g of each single shade) and 40 Tips

33372 A1
33373 A2
33374 A3
33375 A3.5
33376 B1
33385 B2
33377 C2
33378 Univ. Opaque
33379 XL1
33380 XL2

Accessories

33370 Tip Refills (pack of 50)