

**Veneer cement so good
we don't even need to
show pictures of pretty
people smiling.**

FEATURES & BENEFITS



UNDETECTABLE SHIFTING

The restoration visually stays the same shade as on the first day you cured the cement.



VERSATILE SHADES

Different shades allow you to warm or brighten the restoration easily for a natural look; shades can be mixed and customized.



LIGHT CURE

No hand mixing, compatible with LED and Halogen curing lights.



EASY STORAGE

No refrigeration necessary for the Mojo cement and Try-in gels.



EASY CLEAN-UP

Simply light cure for 1-2 seconds, and remove excess cement.



VERSATILE USAGE

Designed for use with porcelain, ceramic and composite veneer restorations.



CLEAR SHADE

This truly clear cement will not alter the shade of the restoration or underlying tooth structure. Over 80% translucency!



UNIVERSAL VISCOSITY

Mojo Cement is designed for optimal consistency.



The try-in went smoothly and the cementation even better. Her dreams really did come true that day, and she was smiling inside and out. I gave her a thank you card for letting me be a part of this incredible transition. Inside, I slipped a gift card for her to go get a

beautiful white blouse. I love it whenever she comes in to the office, she can't stop smiling. Her veneers look as beautiful as the day they were cemented.

— DR. GREGORY GILLESPIE




Try-in seating

Place cement into restoration

Easy clean-up

Final result


4



5



6



7



Order today!

N06K - Mojo Veneer Cement Kit

- 8 Syringes: Cement (2.3gm) & Try-in (1.6gm) 1 each Shade: Clear, Dark, Light, Movie Star White
- 1 4ml bottle Bond-1™ Primer/Adhesive
- 1 1ml syringe 37% Phosphoric Acid Etching Gel
- 1 3ml bottle Silane
- 40 Elephant Tips, Instructions

Mojo Cement Refills

- 1 2.3gm syringe Cement (all refills include 5 tips per syringe)
- N06A** Mojo Light Cement
- N06B** Mojo Dark Cement
- N06C** Mojo Clear Cement
- N06D** Mojo Movie Star Cement

Mojo Try-in Gel Refills

- 1 1.6gm syringe Try-in Gel (all refills include 5 tips per syringe)
- N06AA** Mojo Light Try-in Gel
- N06AB** Mojo Dark Try-in Gel
- N06AC** Mojo Clear Try-in Gel
- N06AD** Mojo Movie Star Try-in Gel

Accessories

- N06E** 100 Elephant Tips



Make the call!

Want Mojo Cement?

Give your Pentron Sales Team a ring at:

855-748-2964

Have questions?

Chat with our Customer Service at:

800-551-0283

or visit:

www.pentron.com

Water absorption/solubility data:
Property Evaluation of Several Light Curable Veneer Cements
S. Jin, W. Jia, and J. Turcotte

Shade shift/flex strength/flex modulus data:
Mechanical Behavior of a New Dental Resin Veneer Cement
Jeffery Thompson, Ph.D.

Step-by-step case study courtesy of Greg Gillespie, DDS.

All trademarks are property of their respective owners.
Copyright © 2016 Pentron Clinical MKT - 1663

PENTRON™

1717 West Collins Avenue
Orange, CA 92867
800.551.0283
855.748.2964
www.pentron.com