


X-Prototype : F1, F2, F3, F4, F5



SIZE	LENGTH (mm)	D2= tip 3mm (Radius)	D3=tip16mm (Radius)	Color	TOLERANCE (mm)
#F1	27 ± 1	0.38	0.80	Yellow	± 0.5
#F2	27 ± 1	0.43	1.00	Red	± 0.5
#F3	27 ± 1	0.52	1.10	Blue	± 0.5
#F4	27 ± 1	0.55	1.10	Black	± 0.5
#F5	27 ± 1	0.65	1.00	Yellow	± 0.5